




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Surgical Gown

Reference(s)

Reference	Length (cm)	Size	Finger loop	Mask	Hand towels	Unit/case
S3505CEC	115	S/M	No	No	2	36
S3506CEC	120	L	No	No	2	36
S3507CEC	140	XL	Yes	No	2	32
S3516CEC	150	XLL	No	No	2	32
S3545CEC	130	XL	No	No	2	36
FS3545CEA	140	XL	No	No	2	32
S3585CEC	170	2XLL	No	No	2	32

Description

A SMS wrap 60x60 cm contains:

- 1 SMMMS sterile back gown, hook and loop adjustable neck
- 2 absorbent hand towels, 39x30 cm unfolded

Main materials

SMMMS (polypropylene)

Packaging

Paper and film peel pouch
2 polyliners contain all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
In compliance with the standard EN13795
MDD Class: I sterile
AAMI level : 2



Latex free



0086



STERILE EO



© 2018 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_GownEssentialStandard01_EN03

01.INT.03.30