

Convertors®

Medical Products and Services

Standard surgical gown

References	Length cm	Size	Pack factor	Hand towels
S3506	120	L	28	2
S3545	140	XL	28	2
S3516	150	XLL	24	2

Description: A sterile wrap 60 cm x 60 cm contains:
1 sterile back gown,
Velcro® adjustable neck
2 hand towels


Materials: SMS (polypropylene) 35g/m²

Packaging: Vented sterilization pouch for the units
2 polyliners contain the units
Carton: 60 cm x 40 cm x 27.5 cm


Manufacturer: Cardinal Health France 205 S.A.S.
Manufacturing site: EN ISO 13485 certified by TÜV Product Service

Sterilisation: Ethylene oxide

Regulation: Compliant with Council Directive 93/42/EEC
Class: I sterile
Standard performance according to EN 13795 standard



ORGANIZATION OF THE BOX (Top View)



0123

© Copyright 2007 Cardinal Health, Inc. or one of its subsidiaries.
All rights reserved.
LIT07088ENG-ED03

Convertors® is a registered trademark of Cardinal Health, Inc. or one of its subsidiaries.

Cardinal Health France 205 S.A.S.
5 rue Charles Lindbergh
F-44110 Châteaubriant
Tél. : +33 (0)2 40 81 56 31
Fax : +33 (0)2 40 81 56 34

www.cardinal.com

