




Medline International France SAS

5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com/en/](http://www.medline.com/en/)



## OPS™ Ultimate

## Thyroid drape pack, two parts

### Reference(s)

**TB277041CE**

### Description

A back table cover 140x190 cm, 50 µm contains:

- 1 tri-laminate Head split sheet 240x145 cm, 20x70 cm adhesive split, 17x34 cm reinforcement
- 1 tri-laminate Foot split sheet, 183x240 cm, 20x25 cm adhesive split, 17x34 cm reinforcement, 1 hook and loop tube holder
- 1 Mayo stand cover, 80x142 cm

### Main materials

Absorbent and impervious tri-laminated material (polyethylene + polypropylene)  
HSPP absorbent reinforcement (polypropylene + polyethylene)

### Pack factor

6 units per carton

### Packaging

Vented sterilization pouch  
1 polyliner contains all the units in a carton box

### Manufacturer

Medline International France SAS  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile  
HP: High performance according to EN 13795 standard



Latex free



0086



STERILE EO



© 2014 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

LIT-TB277041CE-EN02

01.01.INT.12