

Technical sheet

By Convertors®

Laparoscopy drape

Reference(s)	TB29496CE
Description	<p>A sterile wrap 112 x 112 cm contains:</p> <p>A table cover (HP) 140 x 190 cm which contains</p> <p>2 hand towels</p> <p>1 laparoscopy drape Tiburon® (HP) 255 x 282 x 315 cm, reinforcement HSPP</p> <ul style="list-style-type: none"> - abdominal fenestration adhesive 34 x 38 cm with incise film - perineal fenestration adhesive 10 x 14 cm - 2 attached instrument pouches 40 x 46 cm - Attached leggings
Main materials	<p>Tiburon® (Polypropylene + polyolefin film)</p> <p>HSPP (Absorbent spunbond polypropylene)</p>
Pack factor	5 packs per carton
Packaging	<p>Sterilization pouch</p> <p>1 poly-liner contains all the units</p> <p>Carton</p>
Manufacturer	<p>CareFusion Switzerland 317 Sàrl</p> <p>EN ISO 13485 certified by BSI</p>
Sterilisation	Ethylene oxide
Regulation	<p>In compliance with Medical Devices Directive 93/42/EEC</p> <p>MDD Class: I sterile</p> <p>HP: High performance according to EN 13795 standard</p>

