




Medline International Germany GmbH  
Wilhelm-Sinsteden-Str. 5-7  
47533 Kleve  
www.medline.com/en



## Technical sheet

# Kidney Basin, bulk non-sterile

### Reference(s)

9005-EU

### Description

Kidney Basin shaped, clear  
Volume: 800 ml

### Main materials

Polypropylene

### Pack factor

288 units per carton

### Packaging

PE bag  
Carton box

### Manufacturer

Medline International France SAS  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I non-sterile  
(to be sterilized using an approved sterilization process)



© 2015 Medline Industries Inc. or one of his subsidiaries.  
All rights reserved.

LIT-BNSEmesis9005-EU-EN01

01.01.INT.12