





Angiography Drape

OPST™ ESSENTIAL



Reference	ESAG3541CE					
Description	<p>A SMS wrap contains:</p> <p>- 1 Bi-laminate Angiography drape 244x370 cm, 112/244x145cm reinforcement, 2 circular fenestrations diameter 15 cm with fenestrated incise film diameter 9 cm and 2 radial fenestrations 9x12cm with fenestrated incise film (open arms), 2 clear PE panels*</p>					
Main Materials	<p>Absorbent and impervious bi-lamited material (polyethylene + polypropylene)</p> <p>Polypropylene reinforcement</p> <p>Clear polyethylene</p>					
Pack Factor	10 Units per case					
Packaging	All the units are packaged in a case					
Logistic Information	Packaging Unit	GTIN	Length(mm)	Width(mm)	Height(mm)	Weight(gr)
	Case	50888277428841	600	400	278	9400
	Each	00888277428846	-	-	-	940
Legal Manufacturer	<p>Medline International France Sas</p> <p>In compliance with EN ISO 13485</p>					
Regulation	<p>In compliance with Medical Devices Regulation 2017/745 and UK MDR</p> <p>MD Class: I sterile</p> <p>*High Performance per EN 13795-1</p> <p>Please refer to complete labeling before using the product</p>					
Other Information	NA					



Medline Industries Ltd
 3rd Floor, Quayside Wilderspool Business Park
 Greenalls Avenue, Warrington WA4 6HL
 United Kingdom

Tel.: +44 844 334 5237
 Fax: +44 844 334 5238
 www.medline.eu/uk