




Emesis/Sickness Bag

Blue, twist-and-seal



Reference	EMB01B					
Description	Emesis bag used to collect vomit from patients. The bag is plastic and it is easy to handle. Can be easily closed and locked with the plastic ring for safe and clean disposal.					
Main Materials	Film: 100% blue polyethylene Ring: 100% polypropylene					
Pack Factor	500 Units per case					
Packaging	All the units are packaged in a case					
Logistic Information	Packaging Unit	GTIN	Length(mm)	Width(mm)	Height(mm)	Weight(gr)
	Case	40888277447494	675	350	296	8540
	Box	30888277447497	132	334	134	780
Legal Manufacturer	Medline International France Sas In compliance with EN ISO 13485					
Regulation	In compliance with Medical Devices Regulation 2017/745 and UK MDR MD Class: I non-sterile Please refer to complete labeling before using the product					
Other Information	Bag physical properties: Total capacity 2L with graduated volume of 1L (approximate volume) Number of closing points: 2 Disposal: Please refer to the hospital policy for waste disposal.					



Medline Industries Ltd
3rd Floor, Quayside Wilderspool Business Park
Greenalls Avenue, Warrington WA4 6HL
United Kingdom

Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk