




# Emesis Bag, Blue, With Absorbent Pad



<b>Reference</b>	<b>EMB03B</b>					
<b>Description</b>	Emesis bag used to collect vomit from patients. The bag is plastic and it is easy to handle. Can be easily closed and locked with the plastic ring for safe and clean disposal. Absorbent pad inside can transfer the liquid into solid to avoid spill.					
<b>Main Materials</b>	Film: 100% blue polyethylene Ring: 100% polypropylene Absorbent pad: pulp, cellulose, super absorbent polymer					
<b>Pack Factor</b>	300 Units per case					
<b>Packaging</b>	All the units are packaged in a case					
<b>Logistic Information</b>	<b>Packaging Unit</b>	<b>GTIN</b>	<b>Length(mm)</b>	<b>Width(mm)</b>	<b>Height(mm)</b>	<b>Weight(gr)</b>
	Case	40888277449559	416	284	462	10280
	Each	00888277449551	-	-	-	34
	Box	30888277449552	446	134	132	1580
<b>Legal Manufacturer</b>	Medline International France Sas In compliance with EN ISO 13485					
<b>Regulation</b>	In compliance with Medical Devices Regulation 2017/745 and UK MDR MD Class: I non-sterile Please refer to complete labeling before using the product					
<b>Other Information</b>	Bag physical properties: Total capacity 2L with graduated volume of 1L (approximate volume) and absorbent pad capacity 219 ml Number of closing points: 1.					



Medline Industries Ltd  
 3rd Floor, Quayside Wilderspool Business Park  
 Greenalls Avenue, Warrington WA4 6HL  
 United Kingdom

Tel.: +44 844 334 5237  
 Fax: +44 844 334 5238  
[www.medline.eu/uk](http://www.medline.eu/uk)