





# Hip drape for anterior approach – Right hip

## OPST™ ESSENTIAL



Reference	ES29434LCE					
Description	A SMS wrap, 90x90 cm, contains: -1 Bi-laminate drape, 245x310cm, reinforcement 105x150cm, 11x115 cm adhesive split , with integrated stockinette 110cm for left leg*					
Main Materials	Absorbent and impervious bi-laminated material (polyethylene + polypropylene) Absorbent spunbond polypropylene reinforcement					
Pack Factor	12 Units per case					
Packaging	All the units are packaged in a case					
Logistic Information	Packaging Unit	GTIN	Length(mm)	Width(mm)	Height(mm)	Weight(gr)
	Case	50888277433470	600	400	250	9660
	Each	00888277433475	-	-	-	805
Legal Manufacturer	Medline International France Sas In compliance with EN ISO 13485					
Regulation	In compliance with Medical Devices Regulation 2017/745 and UK MDR MD Class: I sterile *High Performance per EN 13795-1 Please refer to complete labeling before using the product					
Other Information	NA					



Medline Industries Ltd  
3rd Floor, Quayside Wilderspool Business Park  
Greenalls Avenue, Warrington WA4 6HL  
United Kingdom

Tel.: +44 844 334 5237  
Fax: +44 844 334 5238  
www.medline.eu/uk