




Instrument Cover  
Invisishield™



|                      |                                                                                                                                                                              |                |            |           |            |            |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|------------|-----------|------------|------------|
| Reference            | ICE2158                                                                                                                                                                      |                |            |           |            |            |
| Description          | Clear PE cover, 61x46cm, adhesive 1.3x61cm, 50 µm*                                                                                                                           |                |            |           |            |            |
| Main Materials       | Clear PE film, 50 µm                                                                                                                                                         |                |            |           |            |            |
| Pack Factor          | 40 Units per case                                                                                                                                                            |                |            |           |            |            |
| Packaging            | All the units are packaged in a case                                                                                                                                         |                |            |           |            |            |
| Logistic Information | Packaging Unit                                                                                                                                                               | GTIN           | Length(mm) | Width(mm) | Height(mm) | Weight(gr) |
|                      | Case                                                                                                                                                                         | 50888277422559 | 203        | 203       | 171        | 1354       |
|                      | Each                                                                                                                                                                         | 00888277422554 | 59         | 59        | 50         | 34         |
|                      | Box                                                                                                                                                                          | 20888277422558 | 128        | 128       | 108        | 339        |
| Legal Manufacturer   | Medline International France Sas<br>In compliance with EN ISO 13485                                                                                                          |                |            |           |            |            |
| Regulation           | In compliance with Medical Devices Regulation 2017/745 and UK MDR<br>MD Class: I sterile<br>*HP per EN 13795-1<br>Please refer to complete labeling before using the product |                |            |           |            |            |
| Other Information    | NA                                                                                                                                                                           |                |            |           |            |            |



Medline Industries Ltd  
3rd Floor, Quayside Wilderspool Business Park  
Greenalls Avenue, Warrington WA4 6HL  
United Kingdom

Tel.: +44 844 334 5237  
Fax: +44 844 334 5238  
www.medline.eu/uk