



Medline Industries Ltd
 No. 5 Booths Park
 Chelford Road Knutsford
 WA16 8GS United Kingdom

www.medline.com/en/



Technical sheet

ULTRASORBS AP DISPOSABLE UNDERPAD

Reference(s)

ULTRASRBE6191, ULTRASRBE4561, ULTRASRBE2540, USAPE6191ES, USAPE6191UF, MSCE4090PUPS, MSCE4090PUPSB

Description

AIR PERMEABLE				
Product No.	Size	Backing Print Color	Case Quantity	Absorbency
ULTRASRBE6191	61x91cm	Blue	70 in 7 bags of 10	1800-2300g
ULTRASRBE4561	45x61cm	Blue	60 in 6 bags of 10	900-1150g
ULTRASRBE2540	25x40cm	Blue	100 in 10 bags of 10	200-300g
USAPE6191ES	61x91cm	Green	60 in 6 bags of 10	1800-2300g
USAPE6191UF	61x91cm	Blue	70 in 6 bags of 10	1800-2300g
MSCE4090PUPS	101 x 228 cm	Blue	30 in 1 bag of 30	4000-4500g
MSCE4090PUPSB	101 x 228 cm	Blue	30 in 30 bags of 1	4000-4500g

- Disposable non-sterile underpads, to be used to absorb fluids, including with patients having burns.

Main materials

Spun-lace polypropylene, Airlaid core with Super absorbent polymer, Air permeable Polyethylene/nonwoven Backsheet

Pack factor

Refer to description item

Packaging

Plastic bags in corrugate
 General Warehouse Conditions, Keep Dry

Manufacturer

Medline Industries Inc, 1 Medline Place, Mundelein, IL 60060 USA
 EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
 MDD Class: I Non Sterile

