



ULTRASORBS AP DISPOSABLE UNDERPAD TECHNICAL DATA

AIR PERMEABLE						
Product No.	Size	Backing Print Color	Packaging	Case Quantity	Absorbency	Box Dimensions
ULTRASRBE6191	60 x 91 cm	Blue	Poly Bag/Corrugate	70	1800-2300g	18x9x18in
ULTRASRBE2540	25 x 40 cm	Blue	Poly Bag/Corrugate	100	200-300g	17.3x10.6x5.5in
ULTRASRBE4561	45 x 61 cm	Blue	Poly Bag/Corrugate	60	900-1150g	17.7x9.8x7.9in
USAPE6191UF	60 x 91 cm	Blue	Poly Bag/Corrugate	70	1800-2300g	18x9x18in
USAPE6191ES	60 x 91 cm	Green	Poly Bag/Corrugate	60	1800-2300g	18x9x18in
ULTRASRBE2436 (North America/ Australia/New Zealand)	24 x 36 in	Blue	Poly Bag/Corrugate	70	1800-2300g	18x9x18in
NON-AIR PERMEABLE						
MSCE282100	60 x 91 cm	N/A	Poly Bag/Corrugate	70	1800-2300g	18x9x18in

- Description :** A package contains : Disposable Underpad
- Material :** Spun-lace polypropylene, Airlaid core with Super absorbent polymer, Air permeable Polyethylene/nonwoven Backsheet
- Storage Conditions:** General Warehouse Conditions, Keep Dry
- Packaging:** Bags of 10 Disposable Underpads in Corrugate
- Sterilization :** NON-STERILE
- Manufacturer :** Conforms to ISO 13485
- Law :** Conforms to EU Directive 93/42/EEC - Class I



Medline Industries, Inc. is committed to reducing environmental impact through increased use of post-consumer recycled materials.



Medline Industries, Inc
The Company of Choice
www.medline.com

Medline
One MedlinePlace
Mundelein, Illinois 60060 USA
Telephone +1 847 643 3292
Fax +1 847 643 3295
International@medline.com