



Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en



Technical Sheet

Suresite® Select

Reference(s)

| References | Width (cm) | Length (cm) | Recommended sites |
|-------------|------------|-------------|---|
| VAD1680-INT | 4,4 | 3,8 | Peripheral lines, PICC lines, Midlines, Central Venous Catheters (CVC), Tunneled CVC, Implantable ports, Arterial catheters, Dialysis catheters |
| VAD1682-INT | 5,7 | 5,1 | |
| VAD1882-INT | 7,0 | 8,9 | |
| VAD1683-INT | 7 | 7,3 | |
| VAD1685-INT | 8,9 | 11,4 | |
| VAD1688-INT | 10,2 | 11,4 | |
| VAD1689-INT | 10,2 | 15,5 | |

Description

Suresite® Select is a highly permeable, transparent film dressing that uses an innovative multi-layer construction.
Transparent film provides a breathable, sterile barrier to external contaminants.
Innovative adhesives hold securely and manage moisture.
Frame delivery ensures accurate and gentle removal.
Clover notch fits easily and confirms around hubs and lumens.
Wide viewing window makes it easy to monitor for signs of infection.

Main materials

Liner: Silicone treated paper
Adhesive: Acrylic based (biocompatible, skin contact grade)
Film: Polyurethane film, 0,02mm.
Fabric: Non-woven fabric, 0,5mm
Strips: Cellulose and Acrylic adhesive.

Pack factor

VAD1680-INT : 400 units per carton, 4 boxes of 100 units
VAD1682-INT : 400 units per carton, 4 boxes of 100 units
VAD1882-INT : 200 units per carton, 4 boxes of 50 units
VAD1683-INT : 400 units per carton, 4 boxes of 100 units
VAD1685-INT : 200 units per carton, 4 boxes of 50 units
VAD1688-INT : 200 units per carton, 4 boxes of 50 units
VAD1689-INT : 100 units per carton, 4 boxes of 25 units

Packaging

Pouch
A carton box contains all the units

Manufacturer

Centurion Medical Product Inc (a Medline company)
EN ISO 13485



Latex free



0086



STERILE EO

© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS-SuresiteSelect-EN02

01.INT.03.30

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile

Shelf Life

5 years



STERILE EO

© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS-SuresiteSelect-EN02

01.INT.03.30