

# Sahara Protective OR Table Sheet


Promote safety. Improve efficiency. Reduce costs.



The Sahara Protective OR Table Sheet is designed for short, low-fluid surgical procedures. The back sheet guarantees effective protection of the operating theatre table and reduces the risk of cross-contamination because it is made of a durable, impervious polyethylene (PE) film.

## Key benefits and features

- Ensures an enhanced surface grip, thanks to the impervious, non-slippery back sheet
- Reduces the risk of cross-contamination because it is single-use
- Eliminates unnecessary linen laundry costs
- Minimises the risk of allergic reactions because it is latex-free
- **Low linting, according to EN ISO 9073-10**



### Outer layer

- Soft, non-woven

### Bottom layer

- Impervious PE backing provides barrier protection and helps prevent fluids from compromising operating table pads over time

## Ordering information

Item number	Description	Quantity
DYND4090NA	Sahara Protective 229 x 102 cm	40/case 10/bag

**For more information** on this product, please contact your Medline account manager or visit our website: [uk.medline.eu](http://uk.medline.eu)