



SKIN HEALTH SOLUTION

BORDERED GAUZE

Island Dressing with Adhesive Border

Product overview

Bordered Gauze, adhesive island dressing with water resistant backing is ideal for the treatment of smaller wounds and incision sites

- It offers a full range of sizes to simplify dressing choice
- Up to 7 days wear time
- Secure fixation, conformability, comfort and absorption properties



Item Details

Product	Item No.	Size	Order Qty	PIPP/ FP10 Code	Bordered Gauze NHSSC code
Dressing Adhesive Island Non-Woven Sterile	MSCE3223	6x7 cm	1 box of 15 dressings	405-1116	EIJ382
	MSCE3278	6.5x8 cm	1 box of 15 dressings	405-1124	EIJ388
	MSCE3244	10x10 cm	1 box of 15 dressings	405-1132	EIJ383
	MSCE3246	10x15 cm	1 box of 15 dressings	405-1140	EIJ384
	MSCE3248	10x20 cm	1 box of 15 dressings	405-1157	EIJ385
	MSCE32410	10x25 cm	1 box of 15 dressings	405-1165	EIJ386
	MSCE32414	10x35 cm	1 box of 15 dressings	405-1173	EIJ387

Indications

- Low to moderately exuding wounds and incision sites

Contact details

Customer Services

Tel.: +44 844 334 5237
uk-customerservice@medline.com

Medline Industries Ltd

3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Fax: +44 844 334 5238
www.medline.eu/uk

FEATURES

- Full range of products and sizes to simplify product choice
- Cost effective solution for dressing low to moderately exuding wounds and incision sites
- Deluxe, anti-adherent soaker pad to prevent maceration and support longer wear time

This product is a medical device of class I sterile. Before use, consult instructions for use and precautions on the corresponding labeling.

© 2018 Medline Industries, Inc. Medline and Medline logo are registered trademark of Medline Industries, Inc. USA.

ML585-EN02/OS 11/2018.



Shanghai International Holding Corp.
Shanghai, China

20137 Hamburg, Germany



0123 (TUV SUD)