

Central Sterile Product Guide

Comprehensive product selection and support



Sterilisation at its best





Meticulous care

Few understand the amount of work that goes into central sterile processing, which aims to eliminate 99.999% of pathogens.

Medline does not take the sterility of surgical instrumentation for granted. We understand the challenge and take a systematic approach to assembling broad lines of essentials for every step of the central sterile process.

Sterilisation

Sterilisation Sterilisation	4
Sterilisation Wraps	5
Sterilisation Trays	8
Sterilisation Wrap Accessories	8

Protective Apparel

Gloves	10
Gowns and Shoe Covers	12
Face Protection	14

Floor Safety

Floor Mat	16
-----------	----

Use this comprehensive guide to work with your Medline representative. Together, you can outfit your central sterile department with the selection, quality and value to facilitate your systematic approach.

Sterilisation

Going from individual components to assembled and sterilised sets is where the central sterile department transitions from preparing tools to delivering outstanding patient care. Medline helps you accomplish your goals with safety, sustainability and value in mind.

Gemini Sterilisation Wraps²

Drapeable and stronger

Gemini wrap is resistant to tears and punctures, performing better than standard SMS wraps. It is very soft and easy to fold, which makes it simple to work with.

Low linting and better adhesion properties

The surface is friendly to sterilisation tape. Naturally low linting, it helps maintain the integrity of sterile packs.

Made with unique SSMMS fabric with proven barrier protection

Gemini's strong multiple layers and 100% polypropylene (EN 11607-1 and EN 868-2 tested) construction provide unique high-level barrier protection against fluid and particles.

Versatile

Medline's Gemini Wrap is indicated for use in all major sterilisation cycles, including:

GEM0-GEM5

- Prevacuum steam cycles
- Gravity steam cycles

GEM1-GEM5

- Ethylene oxide (EtO) sterilisation
- STERRAD® sterilisation
- STERIS V-PRO® low-temperature sterilisation
- STERIZONE® VP4 sterilisation.

WRAP Certification

Gemini wraps are manufactured in a WRAP-certified facility. Promoting and respecting labour and human rights in our supply chain is part of Medline's corporate social responsibility agenda. More detailed descriptions of WRAP's 12 Principles are available on their website: www.wrapcompliance.org.

Gemini Single Sterilisation Wrap

Gemini Single Sterilisation Wrap is a great solution for sterilisation units where two sheets are used to wrap a medical device or an instrument set.



	GEMO 40 gr	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr	GEM5 75 gr
46 x 46 cm	GEMO118-EU 1,000 ea/cs	GEM1118-EU 1,000 ea/cs	GEM2118-EU 1,000 ea/cs	GEM3118-EU 600 ea/cs		
51 x 51 cm	GEMO120-EU 1,000 ea/cs	GEM1120-EU 1,000 ea/cs	GEM2120-EU 500 ea/cs	GEM3120-EU 250 ea/cs		
61 x 61 cm	GEMO124-EU 500 ea/cs	GEM1124-EU 500 ea/cs	GEM2124-EU 500 ea/cs	GEM3124-EU 250 ea/cs	GEM4124-EU 250 ea/cs	
76 x 76 cm	GEMO130-EU 300 ea/cs	GEM1130-EU 300 ea/cs	GEM2130-EU 300 ea/cs	GEM3130-EU 200 ea/cs	GEM4130-EU 250 ea/cs	
91 x 91 cm	GEMO136-EU 300 ea/c	GEM1136-EU 300 ea/cs	GEM2136-EU 300 ea/cs	GEM3136-EU 150 ea/cs	GEM4136-EU 150 ea/cs	GEM5136-EU 144 ea/cs
102 x 102 cm	GEMO140-EU 250 ea/cs	GEM1140-EU 250 ea/cs	GEM2140-EU 250 ea/cs	GEM3140-EU 150 ea/cs	GEM4140-EU 150 ea/cs	GEM5140-EU 96 ea/cs
114 x 114 cm	GEMO145-EU 250 ea/cs	GEM1145-EU 250 ea/cs		GEM3145-EU 100 ea/cs	GEM4145-EU 100 ea/cs	GEM5145-EU 96 ea/cs
122 x 122 cm	GEMO148-EU 250 ea/c	GEM1148-EU 250 ea/cs	GEM2148-EU 100 ea/cs	GEM3148-EU 50 ea/cs	GEM4148-EU 50 ea/cs	GEM5148-EU 48 ea/cs
137 x 137 cm	GEMO154-EU 100 ea/cs	GEM1154-EU 100 ea/cs	GEM2154-EU 100 ea/cs	GEM3154-EU 50 ea/cs	GEM4154-EU 50 ea/cs	GEM5154-EU 48 ea/cs
137 x 183 cm	GEMO172-EU 100 ea/cs			GEM3172-EU 50 ea/cs	GEM4172-EU 50 ea/cs	

Gemini Interleaved Sterilisation Wrap

Gemini Interleaved Sheets offer both our multilayer technology and colour differentiation to sterilisation units where sequential wrapping is in place. This solution combines tear resistance and security. It is indicated for use in all major sterilisation cycles.



	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr
61 x 61 cm	GEM1124INT-EU 500 ea/cs	GEM2124INT-EU 500 ea/cs	GEM3124INT-EU 250 ea/cs	GEM4124INT-EU 250 ea/cs
76 x 76 cm	GEM1130INT-EU 300 ea/cs	GEM2130INT-EU 300 ea/cs	GEM3130INT-EU 200 ea/cs	GEM4130INT-EU 250 ea/cs
91 x 91 cm	GEM1136INT-EU 300 ea/cs	GEM2136INT-EU 300 ea/cs	GEM3136INT-EU 150 ea/cs	GEM4136INT-EU 150 ea/cs
102 x 102 cm	GEM1140INT-EU 250 ea/cs	GEM2140INT-EU 250 ea/cs	GEM3140INT-EU 150 ea/cs	GEM4140INT-EU 150 ea/cs
122 x 122 cm	GEM1148INT-EU 250 ea/cs	GEM2148INT-EU 100 ea/cs	GEM3148INT-EU 50 ea/cs	GEM4148INT-EU 50 ea/cs
137 x 137 cm	GEM1154INT-EU 100 ea/cs	GEM2154INT-EU 100 ea/cs	GEM3154INT-EU 50 ea/cs	GEM4154INT-EU 50 ea/cs
137 x 183 cm				GEM4172INT-EU 50 ea/cs

Gemini Bonded Sterilisation Wrap

Gemini Bonded Wrap can give you greater confidence in your sterilisation process. Two sheets fused together provide exceptional strength and help improve efficiency. The dual colour facilitates hole detection and helps differentiate between instrument sets.



	GEMO 40 gr	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr	GEM5 75 gr
46 x 46 cm	GEMO118T-EU 500 ea/cs	GEM1118T-EU 500 ea/cs	GEM2118T-EU 500 ea/cs	GEM3118T-EU 300 ea/cs		
51 x 51 cm	GEMO120T-EU 500 ea/cs	GEM1120T-EU 500 ea/cs	GEM2120T-EU 250 ea/cs	GEM3120T-EU 125 ea/cs		
61 x 61 cm	GEMO124T-EU 250 ea/cs	GEM1124T-EU 250 ea/cs	GEM2124T-EU 250 ea/cs	GEM3124T-EU 100 ea/cs	GEM4124T-EU 100 ea/cs	
76 x 76 cm	GEMO130T-EU 150 ea/cs	GEM1130T-EU 150 ea/cs	GEM2130T-EU 150 ea/cs	GEM3130T-EU 75 ea/cs	GEM4130T-EU 100 ea/cs	
91 x 91 cm	GEMO136T-EU 150 ea/cs	GEM1136T-EU 150 ea/cs	GEM2136T-EU 150 ea/cs	GEM3136T-EU 75 ea/cs	GEM4136T-EU 75 ea/cs	GEM5136T-EU 72 ea/cs
102 x 102 cm	GEMO140T-EU 100 ea/cs	GEM1140T-EU 100 ea/cs	GEM2140T-EU 100 ea/cs	GEM3140T-EU 75 ea/cs	GEM4140T-EU 75 ea/cs	GEM5140T-EU 48 ea/cs
114 x 114 cm	GEMO145T-EU 100 ea/cs	GEM1145T-EU 100 ea/cs		GEM3145T-EU 50 ea/cs	GEM4145T-EU 50 ea/cs	GEM5145T-EU 48 ea/cs
122 x 122 cm	GEMO148T-EU 100 ea/cs	GEM1148T-EU 100 ea/cs	GEM2148T-EU 50 ea/cs	GEM3148T-EU 30 ea/cs	GEM4148T-EU 30 ea/cs	GEM5148T-EU 24 ea/cs
137 x 137 cm	GEMO154T-EU 50 ea/cs	GEM1154T-EU 50 ea/cs	GEM2154T-EU 50 ea/cs	GEM3154T-EU 30 ea/cs	GEM4154T-EU 30 ea/cs	GEM5154T-EU 24 ea/cs
137 x 183 cm				GEM3172T-EU 30 ea/cs	GEM4172T-EU 30 ea/cs	GEM5172T-EU 24 ea/cs

Wrap Accessories

Dust Cover

These covers protect instrument trays from dust and environmental agents. They are constructed from strong, quality polyethylene and are self-sealable.

Item no.	Description	Packaging
MDS100815	20 x 30 cm	1,000
MDS100825	30 x 50 cm	500
MDS100830	40 x 55 cm	500
MDS100835	40 x 75 cm	250
MDS100840	60 x 75 cm	250
MDS100845	63 x 94 cm	100



Transport Tray⁴

- Rigid plastic is designed to permit air flow
- Helps reduce exposed wire shelving
- Helps protect instrument tray from holes or tears during sterilisation, storage and transport

Item no.	Description	Packaging
GEMTRAY3CE	31 x 56 cm	10
GEMTRAY4CE	43 x 69 cm	8



Tray Liner Sheet³

Paper tray liners facilitate the absorption of condensation during sterilisation.

Item no.	Description	Packaging
TLC009	30 x 30 cm	2,000
TLC010	30 x 45 cm	1,000
TLC011	30 x 60 cm	1,000



Sterilisation Trays

Wire Mesh Stainless Steel Trays

Select the perfect size for the majority of your instrument sets.

For half-sized containers

Item no.	Description	Packaging
SCD 333	24.5 x 24.5 x 3.0 cm	1/ea
SCD 336	24.5 x 24.5 x 6.0 cm	1/ea
SCD 3310	24.5 x 24.5 x 10.0 cm	1/ea

For three-quarter containers

Item no.	Description	Packaging
SCD 433	40.0 x 25.0 x 3.0 cm	1/ea
SCD 436	40.0 x 25.0 x 6.0 cm	1/ea
SCD 4310	40.0 x 25.0 x 10.0 cm	1/ea

For full-sized containers

Item no.	Description	Packaging
SCD 633	51.5 x 24.5 x 3.0 cm	1/ea
SCD 636	51.5 x 24.5 x 6.0 cm	1/ea
SCD 6310	51.5 x 24.5 x 10.0 cm	1/ea
SCDD 636	48.5 x 24.5 x 6.0 cm	1/ea




Sterilisation Accessories

Silicone Mats

Cushion and protect delicate surgical instruments with these silicone mats; sized to fit perfectly into Medline trays.

Item no.	Description	Packaging
MDNSNS315	12.7 x 24.1 cm	1/ea
MDNSNS335	24.1 x 24.1 cm	1/ea
MDNSNS435	39.4 x 24.1 cm	1/ea
MDNSNS635	47 x 24.1 cm	1/ea



Protective Apparel

There is no shortage of hazardous substances and microorganisms in the healthcare environment. Provide an effective barrier to these substances and the highest level of safety with Medline's protective apparel. Each item is designed with protection and comfort in mind.



Exam Gloves⁵

SensiCare® Pro 400 mm Extended Cuff Nitrile Exam Gloves

With **OptiCuff**
anti-cuff-
roll-down
technology

Our SensiCare Pro Nitrile Gloves provide exceptional all-round protection, fit and performance. These gloves are designed for staff who use harsh chemicals for general tasks and procedures; clean, decontaminate and/or disinfect surgical instruments.

- Longer 400 mm cuffs provide extra protection between cuffs and sleeves
- Fully textured grip enhances tactile sensitivity
- Good performance in puncture resistance tests
- Chemo-tested
- PPE-certified for extra protection



Item no.	Description	Packaging
PRO400XS	SensiCare Pro 400 mm nitrile exam gloves, x-small	50/bx, 10 bx/cs
PRO400S	SensiCare Pro 400 mm nitrile exam gloves, small	50/bx, 10 bx/cs
PRO400M	SensiCare Pro 400 mm nitrile exam gloves, medium	50/bx, 10 bx/cs
PRO400L	SensiCare Pro 400 mm nitrile exam gloves, large	50/bx, 10 bx/cs
PRO400XL	SensiCare Pro 400 mm nitrile exam gloves, x-large	50/bx, 10 bx/cs

EN ISO 374-5:2016	Level
Protection against bacteria and fungi	Pass
Protection against viruses	Pass



SensiCare® Extended Cuff Nitrile Exam Gloves with excellent tactile sensitivity

- Durable nitrile protection
- Extended cuff for additional protection
- PPE-certified for extra coverage
- Comfort of latex in a nitrile glove


Item no.	Description	Packaging
MDSEXTXS	SensiCare Extended Cuff Nitrile Exam Gloves, x-small	100/bx, 10 bx/cs
MDSEXTS	SensiCare Extended Cuff Nitrile Exam Gloves, small	100/bx, 10 bx/cs
MDSEXTM	SensiCare Extended Cuff Nitrile Exam Gloves, medium	100/bx, 10 bx/cs
MDSEXTL	SensiCare Extended Cuff Nitrile Exam Gloves, large	100/bx, 10 bx/cs
MDSEXTXL	SensiCare Extended Cuff Nitrile Exam Gloves, x-large	90/bx, 10 bx/cs



Built-in impervious cover gowns compliance


Your staff knows what they need to do to protect themselves and their patients, yet time constraints and inconvenience can lead to cutting corners with proper gown compliance. Medline's protective apparel is quick and easy to don and remove while offering protection and comfort.

Easy donning. Exceptional protection. Breathable comfort.



Donning

1. Begin by placing your arms in the sleeves with thumbs in the thumb loops.
2. Pull the gown over your head.
3. Tie in back at the waist.



Removal

1. Pull at the waist to break ties.
2. Pull gown at chest with arms crossed.
3. Continue to roll gown down the body into a ball and discard after removal.

With no neck ties to undo, you reduce the chance of cross-contamination from gloves to the hair or back of the neck.

Impervious cover gowns

Suitable for protecting the user against prolonged contact with water and cleansing agents; CRIE404xC range and CRIE500x range tested according to EN 16523-1 against cleansing agent glucoprotamin (e.g. Incidin Plus).

Lightweight bilaminate impervious cover gown

- Waist and neck ties
- Knit cuffs

Item no.	Size	Colour	Packaging
CRIE4040C	L	Yellow	10/bag, 100/case
CRIE4041C	XL	Yellow	10/bag, 100/case



CRIE4040C

Polyethylene impervious cover gowns

- Over-the-head design
- Thumb-loop-style wrists
- Individually folded and wrapped in a dispenser box

Item no.	Size	Colour	Packaging
CRIE5010	L	Blue	1/bag, 15/box, 75/case
CRIE5011	XL	Blue	1/bag, 15/box, 75/case



CRIE5010

Polyethylene disposable apron

- Worn over non-fluid-resistant gowns or scrub suits as a convenient solution for users who want to add extra fluid protection at a lower cost

Item no.	Description	Packaging
NONE24274W	Medium weight, 71.1 x 116.8 cm, individually wrapped, dispenser box	50/box, 250/case



NONE24274W

Footwear⁶

Shoe covers

Footwear covers prevent environmental contamination and reduce the transmission of microorganisms. Manufactured from polyethylene or polypropylene materials, some models feature additional fluid protection.

Item no.	Size	Material	Characteristics	Packaging
CRIE2010	One size fits most	Polyethylene film	Impervious	100/bag, 1,000/case
CRIE2020	One size fits most	Polyethylene film and polypropylene (SPP)	Impervious	100/bag, 200/case



CRIE2010



CRIE2020

Boot covers

Medline's boot cover range offers enhanced coverage and increased fluid protection. The products cover the shoe and leg all the way up to the knee and have a convenient one-size-fits-most design, with a textured sole on the bottom. They also feature a hook-and-loop strap or elastic bands around the ankle to help keep them in place and ensure adequate protection.

Item no.	Size	Characteristics	Packaging
NONE27143XL	One size fits most	Fluid-resistant	50/box, 150/case
NONE27144	One size fits most	Impervious	50/box, 150/case

NONE27143XL

- Fluid-resistant material above the ankle
- Fluid protection on the bottom
- Sealed seams on the foot



NONE27143XL

NONE27144

- Impervious PE-coated polypropylene material
- Fluid protection on the bottom
- Hook-and-loop strap for an adjustable fit



NONE27144


Innovative packaging to simplify mask selection

Medline's packaging is designed to make it easy for your staff to instantly identify many important mask features.

Colour-coded packaging

- Purple: Type IIR fluid-resistant masks
- Green: Type II surgical masks
- Orange: Type II procedural masks
- Blue: Respirators

Perforations and a convenient thumb notch
For easy opening and dispensing



Product classification
For proper selection and peace of mind

Face protection

Fluid-resistant facemask with shield⁷



NONE27410EL

Protection against splashes and fluids via this mask with shield and elastic earloops.

- Optically clear polyester shield
- Anti-fog foam strip

Item no.	Description	Packaging
NONE27410EL	Fluid-resistant ear loop facemask with shield	25/bx, 4/cs

Fluid-resistant facemask with shield and ties⁷



NONE27710

Protection against splashes and fluids via this mask with shield and tie attachments.

- Optically clear polyester shield
- Anti-fog foam strip

Item no.	Description	Packaging
NONE27710	Max X fluid-resistant facemask with shield and ties	25/bx, 4/cs
NONE27410	Fluid-resistant ear loop facemask with shield	25/bx, 4/cs

Fluid-resistant facemask with ear loops⁷



NONE27412EL

- Anti-fog foam strip
- Can be worn with additional splash protection, such as face shield or eye wear

Item no.	Description	Packaging
NONE27412EL	Fluid-resistant ear loop facemask	25/bx, 4/cs

Disposable replacement lenses



NONFRAME

Item no.	Description	Packaging
NONFRAME	Reusable frames, available in different colours	100/bx

Item no.	Description	Packaging
NONLENS	Disposable replacement lenses	250/bx

Face shield



NONEFS500

- Anti-fog foam strip
- Elastic headband with foam
- Not made with natural rubber latex

Item no.	Description	Packaging
NONEFS500	Full length disposable face shield with elastic foam strap	96/cs

Upgrade your level of safety in CSSD

Working in wet environments means being exposed to hazards, including slips, trips and falls on a daily basis. In the central sterile department, especially in decontamination areas where fluids accumulate, this represents a risk for staff.

Medline's floor safety solution addresses both the hazard of falls and the need for quick and safe fluid removal.

Liqui-Loc™ Floor Mat



- Easy disposal for hygiene and safety
- Impervious, anti-slippery back for standing or placing tools
- Drip-safe, locking fluid with absorbent material




Item no.	Description	Size	Absorbency	Anti-slip	Packaging
ORLLPBS	Liqui-Loc; small	71.1 x 101.6 cm	2.2 L	Yes	13/cs
ORLLPBM	Liqui-Loc; medium	71.1 x 142.2 cm	3.0 L	Yes	10/cs
ORLLPBL	Liqui-Loc; large	71.1 x 182.8 cm	3.9 L	Yes	8/cs

Double-Sided Absorbant Floor Pad



- Promote hygiene with a high-absorbency pad
- Speed up the cleaning process and absorb fluid spills on the floor
- Do not stand on



Item no.	Description	Size	Absorbency	Anti-slip	Packaging
ORDSFMAT	Double-sided absorbent floor pad	75 x 38 cm	7 L	No	50/cs

For more information, please contact your Medline account manager or visit our website: medline.eu/uk



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
medline.eu/uk
uk-customerservice@medline.com

FOLLOW US  

² Gemini wraps are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labeling.



³ Tray Liner sheets are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labeling.



⁴ GEMTRAY transport trays are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labeling.



⁵ Medline exam gloves are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labeling.



⁶ Shoe covers and boot covers are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labeling.



⁷ NONE27xxx are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labeling.

