

CareSet® Premium

Intravitreal injection (IVT) sets



**Making
intravitreal
injections faster
and safer**



All sets
below are
customisable
upon request

KER70238B

- 1 drape, 100 x 90 cm
- 1 pair of plastic Kocher forceps, 14 cm
- 10 non-woven swabs, 7.5 x 7.5 cm
- 3 eye sticks, 6.6 x 0.5 cm
- 2 cups, 60 ml (clear and red)
- 1 barraquer eye speculum, 4.5 cm
- 1 wooden, cotton-tipped applicator, 15 cm
- 3 adhesive strips, 1.3 x 10.2 cm
- 1 plastic scleral marker, double-ended, 3.5 – 4 mm
- 1 adhesive fenestrated drape, 50 x 60 cm, Ø 5 cm
- 1 peel pouch, 28 x 39 cm
- Pack factor: 24





KER70247B

- 1 drape, 100 x 120 cm
- 5 non-woven swabs, 7.5 x 7.5 cm
- 1 barraquer eye speculum, 4.5 cm
- 1 cup 120 ml, clear
- 3 wooden, cotton-tipped applicators, 15 cm
- 1 plastic scleral marker, double-ended, 3.5 – 4 mm
- 1 adhesive fenestrated drape, 75 x 90 cm, 6 x 8 cm
- 1 large, rigid blister
- Pack factor: 28




RGEYSP01B

- 1 adhesive fenestrated drape, 50 x 60 cm, Ø 5 cm
- 5 eye sticks, 6.6 x 0.5 cm
- 1 plastic scleral marker, double-ended, 3.5 – 4 mm
- 1 barraquer eye speculum, 4.5 cm
- 10 non-woven swabs, 7.5 x 7.5 cm
- 1 back table cover, 100 x 110 cm
- 1 peel pouch, 25 x 39 cm
- Pack factor: 40



KER70179B

- 1 drape, 100 x 90 cm
- 1 plastic scleral marker, double-ended, 3.5 – 4 mm
- 1 cup, 60 ml (clear)
- 10 non-woven swabs, 10 x 10 cm
- 1 barraquer eye speculum, 4.5 cm
- 1 adhesive fenestrated drape, 50 x 60 cm, Ø 5 cm
- 1 peel pouch, 25 x 35 cm
- Pack factor: 32



EYE001

- 1 drape, 75 x 75 cm
- 5 non-woven swabs, 7.5 x 7.5 cm
- 1 adhesive fenestrated drape, 50 x 60 cm, Ø 6 cm
- 1 barraquer eye speculum, 4.5 cm
- 1 peel pouch, 25 x 35 cm
- Pack factor: 40



EYE002A

- 1 drape, 100 x 90 cm
- 3 non-woven swabs, 7.5 x 7.5 cm
- 3 non-woven balls, medium-sized
- 1 cup, 60 ml (clear)
- 1 barraquer eye speculum, 4.5 cm
- 2 eye sticks, 6.6 x 0.5 cm
- 1 pair of plastic Kocher forceps, 14 cm
- 1 plastic scleral marker, double-ended, 3.5 – 4 mm
- 1 adhesive fenestrated drape, 40 x 75 cm, Ø 7 cm
- 1 peel pouch, 28 x 39 cm
- Pack factor: 32



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk
uk-customerservice@medline.com

FOLLOW US  

These sterile packs are procedure packs according to article 12 of the MDD 93/42/EEC. They are intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. ML1104_EN 06/2021.

Procedure pack producer Medline International France SAS
5 Rue Charles Lindbergh,
44110 Châteaubriant, France