

Compass Pressure Monitoring Kits

Measure compartment pressure precisely in just seconds.



**Precise pressure
made simple**



Thoracentesis Pressure Monitoring Kit

Reference: COMPTHORA1

1x	Drape; wrapping, 75 x 75 cm
1x	Compass Thoracentesis
1x	Syringe; 5 ml, L/S, 3 pieces
1x	Needle; 25 g x 2.5 cm
1x	Needle; 22 g x 3.2 cm
1x	Needle; 18 g, 4 cm, with filter
1x	Syringe; 10 ml, L/L, 3 pieces
1x	Syringe; 50 ml, L/L, 3 pieces
1x	Scalpel; n°11
5x	Swabs; non-woven, 7.5 x 7.5 cm
1x	Drape; adhesive, fenestrated, 45 x 75 cm, Ø 7 cm
1x	Dressing; adhesive, 5 x 7 cm
1x	Irrigation set Y
1x	Tray; blue, 25 x 13 x 2.5 cm

Pack factor: 10



Learn more:
Scan the QR code

Lumbar Pressure Monitoring Kit

Reference: COMPLUMB1

1x	Compass Lumbar
1x	Needle; spinal, 22 g, 9 cm
3x	Swabs; non-woven, 7.5 x 7.5 cm

Pack factor: 16



Learn more:
Scan the QR code

Single sterile kits

In addition to the specialised Compass device, each kit contains all the tools you'll need during pressure monitoring procedures.



Intra-Abdominal Pressure Monitoring Kit 1

Reference: COMPABD01

1x	Compass Universal Hg
3x	Syringe; 10 ml, L/L, 3 pieces
1x	Forceps; Kelly, metallic, 14 cm
1x	Tubing extension; 100 cm
5x	Swabs; non-woven, 7.5 x 7.5 cm

Pack factor: 16



Intra-Abdominal Pressure Monitoring Kit 2

Reference: COMPABD02

1x	Compass Universal Hg
1x	Connector LL to Foley
1x	Connector LL to drainage bag
1x	3-way stop cock
1x	Cap LL
1x	3-way stop cock, with extension 30 cm
1x	Tray 25 x 13 x 2.5 cm
1x	Drape 75 x 90 cm

Pack factor: 5



Intra-Abdominal Pressure Monitoring Kit 3

Reference: COMPABD03

1x	Compass Universal Hg
1x	Connector LL to Foley
1x	Connector LL to drainage bag
1x	3-way stop cock
1x	Needleless connector, neutral
1x	Tubing extension 100 cm
1x	Tray 25 x 13 x 2.5 cm
1x	Drape 75 x 90 cm

Pack factor: 10



Learn more about these kits:
Scan the QR code



Acute Compartment Pressure Monitoring Kit

Reference: COMPACUT1

2x	Swabs; non-woven 7.5 x 7.5 cm
1x	Syringe; 10 ml, L/L, 3 pieces
1x	Needle; 18 g, 5 cm
1x	Compass Universal Hg

Pack factor: 16



Learn more:
Scan the QR code



Chronic Compartment Pressure Monitoring Kit

Reference: COMPCHRON1

1x	Compass Universal Hg
6x	Needle; 18 g, 5 cm
1x	Syringe; 20 ml, L/L, 3 pieces
10x	Swabs; non-woven 7.5 x 7.5 cm

Pack factor: 16



Medline Industries Ltd.
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd.
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

These products are class IIa medical devices intended to be used by healthcare professionals.
Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. ML1204_EN 01/2024.



Medline International France SAS
5 Rue Charles Lindbergh,
44110 Châteaubriant, France

