



CENTURION
MEDICAL PRODUCTS

Compass Thoracentesis

Measure pleural pressure instantly


**Precise
pressure
made simple**



Why wonder? Avoid the risks and consequences of too little or too much drainage.

Critical data may be lost with traditional diagnostic tools.

In 2018, around 250,000 cases of pleural effusions were diagnosed in the National Health Service (NHS). Pleural manometry can provide critical information and help prevent complications from insufficient or excessive drainage volumes.



The meaningful target is pressure, not volume

However, many traditional diagnostic tools are time-consuming to set up and prone to inaccuracies. Without the appropriate tools for measurement, clinicians are:

- Left to extract an arbitrary volume for thoracentesis
- Risking the health of the patient
- Potentially performing repeat procedures
- Risking re-expansion pulmonary edema (RPE), which is associated with a 15 – 20% mortality rate

If you're not measuring pleural pressures, how do you know how much pleural effusion is to be drained?

Why just measure volume when your target is measured in pressure? Pleural manometry is evidence-backed as a useful tool during pleural effusion evacuation for diagnostic and therapeutic treatments. **Take the guesswork and compromise out of pleural effusion drainage.**

Compass simplifies pleural manometry.

When attached to your thoracentesis catheter, the sterile, compact, single-use Compass provides instant, precise pleural pressure readings on its digital display, and will continue to provide exact measurements during pleural effusion evacuation, as you reach optimal levels. Small, light, self-contained and cost-effective, Compass removes obstacles that often prevent routine measurements. This can enhance patient outcomes while controlling risks.

Good for you; great for your patients
Compass: Precise pressure made simple

Product features

- **Measures** positive and negative pressures in cmH_2O , eliminating unit conversions
- **Rapid** digital display provides indisputable, quantifiable readings in real time
- **Accurate** measurements
- **Secure**, universal luer lock connections
- **Convenient**, single-item, sterile packaging
- **Compact**, easy-to-use device



Compass Thoracentesis
CTHR002-INT

For more information on this product, please contact your Medline account manager or visit our website: uk.medline.eu



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
uk.medline.eu
uk-customerservice@medline.com


FOLLOW US  

Compass UniversalHG is a class IIa medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. ML953_EN 04/2021

 Centurion Medical Products Corporation
100 Centurion Way
Williamston, MI 48895, USA

 2797 (BSI, The Netherlands)