



Medline offers protective sheets and underpads in a range of sizes and with different levels of absorption.

Keep patients comfortable while protecting their beds.

Protection Plus disposable underpads

Surface protection with high absorbency and odour control

- Soft, non-woven top sheet helps ensure patient comfort
- Absorbent material for fluid containment
- Variety of sizes and weights to help with different budgets and needs
- Impervious polypropylene back sheet


Standard weight

Stock code	Size	Packaging
MSC281230C	43 x 61 cm	25/pk, 300/cs
MSC281241C	58 x 58 cm	25/pk, 200/cs
MSC281242C	58 x 91 cm	25/pk, 150/cs



Deluxe weight

Stock code	Size	Packaging
MSC281264C	58 x 91 cm	25/pk, 150/cs



Deluxe weight (with polymer)

Stock code	Size	Packaging
MSC282046	60 x 60 cm	25/pk, 100/cs

Protective sheets

Cost-effective surface protection

- Polypropylene backing for preventing leakage
- Available in various sizes and in 1- or 2-ply tissue to help you meet your protection needs

Stock code	Size	Packaging
MSCE261212	80 x 140 cm (32 x 55 in.), 1-ply	100/cs
MSCE261220	80 x 210 cm (32 x 83 in.), 1-ply	100/cs
MSCE261420	80 x 210 cm (32 x 83 in.), 2-ply	100/cs



Support for success

We understand that introducing new protocols or products can be a significant investment of time and effort. We can help. Our team of experts offers you clinical and product support, to help ensure successful implementation. They can provide you with:

- Product in-service training
- Specialised education on product usage and best practices
- A comprehensive skin management programme.

To find out more about our clinical and educational support services, contact your local Medline representative.



Medline Industries Ltd.

3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237

www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd.

1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

These products are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. ML560_EN 09/2023



Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France

