

Curad® 3G

Third-generation stretch vinyl powder-free exam gloves



Our proprietary formulation offers a combination of softness, strength and stretchability. The resin in the Curad 3G gloves has much longer polymer chains than other synthetic gloves. The longer polymer chains make the resin **more workable, flexible, and resilient.**

Benefits and features

- Designed to be more conforming to hands for excellent tactile sensitivity; made of material containing polyvinyl chloride (PVC) resin and plasticisers
- Latex-free and accelerator-free
- Acceptable quality level (AQL) for EN455-1: 1.5
- Force at break (FAB) median value: 5.1 newtons
- Personal protective equipment (PPE)-certified


Ordering information

Item number	Size	Packaging
CUR8233	X-small	100 ea/bx; 10 bx/cs
CUR8234	Small	100 ea/bx; 10 bx/cs
CUR8235	Medium	100 ea/bx; 10 bx/cs
CUR8236	Large	100 ea/bx; 10 bx/cs
CUR8237	X-large	90 ea/bx; 10 bx/cs




Sustainable packaging technology

MicroStop™ film



Flat-layered packaging



Colour-coded sizing



100% recyclable



Recommended for use in



Non-critical care unit



Household cleaning



Labour and delivery

Regulations and quality standards

- Class I non-sterile medical device, per Medical Device Regulation (EU) 2017/745 and the UK Medical Device Regulation
- Compliant with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended
- Compliant with standards EN ISO 21420; EN ISO 374-1; EN ISO 374-2; EN 16523-1; EN ISO 374-4 and EN ISO 374-5
- Tested for viral penetration, in accordance with EN ISO 374-5
- Compliant with Regulation (EU) No. 10/2011 for plastic materials and articles intended to come in contact with food, except 'fatty food'



Explore how Medline's gloves range, including Curad® 3G, makes a difference around the world. Contact your local Medline representative, or visit our website: www.medline.eu/uk

This product is a class I non-sterile medical device intended for use by healthcare professionals. Before use, consult the instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML106_EN 06/2024

Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France