



# DERMASSURE™ GREEN

A new alternative. Not made with chemical accelerators or natural rubber latex.

- Made from a next-generation polychloroprene material that provides a soft feel and enhanced flexibility and tactile sensitivity.
- Excellent choice as an underglove in combination with any of our lighter colored gloves.
- When used as an underglove, DermAssure™'s proprietary polymer coating on the outer surface of the glove bonds with the outer glove to feel like you're wearing just one glove for enhanced performance.
- Dark green color allows for rapid detection of perforations to the outer glove.
- Powder-free and coated with a synthetic polymer to assist with damp-hand donning.
- ARC (anti-roll down cuff) was designed to reduce cuff roll-down common with many surgical gloves.
- Improved former design to enhance comfort.

| Item Number | Size | Packaging         |
|-------------|------|-------------------|
| MSG6555I    | 5 ½  | 50 ea/bx, 4 bx/cs |
| MSG6560I    | 6    | 50 ea/bx, 4 bx/cs |
| MSG6565I    | 6 ½  | 50 ea/bx, 4 bx/cs |
| MSG6570I    | 7    | 50 ea/bx, 4 bx/cs |
| MSG6575I    | 7 ½  | 50 ea/bx, 4 bx/cs |
| MSG6580I    | 8    | 50 ea/bx, 4 bx/cs |
| MSG6585I    | 8 ½  | 50 ea/bx, 4 bx/cs |
| MSG6590I    | 9    | 50 ea/bx, 4 bx/cs |

## SPECIFICATIONS:

- Freedom from holes per EN455: AQL 0.65
- Synthetic polychloroprene (not made with natural rubber latex)
- Synthetic polymer coating
- Finger thickness: 0.20 mm
- Colour: dark green
- Grip: smooth
- Accelerator free

**Recommended for use in:**  
General surgery

### Regulations and quality standards

- Certified to EN ISO 13485
- Medical Device Directive 93/42/EEC
- Class IIa sterile Medical Device
- Certified to "Personal Protective Equipment" council directive 89/686/EEC Category III



*To learn more about our products, please contact your Medline representative or visit us at [www.medline.com/uk](http://www.medline.com/uk)*