

# DermAssure™ Green

## A new alternative



- These gloves are made from a next-generation polychloroprene material that provides a soft feel and enhanced flexibility and tactile sensitivity.
- DermAssure Green is an excellent choice as an underglove in combination with any of our lighter-coloured gloves.
- When used as an underglove, DermAssure's proprietary polymer coating, on the outer surface of each glove, bonds with the outer glove to feel like you're wearing just one glove, for enhanced performance.
- Dark-green colour allows for rapid detection of perforations in the outer glove.
- Powder-free, these gloves are coated with a synthetic polymer to assist with damp-hand donning.
- The anti-roll down cuff (ARC) design helps lower the rate of cuff roll-down, which is common with many surgical gloves.
- Through an improved former design, DermAssure gloves offer enhanced comfort.

## Specifications

- Freedom from holes, per EN455-1; Acceptance quality limit (AQL) 0.65
- Synthetic polymer coating
- Accelerator-free
- Synthetic polychloroprene (not made with natural rubber latex)
- Glove length (per EN 455-2, size 7.5): 301 mm
- Finger tip thickness: 0.19 mm
- Dark-green colour
- Smooth texture

## Ordering information

Item code	Size	Packaging
MSG6555	5 ½	50 pairs/box, 4 boxes/case
MSG6560	6	50 pairs/box, 4 boxes/case
MSG6565	6 ½	50 pairs/box, 4 boxes/case
MSG6570	7	50 pairs/box, 4 boxes/case
MSG6575	7 ½	50 pairs/box, 4 boxes/case
MSG6580	8	50 pairs/box, 4 boxes/case
MSG6585	8 ½	50 pairs/box, 4 boxes/case
MSG6590	9	50 pairs/box, 4 boxes/case

## Recommended for use in

- General surgery

## Regulations and quality standards

- Factory is EN ISO 13485-compliant
- Medical Device Directive 93/42/EEC
- Category III personal protective equipment (PPE); compliant with PPE Regulation (EU) 2016/425



This glove is a class IIa sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML423\_EN 07/2023

Learn more. Contact your Medline account manager or visit: [www.medline.eu/uk](http://www.medline.eu/uk)