



**DISPOSABLE ALL-AROUND
PROTECTION.**

Protective Apparel
For Nursing and Research

Why are more healthcare professionals choosing Medline for their all-around protection?

Our products provide high protection and great comfort.

We know you rely on disposable apparel to keep your staff protected and comfortable. For this reason, we use the finest materials available for protection, offering the greatest level of comfort. We apply the highest quality standards to our products and in our manufacturing facilities.



Medline International France is EN ISO 13485: certified. Medline has been manufacturing protective apparel for over 40 years.

Table of Contents:

Guide to Protective Materials	3
Levels of Barrier Protection	3
Isolation Gowns overview	4-5
Cover Gowns	6
Fluid Resistant Gowns	7
Impervious Gowns	8-9
Prevention Plus Procedure Gowns	10
Sterile Isolation Gowns	11
Foot Covers	12
Patient and Visitor Kits	13
Patient gowns	14
Warm up jackets	15
Scrub Suits	16
Clear Air Suits	18
Headwear	20-23



GUIDE TO PROTECTIVE MATERIALS

To help chose the right product, here is a brief description of the features and benefits of the main materials available at Medline. Refer back to this page as you look through the brochure.

Spunbond Polypropylene

The standard fabric for basic protective apparel; formed by bonding fibers together to form a single layer of breathable, woven-like material. Its main advantages are economy and comfort



Multi-Ply Material (SMS)

This fabric consists of inner layers of meltblown polypropylene sandwiched between outer layers of spunbond polypropylene. Strong, fluid-resistant and breathable, this multi-layer material provides an optimal mix of protection and comfort.




Coated Polypropylene

This fabric is created by coating soft, spunbond polypropylene with a layer of impervious polyethylene. This material is ideal for any situation demanding a high level of fluid protection. Available in medium weight and heavy weight materials.



Polyethylene

This single layer of polyethylene film provides high barrier properties at a very attractive price. This film is ideal for situations when impervious properties are needed.



Prevention Plus

This impervious and breathable material consists of a hybrid film laminated between two layers of bi-component spunbond material. This fabric is ideal for long interventions when impervious protection is needed.



Levels of Barrier Protection

European Directives and Norms to classify Protective apparel:


The products listed in the brochure are classified by the European regulation as PPE Category I - Minimum Protection according to the Directive 89/686/EEC or MDD class I non sterile according to the Directive 93/42/EEC. Please consult the Technical Data Sheet of the product in order to find how the product is classified.

Directive 93/42/EEC applies to medical devices and their accessories. Medical devices means any instrument, apparatus, appliance, software, material or other article, whether used alone or in combination, including the software intended by its manufacturer to be used specifically for diagnostic and/or therapeutic purposes and necessary for its proper application.

Directive 89/686/EEC commonly referred to as PPE (Personal Protective Equipment) details the conditions governing the placement and free movement within the European Community Market for the PPE goods and the safety requirements which PPE must satisfy in order to ensure the health protection and safety of users.

ISOLATION GOWNS OVERVIEW

Product family	Item number	Size	Gown length (cm)	Gowns Colour	Material	Wrists type	Neck closure	Resistance to liquid penetration EN20811 cmH ₂ O	Resistance to wet bacterial penetration EN ISO 22610 per EN 14126	Resistance to penetration by contaminated solid particles EN ISO 22612 per EN 14126	Intended Use			
Cover gowns	CRIE4000	L	114		Polypropylene (SPP)	Elastic	Ties	N/a	N/a	N/a	Lightweight gown designed to limit the spread of particulates from visitors. Intended to be used as a visitor gown. Not suitable for fluid or bacterial protection.			
	CRIE4000G	L												
	CRIE4000B	L												
	CRIE4001	XL	127			Knit cuffs								
	CRIE4001G	XL												
	CRIE4001B	XL												
	CRIE4001CF	XL												
	CRIE4001BCF	XL												
Fluid resistant isolation gowns	NONE27SMS3	L	114		Polypropylene (SMS)	Elastic	Ties	20	N/a	N/a	Recommended for procedures with low exposure to non biological fluids.			
	NONE27SMS3XL	XL	127				Ties	20						
	NONELV200	L	114		Polypropylene (SMS)		Tape tab	28.4						
	NONELV200XL	XL	127				Tape tab	28.4						
	NONELV200C	L	114		Polypropylene (SMS)	Knit cuffs	Ties	28.4	Cass 1	Class 3	Recommended for procedures with medium exposure to non biological fluids. Isolation gowns are comfortable, breathable and generously sized for full coverage.			
	NONELV200CXL	XL	127											
	NONE386CF*	L	118		Polypropylene (SMS)	Knit cuffs	Ties	57.3						
	NONE386XLCF*	XL	135											
	NONE326CF	L	118		Polypropylene (SMS)		Velcro							
NONE326XLCF	XL	135												
Impervious isolation gowns,	CRIE4040C*	L	114		Polypropylene/ Polyethylene Bilaminate	Knit cuffs	Ties	124				Class 6	Class 3	Recommended for procedures with high exposure to fluids.
	CRIE4041C*	XL	135					124						
	CRIE4020CF*	L	114					78						
	CRIE4021VCF*	L	114				Velcro	> 149						
	CRIE4022VCF*	XXL	140											
	CRIE4010CF	L	114				Ties	> 120						
	CRIE4011CF	XL	127					> 120						
	CRIE4010RCF	L	114					91						
	NONE27116	L	114					111						
	NONE27116XL	XL	127					111						
	NONE27279OHXL	XL	127					111						
	CRIE5000	L	114		Polyethylene	Thumb loops	Over the head	50.5						
	CRIE5001	XL	134											
	Isolation gowns for chemotherapeutic procedures	NONE27457	L	114						Prevention Plus	Knit cuffs			
NONE27457XL		XL	127		Velcro	> 200								
Isolation gowns sterile	CRIE4020CE*	L	114		Polypropylene/ Polyethylene Bilaminate	Knit cuffs	Ties	78.5	Class 6	Class 3	Recommended for procedures with medium exposure to non biological fluids. Isolation gowns are comfortable, breathable and generously sized for full coverage.			
	NONE326XLCE	XL	135		Polypropylene (SMS)		Velcro	57.3	Class 1	Class 3	Recommended for procedures with high exposure to fluids.			

 * The isolation gown has the back panel made of breathable material

COVER GOWNS

Cover Gowns – Light Weight Polypropylene
Suitable for health care workers and visitors, helps to prevent contamination with particulates.

- Full coverage, breathable and comfortable
- Waist and neck ties, Elastic or Knit Cuffs

Item No.	Size	Color	Pkg
Cover gowns with elastic wrist			
CRIE4000	L	Yellow	10/bag, 50/case
CRIE4000G	L	Green	10/bag, 50/case
CRIE4000B	L	Blue	10/bag, 50/case
CRIE4001	XL	Yellow	10/bag, 50/case
CRIE4001G	XL	Green	10/bag, 50/case
CRIE4001B	XL	Blue	10/bag, 50/case
Cover gowns with knit cuffs			
CRIE4001CF	XL	Yellow	10/bag, 50/case
CRIE4001BCF	XL	Blue	10/bag, 50/case



CRIE4001



CRIE4001G



CRIE4001BCF

FLUID RESISTANT GOWNS

Fluid Resistant Isolation Gowns - Light Weight SMS
Suitable for procedures with low risks of exposure to non-biological fluids.

- Breathable and comfortable
- Waist and neck ties, elastic wrists or cuff

Item No.	Size	Color	Pkg
Blue, neck and waist ties, elastic wrists, flat pack			
NONE27SMS3	L	Blue	10/bag, 50/case
NONE27SMS3XL	XL	Blue	10/bag, 50/case
Elastic wrists, tape tab neck			
NONELV200	L	Yellow	10/bag, 100/case
NONELV200XL	XL	Yellow	10/bag, 100/case
Knit cuffs, neck and waist ties			
NONELV200C	L	Yellow	10/bag, 100/case
NONELV200CXL	XL	Yellow	10/bag, 100/case




NONE27SMS3

Fluid Resistant, Isolation Gown - Medium Weight SMS
Suitable for procedures with moderate risks of exposure to non-biological fluids.

Gowns offer increased protection from non-biological fluids without sacrificing the level of comfort and breathability.

Item No.	Size	Color	Pkg
Knit cuffs, ties at the neck and waist, breathable back			
NONE386CF	L	Yellow	10/bag, 100/case
NONE386XLCF	XL	Yellow	10/bag, 100/case
Velcro neck closure, ties at the waist, knit cuffs			
NONE326CF	L	Yellow	10/bag, 100/case
NONE326XLCF	XL	Yellow	10/bag, 100/case

! NONE326XLCF is available in sterile- NONE326XLC



NONE386CF and NONE386XLCF has SMS front and sleeves and Spunbond Polypropylene back



NONE326CF

IMPERVIOUS ISOLATION GOWNS

Impervious Isolation Gowns – Coated Polypropylene

The Single Use Bilaminate Isolation Gowns are impervious and suitable for procedures with high risks of exposure to fluids.

The isolation gowns are equipped with knit cuffs, ties at the neck and waist. NONE2727900HXL is designed with an over the head style as well as open back. The design ensures that isolation gowns remain securely in place while in use.

- Class 6 for Resistance to wet bacterial penetration (EN ISO 22610 per EN14126)
- Class 3 for Resistance to penetration by contaminated solid particles (EN ISO 22612 per EN 14126)



CRIE4011CF



CRIE4010RCF

Item No.	Size	Color	Pkg
CRIE4010CF	L	Yellow	10/bag, 50/case
CRIE4011CF	XL	Yellow	10/bag, 50/case
CRIE4010RCF	L	Red	10/bag, 50/case
NONE27116	L	Blue	10/bag, 50/case
NONE27116XL	XL	Blue	10/bag, 50/case
NONE272790HXL	XL	Blue	10/bag, 50/case

Impervious Isolation Gowns with breathable back - Coated Polypropylene

Suitable for procedures with high risks of exposure to fluids. The material is tested against viral penetration according to standard ASTM F1671 or ISO 16604

The Single Use Isolation Gown's front and sleeves are made from impervious bilaminate material and the back of breathable material. Due to this design it offers a great balance between comfort and protection. The Isolation Gown are equipped with knit cuffs,ties or hook and loop at the neck, ties at the waist which ensure that gown remain securely in place.

- Class 6 for Resistance to wet bacterial penetration (EN ISO 22610 per EN 14126)
- Class 3 for Resistance to penetration by contaminated solid particles (EN ISO 22612 per EN 14126)

Item No.	Size	Color	Pkg	Material Viral penetration performance
CRIE4040C	L	Yellow	10/bag; 100/cs	Tested per ISO 16604 Class 6- EN 14126
CRIE4041C	XL	Yellow	10/bag; 100/cs	Tested per ISO 16604 Class 6- EN 14126
CRIE4020CF	L	Yellow	10/bag, 50/case	Tested per ISO 16604 Class 6- EN 14126
CRIE4021VCF	L	Yellow	10/bag, 50/case	Material passes the ASTM F1671 test. No penetration of bacteriophages
CRIE4022VCF	XL	Yellow	10/bag, 50/case	

⚠ CRIE4020CF is available in sterile version- CRIE4020CE



CRIE4040C, CRIE4041C, CRIE4020CF, CRIE4021VCF, CRIE4022VCF have impervious front and sleeves, breathable back



CRIE4022VCF

Impervious Isolation Gown – Polyethylene

Recommended against fluid exposure. Suitable for Infective wards, ICU, Emergency.

Impervious Polyethylene isolation gowns with over the head design and thumb loop style wrists. Individually folded and wrapped gowns in dispenser box.

- Class 6 for Resistance to wet bacterial penetration (EN ISO 22610 per EN 14126)
- Class 3 for Resistance to penetration by contaminated solid particles (EN ISO 22612 per EN 14126)
- Resistance to liquid penetration up to 50.5 cmH₂O

Item No.	Size	Color	Pkg
CRIE5000	L	Blue	1/bag, 15/box, 75/case
CRIE5001	XL	Blue	1/bag, 15/box, 75/case



CRIE5001

Poly Disposable White Apron

- Effective protection against liquids and worn over non-fluid resistant gowns or scrub suits to add additional protection at a convenient price.
- Made from polyethylene material

Item No.	Description	Pkg
NONE24274W	Medium weight, 71.1x116,8 cm individually wrapped, dispenser box	50/bx, 250/ cs



NONE24274W

PREVENTION PLUS PROCEDURAL GOWNS

High performance 3 layer material tested to resist against chemotherapy drugs*.

Isolation gown is equipped with hook-and-loop neck and knit cuffs at the wrist and tiest at the waist.

- Class 6 for Resistance to wet bacterial penetration (EN ISO 22610 per EN 14126)
- Class 3 for Resistance to penetration by contaminated solid particles (EN ISO 22612 per EN 14126)
- Resistance to liquid penetration >200cmH₂O

Chemotherapy Drug Tested*

The Prevention Plus Material has been tested against the following chemotherapy drugs per ASTM F739-99a:

- Etoposide (Toposar)

· Paclitaxel (Taxol)

· Thio-Tepa

· Cisplatin
- Mitomycin C

· Vincristine Sulfate

· Dacarbazine

· 5-Fluorouracil
- Carmustine (BICNU)

· Cyclophosphamide (Cytoxan)

· Doxorubicin Hydrochloride

(Test data available upon request)

Item No.	Size	Color	Pkg
NONE27457	L	Blue	10/bag, 30/case
NONE27457XL	XL	Blue	10/bag, 30/case



NONE27457- Gown
NONE27348PXL- Boot covers

Complete protection in combination with the boot covers **NONE27348PXL**

STERILE ISOLATION GOWNS

Medium Weight SMS Sterile Isolation Gown

Suitable for procedures with moderate risks of exposure to non-biological fluids and where sterile protection is required.

The gown is manufactured from SMS material. It offers increased protection from contaminated non-biological fluids and microbial penetration without sacrificing the level of comfort and breathability.

Item No.	Size	Color	Pkg
NONE326XLCE	XL	Yellow	50/case

Coated Polypropylene Sterile Isolation Gown

Suitable for patient care with high exposure to fluids that requires additional sterile protection.

The gown is impervious in the critical areas - sleeves and front. It has a breathable back manufactured from spunbond polypropylene. The design provides protection and comfort to the user.

- Class 6 for Determination of resistance of protective clothing materials to penetration by blood-borne pathogens (EN ISO 16604 per EN 14126)
- Class 6 for Resistance to wet bacterial penetration (EN ISO 22610 per EN14126)
- Class 3 for Resistance to penetration by contaminated solid particles (EN ISO 22612 per EN 14126)

Item No.	Size	Color	Pkg
CRIE4020CE	L	Yellow	50/case

SHOE COVERS

Medline shoe covers are manufactured from polyethylene or polypropylene. Choose the best option based on comfort, protection and price. The shoe covers have specific patterns with non-skid properties.

Item No.	Size	Material	Characteristics	Pkg
CRIE2002	One Size Fits Most	Polypropylene non woven	Protection	100/bx, 300/cs
CRIE2010	One Size Fits Most	Polyethylene film	Impervious	1000cs
CRIE2102W	One Size Fits Most	Polypropylene non woven	Protection	500/bx, 1000/cs
CRIE2020	One Size Fits Most	Polyethylene film and Polypropylene non woven	Impervious	200/cs
NONE28858	L	Polypropylene non woven	Fluid resistant	100/bx, 300/cs
NONE28859	XL	Polypropylene non woven	Fluid resistant	100/bx, 200/cs



CRIE2002



CRIE2010



CRIE2102W



CRIE2020



NONE28858

BOOT COVERS

Medline knee-high boot covers have a convenient one size fits all design. The range offers high breathability, fluid protection and comfort. To ensure better securement the boot cover have a hook and loop strap or elastic bands around the ankle.

Item No.	Size	Characteristics	Pkg
NONE27143XL	One Size Fits Most	Fluid resistant	50/bx, 150/cs
NONE27144	One Size Fits Most	Impervious	50/bx, 150/cs
NONE27348PXL	One Size Fits Most	Breathable and impervious	50/bx, 150/cs

NONE27143XL

- Breathable, fluid-resistant material above the ankle
- High fluid protection on the bottom
- Sealed seams on foot



NONE27143XL

NONE27144

- Impervious PE coated polypropylene material
- High fluid protection on the bottom
- Hook-and-loop strap for adjustable fit



NONE27144

NONE27348PXL

- Made from blue Prevention Plus material
- Hook-and-loop closure
- Foam strip, non-skid bottom



NONE27348PXL

PATIENT AND VISITOR KITS

The kit is designed as a convenient solution that minimizes patient insecurities that may appear during daily procedures. It provides all patient needs, in one easy to use bag.

Item No.	Each PE bag contains	Pkg
Patient kit NONE27046KIT	- 1 patient gown, short sleeves, waist ties, tape tab on neck, blue, L, 114x139 cm - 1 underwear, dark blue - 1 clip cap, white, 53 cm - 2 shoe covers CRIE2102W, white, size L	50 PE bags/cs
Patient kit for child (4 to 9 years) NONE27KITK	- 1 patient gown for child, short sleeves, neck and waist ties, blue, 70x105 cm - 1 underwear for child, dark blue - 1 clip cap for child, white, 40 cm - 2 shoe covers for child, white, size L, 32 cm	50 PE bags/cs
Patient kit NONE35LKIT	- 1 patient gown, short sleeves, waist ties, neck ties, blue, L size - pants, blue, SMS, 35 gsm, L size - 1 clip cap, white, 53 cm - 2 shoe covers CRIE2020	50 PE bags/cs
Patient kit NONE352XLKIT	- 1 patient gown, short sleeves, waist ties, neck ties, blue, 2XL size - pants, blue, SMS, 35 gsm, 2XL size - 1 clip cap, white, 53 cm - 2 shoe covers CRIE2020	50 PE bags/cs
Patient kit KIT1030	- 1 patient gown, short sleeves, neck and waist ties, blue, L, 114x139 cm - 1 kimono gown, long sleeves, waist ties, 120x185 cm, blue - 1 clip cap, blue, 53 cm - 1 underwear, dark blue - 1 pair of non-slip socks, red, XL	40 PE bags/cs

Item No.	Each PE bag contains	Pkg
Visitor kit CRIE4000BKIT	- 1 cover gown, blue, waist and neck ties, elastic L size, - 1 face mask with earloops, type II - 1 clip cap, white, 53 cm - 2 shoe covers, white, size L	50 PE bags/cs

PATIENT GOWNS

The products are manufactured from soft SMS with high opacity to provide the right intimacy and comfort for the patient during the procedures.

Item No.	Description	Size	Pkg
NONE27046	Sleeveless patient gown, waist ties, blue	L	50/cs
NONE27046XL	Sleeveless patient gown, waist ties, blue	XL	50/cs
NONE27246SLV	Patient shirt, waist ties, tape tab, blue	L	50/cs
NONE27046SL	Short sleeved patient gown, neck and waist ties, blue	L	50/cs
NONE27SLK	Short sleeved patient gown for child, neck and waist ties, blue	suitable for child between 4 to 9 years	100/cs
CRIE4060K	Kimono gown, long sleeves, waist ties, 120x185 cm, blue	L	60/cs



DISPOSABLE PATIENT CAPE

- Disposable patient cape is manufactured from SMS
- Material provides warmth for extended wear.
- Tape tab feature provides a secure fit and added comfort for patient.

Item No.	Colour	Material	Pkg
NONE27036	Blue	SMS	100/cs



COLLAR COVER

For additional warmth and protection around the neck.

Item No.	Colour	Size	Material	Pkg
SPTE-UG1317	Blue	Standard	Spunlace 70gsm	50/bg, 100/cs



SPTE-UG1317

CLASSIC WARM UP JACKETS

A comfortable range to wear in the OR and cooler open spaces in the hospital.

- Made of comfortable SMS material
- Five snap closure
- Reinforced pockets
- Gentle knit collar and knit cuffs

Item No.	Description	Sleeves	Pkg
NONERP600S	Blue S	Long	10/bg, 30/cs
NONERP600M	Blue M	Long	10/bg, 30/cs
NONERP600L	Blue L	Long	10/bg, 30/cs
NONERP600XL	Blue XL	Long	10/bg, 30/cs
NONERP700M	Blue, M	Short – no cuffs	10/bg, 30/cs
NONERP700L	Blue, L	Short – no cuffs	10/bg, 30/cs
NONERP700XL	Blue, XL	Short – no cuffs	10/bg, 30/cs
NONERP800S	Blue, S	3/4	10/bg, 30/cs
NONERP800M	Blue, M	3/4	10/bg, 30/cs
NONERP800L	Blue, L	3/4	10/bg, 30/cs
NONERP800XL	Blue, XL	3/4	10/bg, 30/cs



NONERP600L



NONERP700L



NONERP800L

MEDLINE SINGLE USE SCRUB SUITS

Medline single use scrub suits are manufactured from breathable SMS material. The products are intended to be used by healthcare professionals in hospitals and clinics to maintain a high level of hygiene and contribute to infection control. Medline portfolio has two ranges of scrub suits the Advanced range and the Soft range. Both ranges are designed in collaboration with healthcare professionals.

Advanced Scrub suits



- Manufactured from 42 gsm SMS material with an optimal opacity and softness
- Low linting, with an antistatic treatment and non-biological fluid resistance
- Available in 7 sizes from XS to 3XL
- Available in: **blue**, **green** and **purple**
- Blue and green range is available with shirts and pants packed separately or in sets
- Purple range is available with shirts and pants packed separately
- Latex free

Soft Scrub suits

- Manufactured from 35 gsm SMS material with an optimal opacity and softness
- Low linting, with an antistatic treatment and fluid resistance
- Available in 7 sizes from XS to 3XL
- Available in **blue** and **green** colour
- Blue range is available with shirt and pants packed separately or in sets
- Green range is available with shirts and pants packed in sets
- Latex free



Soft Scrub Suit Set

Size	Blue scrub suits	Units per case	Green scrub suits	Units per case
XS	P35PBXS	25	P35PGXS	25
S	P35PBS	25	P35PGS	25
M	P35PBM	25	P35PGM	25
L	P35PBL	25	P35PGL	25
XL	P35PBXL	25	P35PGXL	25
2XL	P35PB2XL	25	P35PG2XL	25
3XL	P35PB3XL	25	P35PG3XL	25

Soft Scrub Suits Packed Separately

Size	Blue shirts	Units per case	Blue pants	Units per case
XS	PXS9110S	50	PXS9110P	50
S	PS9110S	50	PS9110P	50
M	PM9110S	50	PM9110P	50
L	PL9110S	50	PL9110P	50
XL	PXL9110S	50	PXL9110P	50
2XL	P2XL9110S	50	P2XL9110P	50

Advanced Scrub Suits Set

Size	Blue scrub suits	Units per case	Green scrub suits	Units per case
XS	P42PBXS	25	P42PGXS	25
S	P42PBS	25	P42PGS	25
M	P42PBM	25	P42PGM	25
L	P42PBL	25	P42PGL	25
XL	P42PBXL	25	P42PGXL	25
2XL	P42PB2XL	25	P42PG2XL	25
3XL	P42PB3XL	25	P42PG3XL	25

Advanced Scrub Suits Packed Separately

Size	Blue shirts	Units per case	Blue pants	Units per case	Size	Green shirts	Units per case	Green pants	Units per case
XS	P42PBXS-S	50	P42PBXS-P	50	XS	P42PGXS-S	50	P42PGXS-P	50
S	P42PBS-S	50	P42PBS-P	50	S	P42PGS-S	50	P42PGS-P	50
M	P42PBM-S	50	P42PBM-P	50	M	P42PGM-S	50	P42PGM-P	50
L	P42PBL-S	50	P42PBL-P	50	L	P42PGL-S	50	P42PGL-P	50
XL	P42PBXL-S	50	P42PBXL-P	50	XL	P42PGXL-S	50	P42PGXL-P	50
2XL	P42PB2XL-S	50	P42PB2XL-P	50	2XL	P42PG2XL-S	50	P42PG2XL-P	50
3XL	P42PB3XL-S	50	P42PB3XL-P	50	3XL	P42PG3XL-S	50	P42PG3XL-P	50

Size	Purple shirts	Units per case	Purple pants	Units per case
XS	P42PPXS-S	50	P42PPXS-P	50
S	P42PPS-S	50	P42PPS-P	50
M	P42PPM-S	50	P42PPM-P	50
L	P42PPL-S	50	P42PPL-P	50
XL	P42PPXL-S	50	P42PPXL-P	50
2XL	P42PP2XL-S	50	P42PP2XL-P	50
3XL	P42PP3XL-S	50	P42PP3XL-P	50

MEDLINE CLEAN AIR SUITS

Medline single use Clean Air Suit is a unique, innovative product that provides high comfort for healthcare professionals and increased safety for patients. Medline Clean Air Suits are specifically designed for use in hospitals and clinics to maintain a high level of hygiene and contribute to infection control. For more details about clean air suits please check the clean air suit brochure.



Designed with Attention to Details

Medline Clean Air Suits have been designed in collaboration with healthcare professionals through multiple trials, adjustments and validations. This high attention to detail provides all features and benefits required by the intensive activities of the operating room staff.



Knit cuffs at the sleeves



Convenient chest pocket with a badge holder



Round neck



Waist cuffs with 2 pockets on the sides



3 buttons at the neck



2 big pockets on the pants




Kitted cuffs at ankles



Adjustable waist ties



Shirt and pants are in a PE bag



XS to 3XL

Available in 7 sizes

Contributes to Reducing Air Contamination in the Ward and Surgical Theatre

Medline Clean Air Suits comply with the EN 13795 standard. All the tests performed demonstrate how the Clean Air Suit contributes to reduce airborne particles and preserve a clean air and environment in the ward and surgical theatre.*

	EN 13795 Characteristic	Test Method	EN 13795 Clean air suits requirements	Test results**
Material test Performance*	Resistance to microbial penetration - Dry (CFU)	EN ISO 22612	≤300	0
	Cleanliness - Microbial (CFU/100cm²)	EN ISO 11737-1	≤300	<9
	Cleanliness - Particulate matter (IPM)	EN ISO 9073-10	≤3,5	1,9
	Linting (Log10 (lint count))	EN ISO 9073-10	≤4,0	1,9

** Technical data are based on average.

Clean Air Suit Shirts

Size	Item Number	Units per case
XS	PCLBXS-S	50
S	PCLBS-S	50
M	PCLBM-S	50
L	PCLBL-S	50
XL	PCLBXL-S	50
2XL	PCLB2XL-S	50
3XL	PCLB3XL-S	50

Clean Air Suit Pants


Size	Item Number	Units per case
XS	PCLBXS-P	50
S	PCLBS-P	50
M	PCLBM-P	50
L	PCLBL-P	50
XL	PCLBXL-P	50
2XL	PCLB2XL-P	50
3XL	PCLB3XL-P	50

MEDLINE HEADWEAR

Comfort Line

The Comfort Line is Medline’s preferred headwear product range. The line offers a high level of comfort and comes in a variety of shapes and designs that are in accordance with the aesthetics and protection of the hospital. The Comfort Line is created from soft and breathable materials.

The most preferred Medline range for the OR!




Surgical Hood

FS61100

100/bx
600/cs

Surgical hood • Green




Surgical Hood

FS61150A

100/bx
600/cs

Surgical hood • Blue




Surgical Hood

FS61250

50/bx
300/cs

Surgical hood with sweat band • Blue



Bouffant Cap

FS61300

100/bx
600/cs

Bouffant cap with a latex free elastic band • Green • Size 53 cm




Bouffant Cap

FS61350

100/bx
600/cs

Bouffant cap with a latex free elastic band • Blue • Size 53 cm




Surgical Cap

FS61450

100/bx
600/cs

Surgical cap • sailor style • Blue




Surgical Cap

FS61500

100/bx
600/cs

Surgical hood • Elasticated back • Contoured front • Green




Surgical Cap

FS61550A

100/bx
600/cs

Surgical hood • Elasticated back • Contoured front • Blue




Surgical Cap

FS61600

75/bx
450/cs

Surgical cap for long hair • Elasticated • Green




Surgical Cap

FS61650

75/bx
450/cs

Surgical cap for long hair • Elasticated • Blue




Surgical Cap

FS61750A

100/bx
600/cs

Surgical cap with ties • Long hair • Blue




Surgical Cap

FS61850A

100/bx
600/cs

Surgical cap with ties • Short hair • Blue



Surgical Cap


FS61950

100/bx
600/cs

Surgical cap • Elasticated back • Flat front • Blue

Colorful line

The Colorful Line is made from spunlace viscose; a soft and breathable material. The range comes in a variety of colors and has a good comfort.




Surgical Cap

FS60950B

100/bx
600/cs

Surgical cap • Elasticated back • Flat front • Blue




Surgical Cap

FS60950G

100/bx
600/cs


Surgical cap • lasticated back • Flat front • Green



Surgical Cap

FS60950P100/bx600/cs

Surgical cap • Elasticated back • Flat front • Purple




Surgical Cap

FS60950W100/bx600/cs

Surgical cap • Elasticated back • Flat front • White

Basic line


The Basic Line was designed to satisfy the basic needs of the ward at a convenient price. The line is created from spunbond material with a basis weight from 12 gsm to 20 gsm and comes in blue, green, white or printed patterns.



Surgical Hood

FS63150100/bx600/cs


Light surgical hood • Width 17cm • Blue



Bouffant Cap

FS63360100/bx600/cs


Bouffant cap • size 53 cm • White



Bouffant Cap

FS64360200/bx2'000/cs


Bouffant cap • size 53 cm • White



Surgical Cap

FS60950PK100/bx600/cs


Surgical cap • Elasticated back • Flat front • Pink



Bouffant Cap

FS63300100/bx600/cs

Bouffant cap • Size 53 cm • Green



Bouffant Cap

FS63350100/bx600/cs


Bouffant cap • size 53 cm • Blue



Bouffant Cap

CRIE1000100/bx500/cs


Bouffant cap • size 53 cm • Green



Bouffant Cap

CRIE1003100/bx500/cs


Bouffant cap • size 60 cm • White



Bouffant Cap

CRIE1005100/bx500/cs


Bouffant cap • size 60 cm • Green



Bouffant Cap

NONE28234P100/bx500/cs

Bouffant cap with Multi-color print • size 60 cm • SMS



Bonnet

FS67050100/bx600/cs

Spunlace Bonnet • One size fit most • Blue



Bouffant Cap

CRIE1002100/bx500/cs


Bouffant cap • size 53 cm • White



Bouffant Cap

CRIE1004100/bx500/cs

Bouffant cap • size 60 cm • Blue



Bouffant Cap

NONE28233F100/bx500/cs


Bouffant cap with butterfly print • size 60 cm • SMS • Blue



Bouffant Cap

FS653601'000/cs

Clip cap • One size fits most • White



Surgical Cap

NONE28625100/bx500/cs

Surgical cap with ties • Top SMS; Side-absorbent scrim reinforced material • One size fits most •Light Blue

Medline also propose a large range of surgical, examination and respiratory masks.



MEDLINE

**SAFETY AND COMFORT
WITHOUT COMPROMISE.**

FACE MASKS
Surgical, Procedure & Respirator



**ALWAYS
ON**

Medline Industries Ltd
No 5 Booths Park
Chelford Road
Knutsford, Cheshire WA16 8GS
England
Tel: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk
uk-customerservice@medline.com

We reserve the right to correct any errors that may occur within this brochure.

©2017 Medline Industries, Inc. Medline, Boundary and Thumbs Up are registered trademarks of Medline Industries, Inc.

ML179-EN05/FL - 10/2017.