




**REDUCE CROSS-CONTAMINATION  
AND PROTECT YOUR STAFF.**

**Protective Apparel**  
For nursing and research

# Why are more healthcare professionals choosing Medline for their **all-around protection needs**?



**High-quality products**



**Comprehensive product range**



**Reliable partner**

Your staff are exposed to bacteria and bodily fluids on a daily basis, putting them at risk of cross-contamination. We know you rely on disposable apparel to keep them protected, but also comfortable. Therefore, we apply the highest quality standards to our products and in our manufacturing facilities to ensure that healthcare professionals get the protection they deserve.

Medline International France is EN ISO 13485 and ISO 9001 certified and complies with regulatory standards on the market. In addition, we are conscious about social responsibility and, as a result, some of our products are WRAP-certified. WRAP (Worldwide Responsible Accredited Production) is an internationally recognised social compliance certification body that verifies whether factories respect workers' labour and human rights.

## Table of Contents

European Regulations and Directives .....	3
Guide to Protective Materials .....	3
Isolation Gowns Overview .....	4-5
Cover Gowns .....	6
Examination Gowns .....	7
Impervious Cover Gowns .....	8-9
Impervious Isolation Gowns .....	9
Chemo Protection Isolation Gowns .....	10
Footwear .....	11
Disposable Scrub Suits .....	12-14
Warm-Up Jackets .....	15
Headwear .....	16-21
Patient Gowns .....	22
Patient Apparel .....	22
Patient Kits .....	23

## EUROPEAN REGULATIONS AND DIRECTIVES TO **CLASSIFY PROTECTIVE APPAREL**

Regulation (EU) 2016/425 on personal protective equipment (PPE)

The personal protective equipment included in Medline's range provides protection against minimal risks (contact with weak-action cleaning materials or prolonged contact with water) or more serious risks, such as exposure to substances and mixtures that are hazardous to health. Products listed in the brochure that are certified as PPE are compliant with PPE Category I or Category III requirements.

Regulation (EU) 2017/745 on medical devices (MDR)

Medical devices include any instrument, apparatus, appliance, software, implant, reagent, material or other article intended by the manufacturer to be used, alone or in combination, for human beings for a specific medical purpose. The medical devices included in Medline's range are intended to protect the patient from potential contamination. They are Class I devices.

Directive 2001/95/EC on general product safety

Products included in Medline's protective apparel range that do not fall under the definition of PPE or MDR are certified under the general product safety directive. These are general commodity products (e.g. underwear).

Please consult the product's technical data sheet to see its classification.

## GUIDE TO PROTECTIVE MATERIALS

To help you chose the right product, here are brief descriptions of the features and benefits of the main materials available from Medline. You can refer back to this page as you look through the brochure.

### Spunbond polypropylene (SPP)

This standard fabric is made by bonding fibres together to form a single layer of breathable woven-like material. It is suitable for protecting users against mechanical particles and hair.

### Multi-ply polypropylene (SMS)

This fabric consists of inner layers of meltblown polypropylene sandwiched between outer layers of spunbond polypropylene. Strong, fluid-resistant and breathable, this multilayer material provides an optimal mix of fluid protection and comfort.

### Polyethylene film

This single layer of polyethylene film provides high fluid-resistant properties at a very attractive price.

### Bilaminate polypropylene-polyethylene

This fabric is created by laminating a layer of polypropylene (SPP or SMS) with a layer of impervious polyethylene. This fabric is ideal for any situation demanding a material with high tensile strength and fluid protection.

### Trilaminate polypropylene-polyethylene

This impervious and breathable material consists of a hybrid film laminated between two layers of a bicomponent spunbond material. This fabric is ideal for long interventions when impervious protection is needed against fluids and chemotherapy drugs.

ISOLATION GOWNS OVERVIEW

Product Family	Code	Size	PF	Colour	Material	Length (cm)	Width (cm)	Wrist Type	Neck Closure	EN 20811: Resistance to Liquid Penetration (cmH2O)	Product Classification	Intended Use	Sterile/ Non-sterile			
Cover gowns	CRIE4000	L	50		Polypropylene (SPP)	114	135	Elastic	Ties	n/a	MDR Class I	Protect the patient against cross-contamination	Non-sterile			
	CRIE4000G															
	CRIE4000B															
	CRIE4001	XL				127	147							Knit cuffs		
	CRIE4001G															
	CRIE4001B															
	CRIE4001W															
	CRIE4001CF															
	CRIE4001BCF															
Examination gowns	NONE27SMS3	L	100		Multi-ply polypropylene (SMS)	114	142	Elastic	Ties	10	MDR Class I + PPE Category I	MDR: Protect the patient against cross-contamination	Non-sterile			
	NONE27SMS3XL	XL				127	152	Elastic	Tape tab							
	NONELV200	L				114	142									
	NONELV200XL	XL				127	152	Knit cuffs	Ties							
	NONELV200C	L				118	142									
	NONELV200CXL	XL				135	152	Knit cuffs	Ties							
	NONE386CF*	L				106	132									
	NONE386XLCF*	XL				118	142	Knit cuffs	Hook and loop							
	NONE326MCF	M				135	152									
	NONE326CF	L				135	152									
	NONE326XLCF	XL				114	142									
	NONE326XLCE	XL		50			135	152	Knit cuffs	Ties				52	PPE Category I: Protect the user against prolonged contact with water	Sterile
	NONE27280	L	25		114	142										
	NONE27280XL	XL			127	152										
	Impervious cover gowns	CRIE4040C*	L	100		Bilaminate polypropylene-polyethylene	114	135	Knit cuffs	Ties		109	PPE Category I	Protect the user against prolonged contact with water and cleansing agents	Non-sterile	
		CRIE4041C*	XL				135	152							Non-sterile	
CRIE4041CE*		XL	50				135	152			Sterile					
CRIE4040BC*		L	100		114		135	Knit cuffs			Ties			181	[CRIE4040xC range and CRIE500x are tested according to EN 16523-1 against cleansing agent glucoprotamin (e.g. Incidin® Plus) – breakthrough time of >480 minutes: performance level Class 6 according to EN 14325]	Non-sterile
CRIE4041BC*		XL			135		152									
CRIE4040WC*		L			114		135									
CRIE4041WC*		XL		135	152		Knit cuffs	Ties			50					
NONELV350		L	50		Microporous laminate	114			142							
CRIE5000		L	75		Polyethylene	114			94	Thumb loops		Over the head				
CRIE5001		XL				134	113									
Impervious isolation gowns	CRIE4020CF*	L	50		Bilaminate polypropylene-polyethylene	114	135	Knit cuffs	Ties	78	PPE Category I	Protect the user against prolonged contact with water	Non-sterile			
	CRIE4020CE*	L				114	135		Ties				Sterile			
	CRIE4021VCF*	L				114	135		Hook and loop	149						
	CRIE4022VCF*	XXL				140	147									
	CRIE4010CF	L				114	135	Knit cuffs	Ties	120						
	CRIE4011CF	XL				127	147		Ties							
	NONE27116	L				114	142	Knit cuffs	Ties	111						
	NONE27116XL	XL				127	142		Over the head							
	NONE272790HXL	XL				127	86									
Impervious isolation gowns against chemo agents	NONE27457	L	30		Trilaminate polypropylene-polyethylene	114	142	Knit cuffs	Hook and loop	200	PPE Category III	Protect the user against substances and mixtures that are hazardous to health (chemotherapy drugs)	Non-sterile			
	NONE27457XL	XL				127	152									

\*Isolation gown with back panel made of a breathable material

# COVER GOWNS

Suitable for preventing particulate contamination

## Lightweight polypropylene cover gowns

- » Breathable and comfortable
- » Waist and neck ties
- » Available with elastic wrists or knit cuffs

Item No.	Size	Colour	Packaging
Cover gowns with elastic wrists			
CRIE4000	L	Yellow	10/bag, 50/case
CRIE4000G	L	Green	10/bag, 50/case
CRIE4000B	L	Blue	10/bag, 50/case
CRIE4001	XL	Yellow	10/bag, 50/case
CRIE4001G	XL	Green	10/bag, 50/case
CRIE4001B	XL	Blue	10/bag, 50/case
CRIE4001W	XL	White	10/bag, 50/case
Cover gowns with knit cuffs			
CRIE4001CF	XL	Yellow	10/bag, 50/case
CRIE4001BCF	XL	Blue	10/bag, 50/case



CRIE4001G

# EXAMINATION GOWNS

Suitable for protection against minimal risks from prolonged contact with water

## Lightweight SMS examination gowns

- » Breathable and comfortable
- » Waist and neck ties
- » Available with elastic wrists or knit cuffs

Item No.	Size	Colour	Packaging
Neck and waist ties, elastic wrists			
NONE27SMS3	L	Blue	10/bag, 50/case
NONE27SMS3XL	XL	Blue	10/bag, 50/case
Tape tab neck, waist ties, elastic wrists			
NONELV200	L	Yellow	10/bag, 100/case
NONELV200XL	XL	Yellow	10/bag, 100/case
Neck and waist ties, knit cuffs			
NONELV200C	L	Yellow	10/bag, 100/case
NONELV200CXL	XL	Yellow	10/bag, 100/case



NONE27SMS3

## Medium-weight SMS examination gowns

- » The material used in this range of examination gowns has a hydrostatic head performance higher than 50 cmH<sub>2</sub>O according to EN 20811.

Item No.	Size	Colour	Sterile/ Non-sterile	Packaging
Neck and waist ties, knit cuffs, breathable back panel				
NONE386CF	L	Yellow	Non-sterile	10/bag, 100/case
NONE386XLCF	XL	Yellow	Non-sterile	10/bag, 100/case
Hook and loop at the neck, waist ties, knit cuffs				
NONE326MCF	M	Yellow	Non-sterile	10/bag, 100/case
NONE326CF	L	Yellow	Non-sterile	10/bag, 100/case
NONE326XLCF	XL	Yellow	Non-sterile	10/bag, 100/case
NONE326XLCE	XL	Yellow	Sterile	1/bag, 50/case



NONE326CF



# IMPERVIOUS COVER GOWNS

Suitable for protecting the user against prolonged contact with water and cleansing agents; CRIE404xC range and CRIE500X range tested according to EN 16523-1 against cleansing agent glucoprotamin (e.g. Incidin® Plus)

## Lightweight bilaminate impervious cover gowns

- » Waist and neck ties
- » Knit cuffs

Item No.	Size	Colour	Sterile/ Non-sterile	Packaging
CRIE4040C	L	Yellow	Non-sterile	10/bag, 100/case
CRIE4041C	XL	Yellow	Non-sterile	10/bag, 100/case
CRIE4041CE	XL	Yellow	Sterile	1/bag, 50/case
CRIE4040BC	L	Blue	Non-sterile	10/bag, 100/case
CRIE4041BC	XL	Blue	Non-sterile	10/bag, 100/case
CRIE4040WC	L	White	Non-sterile	10/bag, 100/case
CRIE4041WC	XL	White	Non-sterile	10/bag, 100/case



CRIE4040C

## Polyethylene impervious cover gowns

- » Over-the-head design
- » Thumb-loop-style wrists
- » Individually folded and wrapped in a dispenser box

Item No.	Size	Colour	Packaging
CRIE5000	L	Blue	1/bag, 15/box, 75/case
CRIE5001	XL	Blue	1/bag, 15/box, 75/case



CRIE5001

## Polyethylene disposable apron

- » Worn over non-fluid-resistant gowns or scrub suits as a convenient solution for users that want to add extra fluid protection at a lower cost

Item No.	Description	Packaging
NONE24274W	Medium weight, 71.1 x 116.8 cm, individually wrapped, dispenser box	50/box, 250/case



NONE24274W

## Microporous bilaminate cover gown

- » Waist and neck ties
- » Knit cuffs

Item No.	Size	Colour	Packaging
NONELV350	L	White	10/bag, 50/case



NONELV350

# IMPERVIOUS ISOLATION GOWNS

Suitable for protection against prolonged contact with water

## Medium-weight bilaminate impervious isolation gowns with a breathable back

Item No.	Size	Colour	Sterile/ Non-sterile	Packaging
Waist and neck ties, knit cuffs				
CRIE4020CF	L	Yellow	Non-sterile	10/bag, 50/case
CRIE4020CE	L	Yellow	Sterile	1/bag, 50/case
Hook and look at the neck, waist ties, knit cuffs				
CRIE4021VCF	L	Yellow	Non-sterile	10/bag, 50/case
CRIE4022VCF	XL	Yellow	Non-sterile	10/bag, 50/case



CRIE4021VCF

## Medium-weight bilaminate impervious isolation gowns with a breathable back

- » Waist and neck ties
- » Knit cuffs

Item No.	Size	Colour	Packaging
CRIE4010CF	L	Yellow	10/bag, 50/case
CRIE4011CF	XL	Yellow	10/bag, 50/case
NONE27116	L	Blue	10/bag, 50/case
NONE27116XL	XL	Blue	10/bag, 50/case
NONE27279OHL*	XL	Blue	10/bag, 50/case



NONE27116

\*Over-the-head design and waist ties

# CHEMO PROTECTION ISOLATION GOWNS

Trilaminare polypropylene-polyethylene (Prevention Plus) gowns are certified as PPE Category III Type PB-[4] for partial body protection according to PPE Regulation 2016: 425 and are suitable for protection against mixtures that are hazardous to health. The Prevention Plus material is very comfortable due to its breathability properties.

## Tested against 14 chemotherapy drugs per ASTM F73

- Carmustine (BiCNU)

Dacarbazine

Fluorouracil (5FU)

Thiotepa

Cisplatin

Mitomycin C

Cyclophosphamide (Cytosan)
- Etoposide (Toposar)

Vincristine Sulphate

Doxorubicin hydrochloride

Paclitaxel (Taxol)

Leucovorin

Docetaxel

Oxaliplatin

Test data is available on request.

- » Hook-and-loop closure at the neck
- » Waist ties
- » Knit cuffs at the wrists

Item No.	Size	Colour	Packaging
NONE27457	L	Blue	10/bag, 30/case
NONE27457XL	XL	Blue	10/bag, 30/case

Complete protection is available in combination with the trilaminare polypropylene-polyethylene boot covers (NONE27348PXL).



NONE27457: gown  
NONE27348PXL: boot covers

# FOOTWEAR<sup>1</sup>

## Shoe covers

Footwear covers prevent environmental contamination and reduce the transmission of microorganisms. Manufactured from polyethylene or polypropylene materials, some models feature additional fluid protection.

Item No.	Size	Material	Characteristics	Packaging
CRIE2002	One size fits most	Polypropylene (SPP)	Protection	100/box, 300/case
CRIE2010	One size fits most	Polyethylene film	Impervious	100/bag, 1,000/case
CRIE2102W	One size fits most	Polypropylene (SPP)	Protection	500/box, 1,000/case
CRIE2020	One size fits most	Polyethylene film and polypropylene (SPP)	Impervious	100/bag, 200/case
NONE28858	L	Multi-ply polypropylene (SMS)	Fluid-resistant	100/box, 300/case
NONE28859	XL		Fluid-resistant	100/box, 200/case



## Boot covers

Medline's boot cover range offers enhanced coverage and increased fluid protection. The products cover the shoe and leg all the way up to the knee and have a convenient one-size-fits-most design with a textured sole on the bottom. They also feature a hook-and-loop strap or elastic bands around the ankle to help keep them in place and ensure adequate protection.

Item No.	Size	Characteristics	Packaging
NONE27143XL	One size fits most	Fluid-resistant	50/box, 150/case
NONE27144	One size fits most	Impervious	50/box, 150/case
NONE27348PXL	One size fits most	Breathable and impervious	50/box, 150/case

### NONE27143XL

- » Breathable, fluid-resistant material above the ankle
- » High fluid protection on the bottom
- » Sealed seams on the foot



NONE27143XL

### NONE27144

- » Impervious PE-coated polypropylene material
- » High fluid protection on the bottom
- » Hook-and-loop strap for an adjustable fit



NONE27144

### NONE27348PXL

- » Made of blue Prevention Plus material
- » Hook-and-loop closure
- » Foam strip



NONE27348PXL

# DISPOSABLE SCRUB SUITS<sup>2</sup>

Medline's disposable scrub suits are intended to be used by healthcare professionals in hospitals and clinics in order to maintain a high level of hygiene and contribute to infection control. These scrub suits offer your staff comfort because they are manufactured from a breathable and soft SMS material. We offer two ranges of scrub suits: Soft and Advanced. Both ranges have been designed in collaboration with healthcare professionals for your comfort and safety.

## Features of Medline scrub suits



Bell-shaped design  
for a perfect fit



3 convenient pockets  
on the shirt



V-necked



Adjustable waist ties



Shirt and pants packed  
separately or in sets

## Soft Scrub Suits

- » Manufactured from 35-gsm SMS material with optimal opacity and softness
- » Low linting, with an anti-static treatment and fluid resistance
- » Available in 7 sizes from XS to 3XL
- » Available in **blue** and **green** with shirt and pants packed separately or in sets
- » Latex free



## Advanced Scrub Suits

- » Manufactured from 42-gsm SMS material with optimal opacity and softness
- » Low linting, with an anti-static treatment and fluid resistance
- » Available in 7 sizes from XS to 3XL
- » Available in **blue**, **green** and **purple** with shirt and pants packed separately or in sets
- » Latex free





Soft Scrub Suits packed as a set

Size	Blue Scrub Suits	Green Scrub Suits	Packaging
XS	P35PBXS	P35PGXS	1/bag, 25/case
S	P35PBS	P35PGS	1/bag, 25/case
M	P35PBM	P35PGM	1/bag, 25/case
L	P35PBL	P35PGL	1/bag, 25/case
XL	P35PBXL	P35PGXL	1/bag, 25/case
2XL	P35PB2XL	P35PG2XL	1/bag, 25/case
3XL	P35PB3XL	P35PG3XL	1/bag, 25/case

Soft Scrub Suits packed separately

Size	Blue Shirts	Blue Pants	Green Shirts	Green Pants	Packaging
XS	PXS9110S	PXS9110P	P35PGXS-S	P35PGXS-P	1/bag, 50/case
S	PS9110S	PS9110P	P35PGS-S	P35PGS-P	1/bag, 50/case
M	PM9110S	PM9110P	P35PGM-S	P35PGM-P	1/bag, 50/case
L	PL9110S	PL9110P	P35PGL-S	P35PGL-P	1/bag, 50/case
XL	PXL9110S	PXL9110P	P35PGXL-S	P35PGXL-P	1/bag, 50/case
2XL	P2XL9110S	P2XL9110P	P35PG2XL-S	P35PG2XL-P	1/bag, 50/case
3XL			P35PG3XL-S	P35PG3XL-P	1/bag, 50/case

Advanced Scrub Suits packed as a set

Size	Blue Scrub Suits	Green Scrub Suits	Purple Scrub Suits	Packaging
XS	P42PBXS	P42PGXS	P42PPXS	1/bag, 25/case
S	P42PBS	P42PGS	P42PPS	1/bag, 25/case
M	P42PBM	P42PGM	P42PPM	1/bag, 25/case
L	P42PBL	P42PGL	P42PPL	1/bag, 25/case
XL	P42PBXL	P42PGXL	P42PPXL	1/bag, 25/case
2XL	P42PB2XL	P42PG2XL	P42PP2XL	1/bag, 25/case
3XL	P42PB3XL	P42PG3XL	P42PP3XL	1/bag, 25/case

Advanced Scrub Suits packed separately

Size	Blue Shirts	Blue Pants	Green Shirts	Green Pants	Purple Shirts	Purple Pants	Packaging
XS	P42PBXS-S	P42PBXS-P	P42PGXS-S	P42PGXS-P	P42PPXS-S	P42PPXS-P	1/bag, 50/case
S	P42PBS-S	P42PBS-P	P42PGS-S	P42PGS-P	P42PPS-S	P42PPS-P	1/bag, 50/case
M	P42PBM-S	P42PBM-P	P42PGM-S	P42PGM-P	P42PPM-S	P42PPM-P	1/bag, 50/case
L	P42PBL-S	P42PBL-P	P42PGL-S	P42PGL-P	P42PPL-S	P42PPL-P	1/bag, 50/case
XL	P42PBXL-S	P42PBXL-P	P42PGXL-S	P42PGXL-P	P42PPXL-S	P42PPXL-P	1/bag, 50/case
2XL	P42PB2XL-S	P42PB2XL-P	P42PG2XL-S	P42PG2XL-P	P42PP2XL-S	P42PP2XL-P	1/bag, 50/case
3XL	P42PB3XL-S	P42PB3XL-P	P42PG3XL-S	P42PG3XL-P	P42PP3XL-S	P42PP3XL-P	1/bag, 50/case

WARM-UP JACKETS

Keep your staff warm in the cooler open areas of the hospital with Medline’s comfortable range of warm-up jackets. These lab jackets are designed to have a professional look. And the range includes jackets with long, three-quarter length and short sleeves to meet your hygiene requirements.

- » Made of a comfortable SMS material
- » 5-snap closure
- » 3 reinforced convenient pockets
- » Gentle white knit collar and knit cuffs for a professional look

Long-sleeved range

Item No.	Description	Packaging
NONERP600S	Blue, S	10/bag, 30/case
NONERP600M	Blue, M	10/bag, 30/case
NONERP600L	Blue, L	10/bag, 30/case
NONERP600XL	Blue, XL	10/bag, 30/case



NONERP600S

Long-sleeved range with elastic at the waist

Item No.	Description	Packaging
NONERP610S	Blue, S	10/bag, 50/case
NONERP610M	Blue, M	10/bag, 50/case
NONERP610L	Blue, L	10/bag, 50/case
NONERP610XL	Blue, XL	10/bag, 50/case



NONERP800S

NONERP860M

Three-quarter length sleeves range

Item No.	Description	Packaging
NONERP800S	Blue, S	10/bag, 30/case
NONERP800M	Blue, M	10/bag, 30/case
NONERP800L	Blue, L	10/bag, 30/case
NONERP800XL	Blue, XL	10/bag, 30/case
NONERP860M	White, M	10/bag, 30/case
NONERP860L	White, L	10/bag, 30/case
NONERP860XL	White, XL	10/bag, 30/case



NONERP700S

Short-sleeved range

Item No.	Description	Packaging
NONERP700M	Blue, M	10/bag, 30/case
NONERP700L	Blue, L	10/bag, 30/case
NONERP700XL	Blue, XL	10/bag, 30/case



# HEADWEAR<sup>3</sup>

Covering your head and hair helps reduce the risk of foreign particles getting into surgical wounds and protects patients from potential contamination in operating theatres or ICUs. Medline's disposable headwear products offer a high level of comfort and come in a variety of shapes, sizes and designs that address hospitals' protective and aesthetical needs.

Our broad product range, including surgical hoods and caps, bonnets and bouffant caps, is suitable for a variety of operations and hairstyles. To offer you even more choice and meet your budget, Medline's headwear is available in cluster keyback, viscose and spunbond polypropylene materials that have passed biocompatibility requirements, are low linting and do not contain latex.

## Comfort Line

The Comfort Line is Medline's preferred headwear product range for the operating theatre. Made of soft and breathable cluster keyback fabric, this line offers both a high level of comfort and a high level of protection. Comfort Line products are available in blue and green and a large variety of designs. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why our Comfort product range is produced in a WRAP-certified facility.



### Surgical Hood

FS61100

100/box  
600/case

Surgical hood · green



### Surgical Hood

FS61150A

100/box  
600/case

Surgical hood · blue



### Surgical Hood

FS61250

50/box  
300/case

Surgical hood with sweatband · blue



### Bouffant Cap

FS61300

100/box  
600/case

Bouffant cap with a latex-free elastic band · green · size 53 cm



### Bouffant Cap

FS61350

100/box  
600/case

Bouffant cap with a latex-free elastic band · blue · size 53 cm



### Surgical Cap

FS61450

100/box  
600/case

Surgical cap · sailor style · blue



### Surgical Cap

FS61500

100/box  
600/case

Surgical cap · elastic at the back · contoured front · green



### Surgical Cap

FS61550A

100/box  
600/case

Surgical cap · elastic at the back · contoured front · blue



### Surgical Cap

FS61750A

100/box  
600/case

Surgical cap with ties · for long hair · blue



### Surgical Cap

FS61850A

100/box  
600/case

Surgical cap with ties · for short hair · blue



### Surgical Cap


FS61950

100/box  
600/case

Surgical cap · elastic at the back · flat front · blue

Colourful Line


The Colourful Line is made of spunlace viscose, which is a soft, breathable and cotton-like material. The range comes in a variety of colours and offers good comfort. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why our Colourful product range is produced in a WRAP-certified facility.



**Surgical Cap**

**FS60750B** 100/box  
600/case


Surgical cap with ties · for long hair · blue



**Surgical Cap**

**FS60850B** 100/box  
600/case


Surgical cap with ties · for short hair · blue



**Surgical Cap**

**FS60950B** 100/box  
600/case


Surgical cap · elastic at the back · flat front · blue



**Surgical Cap**

**FS60950G** 100/box  
600/case


Surgical cap · elastic at the back · flat front · green



**Surgical Cap**

**FS60950P** 100/box  
600/case


Surgical cap · elastic at the back · flat front · purple



**Surgical Cap**

**FS60950PK** 100/box  
600/case

Surgical cap · elastic at the back · flat front · pink




**Surgical Cap**

**FS60950W** 100/box  
600/case

Surgical cap · elastic at the back · flat front · white

Basic Line


Designed to satisfy the essential needs of the ward at a convenient price, the Basic Line is made of a spunbond material with a basis weight from 12 gsm to 30 gsm. Basic Line headwear comes in blue, green, white or printed patterns.



**Surgical Cap**

**NONE28625** 100/box  
500/case

Surgical cap with ties · SMS top, absorbent scrim-reinforced material on sides · one size fits most · light blue



**Surgical Cap**

**FS67050** 100/box  
600/case


Spunlace bonnet · one size fits most · blue



**Surgical Cap**

**NONE28233F** 100/box  
500/case


Bouffant cap with butterfly print · size 60 cm · blue



**Surgical Hood**

**FS63100A** 100/box  
600/case

Surgical hood · green



**Surgical Hood**

**FS63150** 100/box  
600/case

Surgical hood · green



**Surgical Hood**

**FS63150A** 100/box  
600/case

Light surgical hood · blue



**Bouffant Cap**

**FS63300** 100/box  
600/case

Bouffant cap · size 53 cm · green



**Bouffant Cap**

**FS63350** 100/box  
600/case

Bouffant cap · size 53 cm · blue



**Bouffant Cap**

FS63360

100/box  
600/case

Bouffant cap • size 53 cm • white



**Surgical Cap**

FS63450A

100/box  
600/case

Surgical cap • sailor style • blue



**Clip Cap**

FS65360

100/box  
1,000/case

Clip cap • one size fits most • white



**Clip Cap**

FS65450

100/box  
1,200/case

Clip cap • one size fits most • blue



**Clip Cap**

FS65460

100/box  
1,200/case

Clip cap • one size fits most • white



**Bouffant Cap**

CRIE1000

100/box  
500/case

Bouffant cap • size 53 cm • green



**Bouffant Cap**

CRIE1002

100/box  
500/case

Bouffant cap • size 53 cm • white



**Bouffant Cap**

CRIE1003

100/box  
500/case

Bouffant cap • size 60 cm • white



**Bouffant Cap**

CRIE1004

100/box  
500/case

Bouffant cap • size 60 cm • blue



**Bouffant Cap**

CRIE1005

100/box  
500/case

Bouffant cap • size 60 cm • green

# Choose your headwear



## Surgical hoods

Offering a helmet-like cap with ties at the neck, one-size-fits-all surgical hoods provide the most coverage of all products in our headwear range for worry-free performance during surgery.

Some products feature an absorbent sweatband for increased staff comfort, especially for longer procedures.

## Ordering Information

Item No.	Description	Packaging
FS61100	Cluster keyback, green	100/box, 600/case
FS61150A	Cluster keyback, blue	100/box, 600/case
FS61250	Cluster keyback, with sweatband, blue	50/box, 300/case
FS63100A	Spunbond polypropylene 30 gsm, green	100/box, 600/case
FS63150A	Spunbond polypropylene 30 gsm, blue	100/box, 600/case
FS63150	Spunbond polypropylene 20 gsm, blue	100/box, 600/case



## Surgical caps with ties for long and short hair

Featuring ties at the back of the neck, these adjustable surgical caps are designed for people with long or short hair to safely keep the hair covered.

This surgical cap is ideal for use in the operating theatre or for day surgeries.

Item No.	Description	Packaging
<b>For long hair</b>		
FS61750A	Cluster keyback, blue	100/box, 600/case
FS60750B	Viscose, blue	100/box, 600/case
<b>For short hair</b>		
FS61850A	Cluster keyback, blue	100/box, 600/case
FS60850B	Viscose, blue	100/box, 600/case
NONE28625	SMS and scrim-reinforced material, blue	100/box, 500/case



## Surgical caps with elastic at the back

Featuring a broad shape and elastic at the back, these one-size-fits-most surgical caps provide a comfortable and secure fit, meeting the requirements of the most active of days.

This cap design is well-suited for long hairstyles and long periods of wear in the operating theatre, ICU and other laboratory departments.

Item No.	Description	Packaging
FS61500	Cluster keyback, green	100/box, 600/case
FS61550A	Cluster keyback, blue	100/box, 600/case
FS61950	Cluster keyback, blue	100/box, 600/case
FS60950B	Viscose, blue	100/box, 600/case
FS60950G	Viscose, green	100/box, 600/case
FS60950P	Viscose, purple	100/box, 600/case
FS60950PK	Viscose, pink	100/box, 600/case
FS60950W	Viscose, white	100/box, 600/case



# PATIENT GOWNS

Improving patient care is a priority for all healthcare providers. Most hospitals and care professionals take this goal a step further and have the overall objective to achieve a high level of patient-centred care and satisfaction.

At Medline, we think the quality of the product contributes to the quality of care that patients receive during their time in the hospital.

Therefore, our patient gowns provide your patients with opacity, comfort and modesty during hospital stays because they are manufactured from a soft, non-woven SMS material.

Medline offers a wide variety of designs to meet the needs of your department. Choose from a selection of patient gowns with front openings, back openings or hook-and-loop fasteners on the sleeves.

Item No.	Size	Description	Packaging
NONE27046 NONE27046XL	L XL	Sleeveless patient gown, blue	50/case
NONE27046SL	L	Patient gown with short sleeves, blue	50/case
NONE27246SLV	L	Patient shirt with hook and loops on sleeves for removal, blue	50/case
NONE27SLK	Suitable for children (4 to 9 years)	Patient gown for child with short sleeves, blue	100/case
CRIE4060K	L	Kimono gown with long sleeves, waist ties, blue	60/case
NONE27036	Standard size	Patient cape with tape tab at the neck, dark blue	100/case
NONE27100	L	Patient gown with front opening and short sleeves, waist ties, blue	100/case



# PATIENT APPAREL

Our apparel provides patients with comfort and dignity during their day surgery or diagnostic procedure because it is manufactured from a soft and opaque polypropylene material. Discover Medline’s individual components and choose between colonoscopy shorts, underwear or elastic slippers to offer your patients the very best care and comfort.

Item No.	Size	Description	Packaging
NONE27UP NONE27UPXXL	L XXL	Underwear with elastic, dark blue	200/case
NONECSHL NONECSH2XL	L XXL	Colonoscopy shorts with back opening and hook-and-loop closures, blue	100/case
CRIE2210	Standard size (32 cm length)	Elastic slippers with rigid sole, white	100/case
MDTEFP218RXL	XL	Socks, red	48/case




# PATIENT KITS<sup>4</sup>

Our expertise in medical apparel design and manufacturing allows us to offer new solutions for day surgeries. With efficiency in mind, Medline has created a variety of patient kit solutions and a visitor kit to suit your needs and budget.

The kits contain single-use patient apparel, which patients can wear during the diagnostic procedure or day surgery. All components are conveniently packed in a single bag that can be given to the patient upon arrival. Due to its packaging, organising storage and replenishing stock are easy to handle, and you only have to order one item instead of four or five individual components.

As a leading international manufacturer and distributor of medical devices, we understand your focus and your patients’ needs. We apply the highest quality standards to our products along each step of the manufacturing and assembly process.



Item No.	Size	Description	Packaging
KIT1030	L	Kimono, blue Patient gown with short sleeves, blue Underwear, dark blue Clip cap, blue Socks, red	40/case
KIT1035 KIT1035XXL	L XXL	Kimono, purple Patient gown with short sleeves, purple Underwear, dark blue Clip cap, blue Socks, red	40/case
KIT1040 KIT1040XXL	L XXL	Patient gown with short sleeves, front opening, blue Scrub suit pants, blue Clip cap, blue Socks, red	40/case
KIT3020 KIT3020XXL	L XXL	Patient gown with short sleeves, blue Underwear, dark blue Clip cap, white Shoe covers, white	50/case
KIT2020 KIT2020XXL	L XXL	Patient gown with short sleeves, blue Clip cap, white Underwear, dark blue Elastic slippers with rigid sole, white	40/case
NONE27KITK	Suitable for children (4 to 9 years)	Patient gown for child with short sleeves, blue Underwear, dark blue Clip cap, white Shoe covers, white	50/case
NONE35LKIT NONE352XLKIT	L XXL	Patient gown with short sleeves, blue Scrub suit pants, blue Clip cap, white Shoe covers, textured bottom, white and blue	50/case
CRIE4000BKIT	L	Cover gown with elastic wrists, blue Facemask Type II with ear loops, blue Clip cap, white Shoe covers, white	50/case



# MEDLINE ALSO OFFERS A LARGE RANGE OF SURGICAL, EXAMINATION AND RESPIRATORY MASKS.

CHECK OUT THE FACEMASK BROCHURE FOR MORE DETAILS.

## Facemasks



For more information on our products, please contact your account manager or visit us at: [www.medline.eu/uk](http://www.medline.eu/uk)



**ALWAYS  
ON.**

**Medline Industries Ltd**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
Fax: +44 844 334 5238  
[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

<sup>1)</sup> Footwear are class I medical devices according to Regulation (EU) 2017/745.

<sup>2)</sup> Disposable scrub suits are class I medical devices according to Regulation (EU) 2017/745.

<sup>3)</sup> Headwear products are class I medical devices according to Regulation (EU) 2017/745.

<sup>4)</sup> Patient kits are procedure packs according to Article 22 of Regulation (EU) 2017/745.

The procedure pack producer is Medline International France SAS, 5 Rue Charles Lindbergh, 44110 Châteaubriant, France.

These products are intended to be used by healthcare professionals.  
Before use, please consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct any errors that may have occurred.

© 2019 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML179-EN 09/2019



Medline International France SAS  
5, rue Charles Lindbergh  
44110 Châteaubriant, France

