

Protective Apparel

Enhance safety and protection




**Reduce cross-
contamination
and protect
your staff.**



Why are more healthcare professionals choosing Medline for their all-around protection needs?



High-quality products



Comprehensive product range



Reliable partner

Your staff are exposed to bacteria and bodily fluids on a daily basis, putting them at risk of cross-contamination. We know you rely on disposable apparel to keep them protected, but also comfortable. Therefore, we apply the highest quality standards to our products and in our manufacturing facilities to ensure that healthcare professionals get the protection they deserve.

Medline International France is EN ISO 13485 and ISO 9001 certified and complies with regulatory standards on the market. In addition, we are conscious about social responsibility and, as a result, some of our products are WRAP-certified. WRAP (Worldwide Responsible Accredited Production) is an internationally recognised social compliance certification body that verifies whether factories respect workers' labour and human rights.

Table of contents

European regulations and directives	3
Guide to protective materials	3
Isolation gowns overview	4-5
Cover gowns	6
Isolation gowns	6-7
Impervious cover gowns	7-8
Medium-weight bilaminate impervious isolation gowns	8
Chemo protection isolation gowns	9
Footwear	10
Disposable scrub suits	11-14
Warm-Up Jackets	15
Headwear	16-21
Patient gowns	22
Patient apparel	22
Patient accessories	22
Patient kits	23



It is critical for our customers to receive protective apparel on time. That is why we have selected certain items where we have built a large safety stock in order to react faster to any unexpected increase in your demand. **These items are marked in bold for easy identification.**

European regulations and directives to classify protective apparel

Regulation (EU) 2016/425 on personal protective equipment (PPE)	Regulation (EU) 2017/745 on medical devices (MDR)	Directive 2001/95/EC on general product safety
The personal protective equipment included in Medline's range provides protection against minimal risks (contact with weak-action cleaning materials or prolonged contact with water) or more serious risks, such as exposure to substances and mixtures that are hazardous to health. Products listed in the brochure that are certified as PPE are compliant with PPE Category I or Category III requirements.	Medical devices include any instrument, apparatus, appliance, software, implant, reagent, material or other article intended by the manufacturer to be used, alone or in combination, for human beings for a specific medical purpose. The medical devices included in Medline's range are intended to protect the patient from potential contamination. They are Class I devices.	Products included in Medline's protective apparel range that do not fall under the definition of PPE or MDR are certified under the general product safety directive. These are general commodity products (e.g. underwear).

Please consult the product's technical data sheet to see its classification.

Guide to protective materials

To help you choose the right product, here are brief descriptions of the features and benefits of the main materials available from Medline. You can refer back to this page as you look through the brochure.

Spunbond polypropylene (SPP)	This standard fabric is made by bonding fibres together to form a single layer of breathable woven-like material. It is suitable for protecting users against mechanical particles and hair.
Multi-ply polypropylene (SMS)	This fabric consists of inner layers of meltblown polypropylene sandwiched between outer layers of spunbond polypropylene. Strong, fluid-resistant and breathable, this multilayer material provides an optimal mix of fluid protection and comfort.
Polyethylene film	This single layer of polyethylene film provides high fluid-resistant properties at a very attractive price.
Bilaminate polypropylene-polyethylene	This fabric is created by laminating a layer of polypropylene (SPP or SMS) with a layer of impervious polyethylene. This fabric is ideal for any situation demanding a material with high tensile strength and fluid protection.
Trilaminate polypropylene-polyethylene	This impervious and breathable material consists of a hybrid film laminated between two layers of a bicomponent spunbond material. This fabric is ideal for long interventions when impervious protection is needed against fluids and chemotherapy drugs.

Isolation gowns overview

Product family	Code	Size	PF	Colour	Material	Length (cm)	Width (cm)	Wrist type	Neck closure	ISO811: Resistance to liquid penetration (cmH2O)	Product classification	Intended use	Sterile/ non-sterile
Cover gowns	CRIE4000	L	50		Polypropylene (SPP)	114	135	Elastic	Ties	n/a	MDR Class I	Protect the patient against cross-contamination	Non-sterile
	CRIE4000G												
	CRIE4000B												
	CRIE4001	XL				127	147	Knit cuffs					
	CRIE4001G												
	CRIE4001B												
	CRIE4001CF												
	CRIE4001BCF												
Isolation gowns	NONE27SMS3	L	100		Multi-ply polypropylene (SMS)	114	142	Elastic	Ties	10	PPE Category I	PPE Category I: Protect the user against prolonged contact with water	Non-sterile
	NONE27SMS3XL	XL				127	152						
	NONELV200C	L				114	142	Knit cuffs	Ties				
	NONELV200CXL	XL				127	152						
	NONE326MCF	M				106	132	Knit cuffs	Hook and loop	50.8			
	NONE326CF	L				118	142						
	NONE326XLCF	XL				135	152						
	NONE27280	L	25			114	142	Knit cuffs	Ties	54.5		Sterile	
	NONE27280XL	XL			127	152							
Impervious cover gowns	CRIE4040C*	L	100		Bilaminate polypropylene-polyethylene	114	135	Knit cuffs	Ties	107	PPE Category I	Protect the user against prolonged contact with water and cleansing agents [CRIE4040xC range and CRIE500x are tested according to EN 16523-1 against the cleansing agent glucoprotamin – breakthrough time of >480 minutes: performance level Class 6 according to EN 14325]	Non-sterile
	CRIE4041C*	XL				135	152						
	CRIE4042C*	2XL				147	162						
	CRIE4040BC*	L	100			114	135						
	CRIE4041BC*	XL				135	152						
	NONELV350	L	50		Microporous laminate	114	142	Knit cuffs	Ties	181			
	CRIE5010	L	75		Polyethylene	114	94	Thumb loops	Over the head	75.5			
	CRIE5011	XL				134	113			99			
	CRIE5021	XL				140	110						
Impervious isolation gowns	NONE27116	L	50		Bilaminate polypropylene-polyethylene	114	142	Knit cuffs	Ties	114	PPE Category I	Protect the user against prolonged contact with water	Non-sterile
	NONE27116XL	XL				127	142		Over the head				
	NONE272790HXL	XL				127	86						
Impervious isolation gowns against chemo agents	NONE27457	L	30		Trilaminate polypropylene-polyethylene	114	142	Knit cuffs	Hook and loop	200	PPE Category III	Protect the user against substances and mixtures that are hazardous to health (chemotherapy drugs)	Non-sterile
	NONE27457XL	XL				127	152						

*Isolation gown with back panel made of a breathable material

Cover gowns

Suitable for preventing particulate contamination

Lightweight polypropylene cover gowns

- Breathable and comfortable
- Waist and neck ties
- Available with elastic wrists or knit cuffs

Item number	Size	Colour	Packaging
Cover gowns with elastic wrists			
CRIE4000	L	Yellow	10/bag, 50/case
CRIE4000G	L	Green	10/bag, 50/case
CRIE4000B	L	Blue	10/bag, 50/case
CRIE4001	XL	Yellow	10/bag, 50/case
CRIE4001G	XL	Green	10/bag, 50/case
CRIE4001B	XL	Blue	10/bag, 50/case
Cover gowns with knit cuffs			
CRIE4001CF	XL	Yellow	10/bag, 50/case
CRIE4001BCF	XL	Blue	10/bag, 50/case



CRIE4001G

Isolation gowns

Suitable for protection against minimal risks from prolonged contact with water

Lightweight SMS isolation gowns

- Breathable and comfortable
- Waist and neck ties
- Available with elastic wrists or knit cuffs

Item number	Size	Colour	Packaging
Neck and waist ties, elastic wrists			
NONE27SMS3	L	Blue	10/bag, 50/case
NONE27SMS3XL	XL	Blue	10/bag, 50/case
Neck and waist ties, knit cuffs			
NONELV200C	L	Yellow	10/bag, 100/case
NONELV200CXL	XL	Yellow	10/bag, 100/case



NONE27SMS3

Medium-weight SMS isolation gowns

The material used in this range of SMS isolation gowns has a hydrostatic head performance higher than 50 cmH2O according to ISO 811.

Item number	Size	Colour	Sterile/ non-sterile	Packaging
Hook and loop at the neck, waist ties, knit cuffs				
NONE326MCF	M	Yellow	Non-sterile	10/bag, 100/case
NONE326CF	L	Yellow	Non-sterile	10/bag, 100/case
NONE326XLCF	XL	Yellow	Non-sterile	10/bag, 100/case



NONE326CF

Impervious cover gowns

Suitable for protecting the user against prolonged contact with water and cleansing agents; CRIE404xC range and CRIE501x range tested according to EN 16523-1 against the cleansing agent glucoprotamin

Lightweight bilaminate impervious cover gowns

- Waist and neck ties
- Knit cuffs

Item number	Size	Colour	Sterile/ non-sterile	Packaging
CRIE4040C	L	Yellow	Non-sterile	10/bag, 100/case
CRIE4041C	XL	Yellow	Non-sterile	10/bag, 100/case
CRIE4042C	2XL	Yellow	Non-sterile	10/bag, 100/case
CRIE4040BC	L	Blue	Non-sterile	10/bag, 100/case
CRIE4041BC	XL	Blue	Non-sterile	10/bag, 100/case



CRIE4040C

Polyethylene impervious cover gowns

- Over-the-head design
- Thumb-loop-style wrists
- Individually folded and wrapped in a dispenser box

Item number	Size	Colour	Packaging
CRIE5010	L	Blue	1/bag, 15/box, 75/case
CRIE5011	XL	Blue	1/bag, 15/box, 75/case
CRIE5021	XL	Blue	1/bag, 15/box, 75/case



CRIE5010

Polyethylene disposable apron

Worn over non-fluid-resistant gowns or scrub suits as a convenient solution for users that want to add extra fluid protection at a lower cost.

Item number	Description	Packaging
NONE24274W	Medium weight, 71.1 x 116.8 cm, individually wrapped, dispenser box	50/box, 250/case

Microporous bilaminate cover gown

- Waist and neck ties
- Knit cuffs

Item number	Size	Colour	Packaging
NONELV350	L	White	10/bag, 50/case



Medium-weight bilaminate impervious isolation gowns

Suitable for protection against prolonged contact with water

- Waist and neck ties
- Knit cuffs

Item number	Size	Colour	Packaging
NONE27116	L	Blue	10/bag, 50/case
NONE27116XL	XL	Blue	10/bag, 50/case
NONE27279OHXL*	XL	Blue	10/bag, 50/case

*Over-the-head design and waist ties



NONE27116

Chemo protection isolation gowns

Trilaminate polypropylene-polyethylene (Prevention Plus) gowns are certified as PPE Category III Type PB-[4] for partial body protection according to PPE Regulation 2016: 425 and are suitable for protection against mixtures that are hazardous to health. The Prevention Plus material is very comfortable due to its breathability properties.

- Hook-and-loop closure at the neck
- Waist ties
- Knit cuffs at the wrists

Tested against 14 chemotherapy drugs per ASTM F739

Carmustine (BiCNU)	Etoposide (Toposar)
Dacarbazine	Vincristine Sulphate
Fluorouracil (5FU)	Doxorubicin hydrochloride
Thiotepa	Paclitaxel (Taxol)
Cisplatin	Leucovorin
Mitomycin C	Docetaxel
Cyclophosphamide (Cytoxan)	Oxaliplatin

Test data is available on request.

Item number	Size	Colour	Packaging
NONE27457	L	Blue	10/bag, 30/case
NONE27457XL	XL	Blue	10/bag, 30/case

Complete protection is available in combination with the trilaminate polypropylene-polyethylene boot covers (NONE27348PXL).



NONE27457: gown
NONE27348PXL: boot covers

Footwear¹

Shoe covers

Footwear covers prevent environmental contamination and reduce the transmission of microorganisms. Manufactured from polyethylene or polypropylene materials, some models feature additional fluid protection.

Item number	Size	Material	Characteristics	Packaging
CRIE2002	One size fits most	Polypropylene (SPP)	Protection	100/box, 300/case
CRIE2010	One size fits most	Polyethylene film	Impervious	100/bag, 1,000/case
CRIE2102W	One size fits most	Polypropylene (SPP)	Protection	500/box, 1,000/case
CRIE2020	One size fits most	Polyethylene film and polypropylene (SPP)	Impervious	100/bag, 200/case
NONE28858	L	Multi-ply polypropylene (SMS)	Fluid-resistant	100/box, 300/case



Boot covers

Medline's boot cover range offers enhanced coverage and increased fluid protection. The products cover the shoe and leg all the way up to the knee and have a convenient one-size-fits-most design with a textured sole on the bottom. They also feature a hook-and-loop strap or elastic bands around the ankle to help keep them in place and ensure adequate protection.

Item number	Size	Characteristics	Packaging
NONE27143XL	One size fits most	Fluid-resistant	50/box, 150/case
NONE27144	One size fits most	Impervious	50/box, 150/case
NONE27348PXL	One size fits most	Breathable and impervious	50/box, 150/case

NONE27143XL

- Breathable, fluid-resistant material above the ankle
- High fluid protection on the bottom
- Sealed seams on the foot



NONE27144

- Impervious coated polyethylene-polypropylene material
- High fluid protection on the bottom
- Hook-and-loop strap for an adjustable fit



NONE27348PXL

- Made of blue Prevention Plus material
- Hook-and-loop closure
- Foam strip



Disposable scrub suits²

A scrub suit designed with you in mind

We have worked with healthcare professionals like you to create our scrub suits with features you want and need to perform at your best. Experience new, clean apparel every day.



Carefully selected materials

Medline's disposable scrub suits offer a high level of comfort thanks to their ultra-soft and breathable SMS material with the following features:



Soft Scrub Suits range

Essential apparel to meet daily challenges

- Manufactured from 35-gsm SMS material with optimal opacity and softness
- Available in 7 sizes from XS to 3XL
- Available in **blue** and **green** with shirt and pants packed separately or in sets



Advanced Scrub Suits range

The choice you deserve

- Manufactured from 42-gsm SMS material with great opacity and optimal softness
- Available in 7 sizes from XS to 3XL
- Available in **blue**, **green** and **purple** with shirt and pants packed separately or in sets



Ultimate Scrub Suits range

Our top-of-the-line ultra-soft range for demanding activities

Medline's Ultimate Scrub Suits offer that extra comfort, breathability and protection you and your patients deserve.

- Manufactured from 45-gsm SMS material with excellent opacity and softness
- Available in 7 sizes from XS to 3XL
- Available in **blue** and **green** with shirt and pants packed separately

Top 3 reasons to choose Medline

- 1 As a worldwide leader in disposable and reusable apparel, Medline provides stable supply and consistent product quality.
- 2 Medline helps you increase patient safety by producing apparel made of low-linting material.
- 3 Medline offers durable and comfortable apparel, designed in partnership with healthcare professionals.



Disposable scrub suit ranges packed as a set

Size	Soft range (blue)	Soft range (green)	Advanced range (blue)	Advanced range (green)	Advanced range (purple)	Packaging
XS	P35PBXS	P35PGXS	P42PBXS	P42PGXS	P42PPXS	1/bag, 25/case
S	P35PBS	P35PGS	P42PBS	P42PGS	P42PPS	1/bag, 25/case
M	P35PBM	P35PGM	P42PBM	P42PGM	P42PPM	1/bag, 25/case
L	P35PBL	P35PGL	P42PBL	P42PGL	P42PPL	1/bag, 25/case
XL	P35PBXL	P35PGXL	P42PBXL	P42PGXL	P42PPXL	1/bag, 25/case
2XL	P35PB2XL	P35PG2XL	P42PB2XL	P42PG2XL	P42PP2XL	1/bag, 25/case
3XL	P35PB3XL	P35PG3XL	P42PB3XL	P42PG3XL	P42PP3XL	1/bag, 25/case

Disposable scrub suit ranges packed separately

Size	Soft range (blue)	Advanced range (blue)	Advanced range (green)	Ultimate range (blue)	Ultimate range (green)	Packaging
	Pants	Pants	Pants	Pants	Pants	
XS	PXS9110P	P42PBXS-P	P42PGXS-P	P45PBXS-P	P45PGXS-P	1/bag, 50/case
S	PS9110P	P42PBS-P	P42PGS-P	P45PBS-P	P45PGS-P	1/bag, 50/case
M	PM9110P	P42PBM-P	P42PGM-P	P45PBM-P	P45PGM-P	1/bag, 50/case
L	PL9110P	P42PBL-P	P42PGL-P	P45PBL-P	P45PGL-P	1/bag, 50/case
XL	PXL9110P	P42PBXL-P	P42PGXL-P	P45PBXL-P	P45PGXL-P	1/bag, 50/case
2XL	P2XL9110P	P42PB2XL-P	P42PG2XL-P	P45PB2XL-P	P45PG2XL-P	1/bag, 50/case
3XL	-	P42PB3XL-P	P42PG3XL-P	P45PB3XL-P	P45PG3XL-P	1/bag, 50/case
	Shirt	Shirt	Shirt	Shirt	Shirt	
XS	PXS9110S	P42PBXS-S	P42PGXS-S	P45PBXS-S	P45PGXS-S	1/bag, 50/case
S	PS9110S	P42PBS-S	P42PGS-S	P45PBS-S	P45PGS-S	1/bag, 50/case
M	PM9110S	P42PBM-S	P42PGM-S	P45PBM-S	P45PGM-S	1/bag, 50/case
L	PL9110S	P42PBL-S	P42PGL-S	P45PBL-S	P45PGL-S	1/bag, 50/case
XL	PXL9110S	P42PBXL-S	P42PGXL-S	P45PBXL-S	P45PGXL-S	1/bag, 50/case
2XL	P2XL9110S	P42PB2XL-S	P42PG2XL-S	P45PB2XL-S	P45PG2XL-S	1/bag, 50/case
3XL	-	P42PB3XL-S	P42PG3XL-S	P45PB3XL-S	P45PG3XL-S	1/bag, 50/case

Warm-Up Jackets

1

Reduce costs

During the winter months, you may see an increase in the use of sterile surgical gowns because staff often wear these for extra warmth. By offering Warm-Up Jackets, you can help reduce costs while helping your staff comply with hygiene protocols.

2

Keep your staff warm and comfortable

Medline's single-use Warm-Up Jacket range is designed with materials that ensure a high level of comfort while keeping your staff warm. Additionally, they feature an intuitive design with a gentle knit collar, cuffs or elastic wrists, and handy pockets.



Soft long-sleeved jacket range



Essential long-sleeved jacket range



Essential three-quarter sleeve range



Essential short-sleeved range

Disposable Warm-Up Jackets

Essential long-sleeved jacket range			Essential three-quarter sleeve range			Essential short-sleeved range without cuffs		
Item number	Description	Packaging	Item number	Description	Packaging	Item number	Description	Packaging
NONERP600S	Blue, S	10/bag, 30/case	NONERP800S	Blue, S	10/bag, 30/case			
NONERP600M	Blue, M	10/bag, 30/case	NONERP800M	Blue, M	10/bag, 30/case	NONERP700M	Blue, M	10/bag, 30/case
NONERP600L	Blue, L	10/bag, 30/case	NONERP800L	Blue, L	10/bag, 30/case	NONERP700L	Blue, L	10/bag, 30/case
NONERP600XL	Blue, XL	10/bag, 30/case	NONERP800XL	Blue, XL	10/bag, 30/case	NONERP700XL	Blue, XL	10/bag, 30/case
Soft long-sleeved jacket range								
NONERP610S	Elastic at the waist, blue, S	10/bag, 50/case						
NONERP610M	Elastic at the waist, blue, M	10/bag, 50/case						
NONERP610L	Elastic at the waist, blue, L	10/bag, 50/case						
NONERP610XL	Elastic at the waist, blue, XL	10/bag, 50/case						

Headwear³

Covering your head and hair helps reduce the risk of foreign particles getting into surgical wounds and protects patients from potential contamination in operating theatres or ICUs. Medline's disposable headwear products offer a high level of comfort and come in a variety of shapes, sizes and designs that address hospitals' protective and aesthetical needs.

Our broad product range, including surgical hoods and caps, bonnets and bouffant caps, is suitable for a variety of operations and hairstyles. To offer you even more choice and meet your budget, Medline's headwear is available in cluster keyback, viscose and spunbond polypropylene materials that have passed biocompatibility requirements, are low linting and do not contain latex.

Comfort Line

The *Comfort Line* is Medline's preferred headwear product range for the operating theatre. Made of soft and breathable cluster keyback fabric, this line offers both a high level of comfort and a high level of protection. Comfort Line products are available in blue and green and a large variety of designs. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why our Comfort product range is produced in a WRAP-certified facility.



Surgical hood

FS61100

100/box, 600/case

Surgical hood, green



Surgical hood

FS61150A

100/box, 600/case

Surgical hood, blue



Surgical hood

FS61250

50/box, 300/case

Surgical hood with sweatband, blue



Bouffant cap

FS61300

100/box, 600/case

Bouffant cap with a latex-free elastic band, green, size 53 cm



Bouffant cap

FS61350

100/box, 600/case

Bouffant cap with a latex-free elastic band, blue, size 53 cm



Surgical cap

FS61450

100/box, 600/case

Surgical cap, sailor style, blue



Surgical cap

FS61500

100/box, 600/case

Surgical cap, elastic at the back, contoured front, green



Surgical cap

FS61550A

100/box, 600/case

Surgical cap, elastic at the back, contoured front, blue



Surgical cap

FS61750A

100/box, 600/case

Surgical cap with ties, for long hair, blue



Surgical cap

FS61950

100/box, 600/case

Surgical cap, elastic at the back, flat front, blue



Surgical cap

FS61850A

100/box, 600/case

Surgical cap with ties, for short hair, blue

Colourful Line

The *Colourful Line* is made of spunlace viscose, which is a soft, breathable and cotton-like material. The range comes in a variety of colours and offers good comfort. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why our Colourful product range is produced in a WRAP-certified facility.



Over-the-shoulders hood

FS60000

50/box, 300/case

Over-the-shoulders hood with ties and sweatband, blue



Surgical hood

FS60100G

100/box, 600/case

Surgical hood, green



Surgical hood

FS60150B

100/box, 600/case

Surgical hood, blue



Surgical hood

FS60250B

50/box, 300/case

Surgical hood with sweatband, blue



Surgical cap

FS60750B

100/box, 600/case

Surgical cap with ties, for long hair, blue



Surgical cap

FS60850B

100/box, 600/case

Surgical cap with ties, for short hair, blue



Surgical cap

FS60950B

100/box, 600/case

Surgical cap, elastic at the back, flat front, blue



Surgical cap

FS60950G

100/box, 600/case

Surgical cap, elastic at the back, flat front, green



Bouffant cap

FS63300

100/box, 600/case

Bouffant cap, size 53 cm, green



Bouffant cap

FS63350

100/box, 600/case

Bouffant cap, size 53 cm, blue



Surgical cap

FS60950P

100/box, 600/case

Surgical cap, elastic at the back, flat front, purple



Surgical cap

FS60950PK

100/box, 600/case

Surgical cap, elastic at the back, flat front, pink



Surgical cap

FS63450A

100/box, 600/case

Surgical cap, sailor style, blue



Clip cap

FS65450

100/box, 1,200/case

Clip cap, one size fits most, blue



Surgical cap

FS60950W

100/box, 600/case

Surgical cap, elastic at the back, flat front, white



Surgical cap

FS60950Y

100/box, 600/case

Surgical cap, elastic at the back, flat front, yellow



Clip cap

FS65460

100/box, 1,200/case

Clip cap, one size fits most, white



Bouffant cap

CRIE1000

100/box, 500/case

Bouffant cap, size 53 cm, green

Basic Line

Designed to satisfy the essential needs of the ward at a convenient price, the *Basic Line* is made of a spunbond material with a basis weight from 12 gsm to 30 gsm. Basic Line headwear comes in blue, green, white or printed patterns.



Bouffant cap

CRIE1002

100/box, 500/case

Bouffant cap, size 53 cm, white



Bouffant cap

CRIE1003

100/box, 500/case

Bouffant cap, size 60 cm, white



Surgical cap

NONE28625

100/box, 500/case

Surgical cap with ties, SMS top, absorbent scrim-reinforced material on sides, one size fits most, light blue



Surgical cap

FS67050

100/box, 600/case

Spunlace bonnet, one size fits most, blue



Bouffant cap

CRIE1004

100/box, 500/case

Bouffant cap, size 60 cm, blue



Surgical hood

FS63100A

100/box, 600/case

Surgical hood, green



Surgical hood

FS63150A

100/box, 600/case

Light surgical hood, blue

Choose your headwear



Surgical hoods

Offering a helmet-like cap with ties at the neck, one-size-fits-all surgical hoods provide the most coverage of all products in our headwear range for worry-free performance during surgery.

Some products feature an absorbent sweatband for increased staff comfort, especially for longer procedures.

Item number	Description	Packaging
FS61100	Cluster keyback, green	100/box, 600/case
FS61150A	Cluster keyback, blue	100/box, 600/case
FS61250	Cluster keyback, with sweatband, blue	50/box, 300/case
FS60100G	Viscose, green	100/box, 600/case
FS60150B	Viscose, blue	100/box, 600/case
FS60250B	Viscose, with sweatband, blue	50/box, 300/case
FS63100A	Spunbond polypropylene 30 gsm, green	100/box, 600/case
FS63150A	Spunbond polypropylene 30 gsm, blue	100/box, 600/case



Surgical caps with ties for long and short hair

Featuring ties at the back of the neck, these adjustable surgical caps are designed for people with long or short hair to safely keep the hair covered.

This surgical cap is ideal for use in the operating theatre or for day surgeries.

Item number	Description	Packaging
For long hair		
FS61750A	Cluster keyback, blue	100/box, 600/case
FS60750B	Viscose, blue	100/box, 600/case
For short hair		
FS61850A	Cluster keyback, blue	100/box, 600/case
FS60850B	Viscose, blue	100/box, 600/case
NONE28625	SMS and scrim-reinforced material, blue	100/box, 500/case



Surgical caps with elastic at the back

Featuring a broad shape and elastic at the back, these one-size-fits-most surgical caps provide a comfortable and secure fit, meeting the requirements of the most active of days.

This cap design is well-suited for long hairstyles and long periods of wear in the operating theatre, ICU and other laboratory departments.

Item number	Description	Packaging
FS61500	Cluster keyback, green	100/box, 600/case
FS61550A	Cluster keyback, blue	100/box, 600/case
FS61950	Cluster keyback, blue	100/box, 600/case
FS60950B	Viscose, blue	100/box, 600/case
FS60950G	Viscose, green	100/box, 600/case
FS60950P	Viscose, purple	100/box, 600/case
FS60950PK	Viscose, pink	100/box, 600/case
FS60950W	Viscose, white	100/box, 600/case
FS60950Y	Viscose, yellow	100/box, 600/case



Surgical caps without ties or elastic

These surgical caps do not contain ties or elastic, which makes them comfortable for longer periods of wear thanks to their lack of pressure areas. The caps are designed to fit well around the head due to the natural elasticity and flexibility of the material.

The caps are ideal for nurses, visitors in the operating theatre and ward staff with short hair. Additionally, the edge of the bonnet (FS67050) can be tuned up to function as a sweatband for even more comfort.

Item number	Description	Packaging
FS61450	Cluster keyback, blue	100/box, 600/case
FS63450A	Spunbond polypropylene 30 gsm, blue	100/box, 600/case
FS67050	Spunlace 35 gsm, blue	100/box, 600/case



Bouffant Caps

These light and roomy caps with a round shape feature elastic all the way around, ensuring all hair is fully and securely covered and that the cap fits comfortably.

Bouffant caps are suitable for long and short hairstyles and provide comfort and breathability for staff in less critical areas in the hospital.

Very similar to the design of the bouffant cap, the compactly packed clip cap is ideal headwear for patients and visitors in less critical areas of the hospital.

Item number	Description	Packaging
FS61300	Cluster keyback, 53 cm, green	100/box, 600/case
FS61350	Cluster keyback, 53 cm, blue	100/box, 600/case
FS63300	SPP 20 gsm, 53 cm, green	100/box, 600/case
FS63350	SPP 20 gsm, 53 cm, blue	100/box, 600/case
CRIE1000	SPP 12 gsm, 53 cm, green	100/box, 500/case
CRIE1002	SPP 12 gsm, 53 cm, white	100/box, 500/case
CRIE1003	SPP 12 gsm, 60 cm, white	100/box, 500/case
CRIE1004	SPP 12 gsm, 60 cm, blue	100/box, 500/case
FS65450	SPP 14 gsm, one size fits most, blue	100/box, 1,200/case
FS65460	SPP 14 gsm, one size fits most, white	100/box, 1,200/case

SPP = spunbond polypropylene

Patient gowns

Improving patient care is a priority for all healthcare providers. Most hospitals and care professionals take this goal a step further and have the overall objective to achieve a high level of patient-centred care and satisfaction.

At Medline, we think the quality of the product contributes to the quality of care that patients receive during their time in the hospital. Therefore, our patient gowns provide your patients with opacity, comfort and modesty during hospital stays because they are manufactured from a soft, non-woven SMS material.

Medline offers a wide variety of designs to meet the needs of your department. Choose from a selection of patient gowns with front openings, back openings or hook-and-loop fasteners on the sleeves.

Item number	Size	Description	Packaging
NONE27046 NONE27046XL	L XL	Sleeveless patient gown, blue	50/case
NONE27046SL NONE27046SLXL	L XL	Patient gown with short sleeves, blue	50/case
NONE27SLK	Suitable for children (4 to 9 years)	Patient gown for child with short sleeves, blue	100/case
CRIE4060K	L	Kimono gown with long sleeves, waist ties, blue	60/case
NONE27100 NONE27100XXL	L XXL	Patient gown with front opening and short sleeves, waist ties, blue	100/case



Patient apparel

Our apparel provides patients with comfort and dignity during their day surgery or diagnostic procedure because it is manufactured from a soft and opaque polypropylene material. Discover Medline's individual components and choose between colonoscopy shorts, underwear or elastic slippers to offer your patients the very best care and comfort.

Item number	Size	Description	Packaging
NONE27UP NONE27UPXXL	L XXL	Underwear with elastic, dark blue	200/case
NONECSHL NONECSH2XL	L XXL	Colonoscopy shorts with back opening and hook-and-loop closures, blue	100/case
CRIE2210	Standard size (32 cm length)	Elastic slippers with rigid sole, white	100/case
MDTEFP218RXL	XL	Socks, red	48/case



Patient accessories

Medline offers convenient accessories to make patient flow and stay in the ambulatory care area more effortless.

Item number	Size	Description	Packaging
PBAG5060	50 x 60 cm	Patient belongings bag, white	250/case
NWP150	106 x 150 cm	Non-woven panel	200/case
TLPC1527	15 x 27 cm	Telemetry pouch for monitoring devices, blue	200/case




Patient kits⁴

Our expertise in medical apparel design and manufacturing allows us to offer new solutions for day surgeries. With efficiency in mind, Medline has created a variety of patient kit solutions and a visitor kit to suit your needs and budget.

The kits contain single-use patient apparel, which patients can wear during the diagnostic procedure or day surgery. All components are conveniently packed in a single bag that can be given to the patient upon arrival. Due to its packaging, organising storage and replenishing stock are easy to handle, and you only have to order one item instead of four or five individual components.

As a leading international manufacturer and distributor of medical devices, we understand your focus and your patients' needs. We apply the highest quality standards to our products along each step of the manufacturing and assembly process.

Item number	Size	Description	Packaging
KIT1030	L	Kimono, blue Patient gown with short sleeves, blue Underwear, dark blue Clip cap, blue Socks, red	40/case
KIT1035 KIT1035XXL	L XXL	Kimono, purple Patient gown with short sleeves, purple Underwear, dark blue Clip cap, blue Socks, red	40/case
KIT3020 KIT3020XXL	L XXL	Patient gown with short sleeves, blue Underwear, dark blue Clip cap, white Shoe covers, white	50/case
NONE27KITK	Suitable for children (4 to 9 years)	Patient gown for child with short sleeves, blue Underwear, dark blue Clip cap, white Shoe covers, white	50/case
NONE35LKIT NONE352XLKIT	L XXL	Patient gown with short sleeves, blue Scrub suit pants, blue Clip cap, white Shoe covers, textured bottom, white and blue	50/case
CRIE400BKIT	L	Cover gown with elastic wrists, blue Facemask Type II with ear loops, blue Clip cap, white Shoe covers, white	50/case
KIT1050E	L	Kimono, blue Patient gown with short sleeves, blue Scrub suit pants, blue Underwear, dark blue Shoe covers, textured bottom, white and blue Socks, red	25/case
KIT1440E KIT14402XE	L XXL	Patient gown with short sleeves, front opening, blue Scrub suit pants, blue Clip cap, blue Shoe covers, textured bottom, white and blue Socks, red	30/case



KIT1030

Medline also offers a large range of surgical, examination and respiratory masks.

Check out the facemask brochure for more details.

Facemasks

Surgical, procedural and respirators



Masks at the ready

For more information, contact your Medline account manager or visit us at: www.medline.eu/uk.



Medline Industries Ltd

3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237

www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd

1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

¹⁾ Footwear products are class I medical devices according to Regulation (EU) 2017/745.

²⁾ Disposable scrub suits are class I medical devices according to Regulation (EU) 2017/745.

³⁾ Headwear products are class I medical devices according to Regulation (EU) 2017/745.

⁴⁾ Patient kits are procedure packs according to Article 22 of Regulation (EU) 2017/745.

The procedure pack producer is Medline International France SAS, 5 Rue Charles Lindbergh, 44110 Châteaubriant, France.

These products are intended to be used by healthcare professionals.

Before use, please consult instructions and precautions on the corresponding labelling.



Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France



All rights reserved. We reserve the right to correct errors that may have occurred.

© 2022 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093, USA. ML179_EN 01/2022