

# Personal Protection

We got you covered!




**Reduce cross-contamination and protect your staff.**

# Why are more healthcare professionals choosing Medline for their all-round protection needs?



**High-quality products**



**Comprehensive product range**



**Reliable partner**

Your staff are exposed to bacteria and bodily fluids on a daily basis, putting them at risk of cross-contamination. We know you rely on disposable apparel to keep them protected, but also comfortable. Therefore, we apply the highest quality standards to our products and in our manufacturing facilities to ensure that healthcare professionals get the protection they deserve.

Medline International France is EN ISO 13485 and ISO 9001 certified and complies with regulatory standards on the market. In addition, we are conscious about social responsibility and, as a result, some of our products are WRAP-certified. WRAP (Worldwide Responsible Accredited Production) is an internationally recognised social compliance certification body that verifies whether factories respect workers' labour and human rights.

## Table of contents

European regulations and directives .....	3
Guide to protective materials .....	3
Isolation gowns overview .....	4-5
Cover gowns .....	6
Isolation gowns .....	6-7
Impervious cover gowns .....	7-8
Medium-weight bilaminate impervious isolation gowns .....	8
Chemo protection isolation gowns .....	9
Footwear .....	10
Disposable scrub suits .....	11-14
Warm-up jackets .....	15
Headwear .....	16-21
Patient gowns .....	22
Patient apparel .....	22
Patient accessories .....	22
Patient kits .....	23



It is critical for our customers to receive protective apparel on time. That is why we have selected certain items where we have built a large safety stock in order to react faster to any unexpected increase in your demand. **These items are marked in bold for easy identification.**

# European regulations and directives to classify protective apparel

## Regulation (EU) 2016/425 on personal protective equipment (PPE)

The personal protective equipment included in Medline's range provides protection against minimal risks (contact with weak-action cleaning materials or prolonged contact with water) or more serious risks, such as exposure to substances and mixtures that are hazardous to health. Products listed in the brochure that are certified as PPE are compliant with PPE Category I or Category III requirements.

## Regulation (EU) 2017/745 on medical devices (MDR)

Medical devices include any instrument, apparatus, appliance, software, implant, reagent, material or other article intended by the manufacturer to be used, alone or in combination, for human beings for a specific medical purpose. The medical devices included in Medline's range are intended to protect the patient from potential contamination. They are Class I devices.

## Directive 2001/95/EC on general product safety

Products included in Medline's protective apparel range that do not fall under the definition of PPE or MDR are certified under the general product safety directive. These are general commodity products (e.g. underwear).

*Please consult the product's technical data sheet to see its classification.*

## Guide to protective materials

To help you choose the right product, here are brief descriptions of the features and benefits of the main materials available from Medline. You can refer back to this page as you look through the brochure.

### Spunbond polypropylene (SPP)

This standard fabric is made by bonding fibres together to form a single layer of breathable woven-like material. It is suitable for protecting users against mechanical particles and hair.

### Multi-ply polypropylene (SMS)

This fabric consists of inner layers of meltblown polypropylene sandwiched between outer layers of spunbond polypropylene. Strong, fluid-resistant and breathable, this multilayer material provides an optimal mix of fluid protection and comfort.

### Polyethylene film

This single layer of polyethylene film provides high fluid-resistant properties at a very attractive price.








### Bilaminate polypropylene-polyethylene

This fabric is created by laminating a layer of polypropylene (SPP or SMS) with a layer of impervious polyethylene. This fabric is ideal for any situation demanding a material with high tensile strength and fluid protection.

### Trilaminate polypropylene-polyethylene

This impervious and breathable material consists of a hybrid film laminated between two layers of a bicomponent spunbond material. This fabric is ideal for long interventions when impervious protection is needed against fluids and chemotherapy drugs.

# Protective gowns overview

Product family	Code	Size	PF	Colour	Material	Length (cm)	Width (cm)	Wrist type		
Cover gowns	CRIE4000	L	50		Polypropylene (SPP)	114	135	Elastic		
	CRIE4000G									
	<b>CRIE4000B</b>									
	CRIE4001	XL				127	147			
	CRIE4001G									
	<b>CRIE4001B</b>									
	CRIE4001BCF									
Isolation gowns	NONE27SMS3	L	100		Multi-ply polypropylene (SMS)	114	142	Elastic		
	NONE27SMS3XL	XL				127	152			
	NONELV200C	L				114	142	Knit cuffs		
	NONELV200CXL	XL				127	152			
	NONE326MCF	M				106	132	Knit cuffs		
	<b>NONE326CF</b>	L				118	142			
	<b>NONE326XLCF</b>	XL				135	152			
	NONE27280	L				25		114	142	Knit cuffs
	NONE27280XL	XL								
Impervious gowns	<b>CRIE4040C*</b>	L	100		Bilaminate polypropylene-polyethylene	114	135	Knit cuffs		
	<b>CRIE4041C*</b>	XL				135	152			
	CRIE4042C*	2XL				147	162			
	CRIE4040BC*	L	100		114	135				
	CRIE4041BC*	XL			135	152				
	NONELV350	L	50		Microporous laminate	114	142	Knit cuffs		
	<b>CRIE5010</b>	L	75		Polyethylene	114	94	Thumb loops		
	<b>CRIE5011</b>	XL				134	113			
	NONE27116	L	50		Bilaminate polypropylene-polyethylene	114	142	Knit cuffs		
	NONE27116XL	XL				127	142			
	NONE272790HXL	XL				127	86			
	760510-N	L	100		130	140				
	760520-N	XL				130	160			
	NONE27457	L	30		Trilaminate polypropylene-polyethylene	114	142			
NONE27457XL	XL	127				152				
One-piece coverall with hood	CN4013DTM	M	25		Microporous laminate	190	61	Knit cuffs		
	CN4013DTL	L				200	66			
	CN4013DTXL	XL				207	71			
	CN4013DT2XL	2XL				212	76			

\*Isolation gown with back panel made of a breathable material

Neck closure	ISO811: Resistance to liquid penetration (cmH2O)	Product classification	Intended use	Sterile/ non-sterile
Ties	n/a	MDR Class I	Protects the patient against cross-contamination	Non-sterile
Ties	10	PPE Category I	PPE Category I: Protects the user against prolonged contact with water	Non-sterile
Ties				
Hook and loop	50.8			
Ties	54.5			Sterile
Ties	107	PPE Category I	Protects the user against prolonged contact with water and cleansing agents  [CRIE4040xC range and CRIE500x are tested according to EN 16523-1 against the cleansing agent glucoprotamin – breakthrough time of >480 minutes: performance level Class 6 according to EN 14325]	Non-sterile
Ties	181			
Over-the-head	75.5			
Ties	114	PPE Category I	Protects the user against prolonged contact with water	
Over-the-head				
Ties	-	PPE Category III and MDR Class I	Protects the user against infective pathogens (bacteria and viruses), also MRSA-MRGN and noroviruses	
	-			
Hook and loop	200	PPE Category III	Protects the user against substances and mixtures that are hazardous to health (e.g. chemotherapy drugs)	
Over-the-head hood	-	PPE Category III	Protects against against hazardous substances and infectious agents (bacteria's and viruses)	Non-sterile
	-			
	-			
	-			

# Cover gowns<sup>1</sup>

Suitable for preventing particulate contamination

## Lightweight polypropylene cover gowns

- Breathable and comfortable
- Waist and neck ties
- Available with elastic wrists or knit cuffs

Item number	Size	Colour	Packaging
<b>Cover gowns with elastic wrists</b>			
CRIE4000	L	Yellow	10/bag, 50/case
CRIE4000G	L	Green	10/bag, 50/case
<b>CRIE4000B</b>	L	Blue	10/bag, 50/case
CRIE4001	XL	Yellow	10/bag, 50/case
CRIE4001G	XL	Green	10/bag, 50/case
<b>CRIE4001B</b>	XL	Blue	10/bag, 50/case



CRIE4001G

## Isolation gowns

Suitable for protection against minimal risks from prolonged contact with water

### Lightweight SMS isolation gowns

- Breathable and comfortable
- Waist and neck ties
- Available with elastic wrists or knit cuffs

Item number	Size	Colour	Packaging
<b>Neck and waist ties, elastic wrists</b>			
NONE27SMS3	L	Blue	10/bag, 50/case
NONE27SMS3XL	XL	Blue	10/bag, 50/case
<b>Neck and waist ties, knit cuffs</b>			
NONELV200C	L	Yellow	10/bag, 100/case
NONELV200CXL	XL	Yellow	10/bag, 100/case



NONE27SMS3



## Medium-weight SMS isolation gowns

The material used in this range of SMS isolation gowns has a hydrostatic head performance higher than 50 cmH<sub>2</sub>O according to ISO 811.

Item number	Size	Colour	Sterile/ non-sterile	Packaging
Hook and loop at the neck, waist ties, knit cuffs				
NONE326MCF	M	Yellow	Non-sterile	10/bag, 100/case
<b>NONE326CF</b>	L	Yellow	Non-sterile	10/bag, 100/case
<b>NONE326XLCF</b>	XL	Yellow	Non-sterile	10/bag, 100/case



NONE326CF

## Impervious cover gowns

Suitable for protecting the user against prolonged contact with water and cleansing agents; CRIE404xC range and CRIE501x range tested according to EN 16523-1 against the cleansing agent glucoprotamin

### Lightweight bilaminate impervious cover gowns

- Waist and neck ties
- Knit cuffs

Item number	Size	Colour	Sterile/ non-sterile	Packaging
<b>CRIE4040C</b>	L	Yellow	Non-sterile	10/bag, 100/case
<b>CRIE4041C</b>	XL	Yellow	Non-sterile	10/bag, 100/case
CRIE4042C	2XL	Yellow	Non-sterile	10/bag, 100/case
CRIE4040BC	L	Blue	Non-sterile	10/bag, 100/case
CRIE4041BC	XL	Blue	Non-sterile	10/bag, 100/case



CRIE4040C

### Polyethylene impervious cover gowns

- Over-the-head design
- Thumb-loop-style wrists
- Individually folded and wrapped in a dispenser box

Item number	Size	Colour	Packaging
<b>CRIE5010</b>	L	Blue	1/bag, 15/box, 75/case
<b>CRIE5011</b>	XL	Blue	1/bag, 15/box, 75/case



CRIE5010

## Polyethylene disposable apron

Worn over non-fluid-resistant gowns or scrub suits as a convenient solution for users that want to add extra fluid protection at a lower cost.

Item number	Description	Packaging
762050	Medium weight, 125 x 80 cm, individually wrapped, dispenser box	100/box



762050

## Microporous bilaminate cover gown

- Waist and neck ties
- Knit cuffs

Item number	Size	Colour	Packaging
NONELV350	L	White	10/bag, 50/case



NONELV350

## Medium-weight bilaminate impervious isolation gowns

Suitable for protection against prolonged contact with water

- Waist and neck ties
- Knit cuffs

Item number	Size	Colour	Packaging
NONE27116	L	Blue	10/bag, 50/case
NONE27116XL	XL	Blue	10/bag, 50/case
NONE272790HXL*	XL	Blue	10/bag, 50/case

\*Over-the-head design and waist ties



NONE27116



# Isolation gown for protection against infective agents and chemotherapy drugs

Trilaminate polypropylene-polyethylene (Prevention Plus) gowns are certified as PPE Category III Type PB-[4] for partial body protection according to PPE Regulation 2016: 425 and are suitable for protection against mixtures that are hazardous to health. It is also certified as EN14126 for protection against infective agents; bacteria and viruses. Each gown is made with Prevention Plus material which makes it very comfortable to wear due to its breathability properties.

- Hook-and-loop closure at the neck
- Waist ties
- Knit cuffs at the wrists

## Tested against 14 chemotherapy drugs per ASTM F739

Carmustine (BiCNU)	Etoposide (Toposar)
Dacarbazine	Vincristine Sulphate
Fluorouracil (5FU)	Doxorubicin hydrochloride
Thiotepa	Paclitaxel (Taxol)
Cisplatin	Leucovorin
Mitomycin C	Docetaxel
Cyclophosphamide (Cytosan)	Oxaliplatin

*Test data is available on request.*



NONE27457: gown  
NONE27348PXL: boot covers

Item number	Size	Colour	Packaging
NONE27457	L	Blue	10/bag, 30/case
NONE27457XL	XL	Blue	10/bag, 30/case
760510-N	L	Yellow	10/bag, 10/case
760520-N	XL	Yellow	10/bag, 10/case

**Complete protection is available in combination with the trilaminate polypropylene-polyethylene boot covers (NONE27348PXL).**

# Footwear<sup>2</sup>

## Shoe covers

Footwear covers serve as a vital barrier, providing an effective solution to address environmental contamination and the transmission of microorganisms. Manufactured from polyethylene or polypropylene materials, some models feature additional fluid protection.

Item number	Size	Material	Characteristics	Packaging
CRIE2002	One size fits most	Polypropylene (SPP)	Protection	100/box, 300/case
<b>CRIE2010</b>	One size fits most	Polyethylene film	Impervious	100/bag, 1,000/case
CRIE2102W	One size fits most	Polypropylene (SPP)	Protection	500/box, 1,000/case
<b>CRIE2020</b>	One size fits most	Polyethylene film and polypropylene (SPP)	Impervious	100/bag, 200/case
NONE28858	L	Multi-ply polypropylene (SMS)	Fluid-resistant	100/box, 300/case



CRIE2002



CRIE2010



CRIE2102W



CRIE2020



NONE28858

## Boot covers

Medline's boot cover range offers enhanced coverage and increased fluid protection. The products cover the shoe and leg all the way up to the knee and have a convenient one-size-fits-most design with a textured sole on the bottom. They also feature a hook-and-loop strap or elastic bands around the ankle to help keep them in place and ensure adequate protection.

Item number	Size	Characteristics	Packaging
NONE27143XL	One-size-fits-most	Fluid-resistant	50/box, 150/case
NONE27144	One-size-fits-most	Impervious	50/box, 150/case
NONE27348PXL	One-size-fits-most	Breathable and impervious	50/box, 150/case

### NONE27143XL

- Breathable, fluid-resistant material above the ankle
- High fluid protection on the bottom
- Sealed seams on the foot



### NONE27144

- Impervious coated polyethylene-polypropylene material
- High fluid protection on the bottom
- Hook-and-loop strap for an adjustable fit



### NONE27348PXL

- Made of blue Prevention Plus material
- Hook-and-loop closure
- Foam strip



# Disposable scrub suits<sup>3</sup>

## A scrub suit designed with you in mind

We have worked with healthcare professionals like you to create our scrub suits with features you want and need to perform at your best. Experience new, clean apparel every day.



Bell-shaped design for a perfect fit



3 convenient pockets on the shirt



V-necked




Adjustable waist ties



With a side pocket



Shirt and pants packed separately



Available  
in **7** sizes

## Carefully selected materials

Medline's disposable scrub suits offer a high level of comfort thanks to their ultra-soft and breathable SMS material with the following features:



Great opacity



Low-linting



Anti-static



Latex-free

## Soft Scrub Suits range

Essential apparel to meet daily challenges

- Manufactured from 35-gsm SMS material with optimal opacity and softness
- Available in 7 sizes from XS to 3XL
- Available in **blue**, **green** and **purple** with shirt and pants packed separately or in sets





## Ultimate Scrub Suits range

### Our top-of-the-line ultra-soft range for demanding activities

Medline's ultimate scrub suits offer that extra comfort, breathability and protection you and your patients deserve.

- Manufactured from 45-gsm SMS material with excellent opacity and softness
- Available in 7 sizes from XS to 3XL
- Available in **blue** and **green** with shirts and pants packed separately or in sets

## Top 3 reasons to choose Medline

1

As a worldwide leader in disposable and reusable apparel, Medline provides stable supply and consistent product quality.

2

Medline helps you increase patient safety by producing apparel made of low-linting material.

3

Medline offers durable and comfortable apparel, designed in partnership with healthcare professionals.



## Disposable scrub suit ranges packed as a set

Size	Soft range (blue)	Soft range (green)	Soft range (purple)	Ultimate range (blue)	Ultimate range (green)	Packaging
XS	<b>P35PBXS</b>	P35PGXS	P35PPXS	P45PBXS	P45PGXS	1/bag, 25/case
S	<b>P35PBS</b>	P35PGS	P35PPS	P45PBS	P45PGS	1/bag, 25/case
M	<b>P35PBM</b>	P35PGM	P35PPM	P45PBM	P45PGM	1/bag, 25/case
L	<b>P35PBL</b>	P35PGL	P35PPL	P45PBL	P45PGL	1/bag, 25/case
XL	<b>P35PBXL</b>	P35PGXL	P35PPXL	P45PBXL	P45PGXL	1/bag, 25/case
2XL	<b>P35PB2XL</b>	P35PG2XL	P35PP2XL	P45PB2XL	P45PG2XL	1/bag, 25/case
3XL	<b>P35PB3XL</b>	P35PG3XL	P35PP3XL	P45PB3XL	P45PG3XL	1/bag, 25/case

## Disposable scrub suit ranges packed separately

Size	Soft range (blue)	Soft range (green)	Ultimate range (blue)	Packaging
	<b>Pants</b>	<b>Pants</b>	<b>Pants</b>	
XS	PXS9110P	P35PGXS-P	P45PBXS-P	1/bag, 50/case
S	PS9110P	P35PGS-P	P45PBS-P	1/bag, 50/case
M	PM9110P	P35PGM-P	P45PBM-P	1/bag, 50/case
L	PL9110P	P35PGL-P	P45PBL-P	1/bag, 50/case
XL	PXL9110P	P35PGXL-P	P45PBXL-P	1/bag, 50/case
2XL	P2XL9110P	P35PG2XL-P	P45PB2XL-P	1/bag, 50/case
3XL	-	P35PG3XL-P	P45PB3XL-P	1/bag, 50/case
	<b>Shirt</b>	<b>Shirt</b>	<b>Shirt</b>	
XS	PXS9110S	P35PGXS-S	P45PBXS-S	1/bag, 50/case
S	PS9110S	P35PGS-S	P45PBS-S	1/bag, 50/case
M	PM9110S	P35PGM-S	P45PBM-S	1/bag, 50/case
L	PL9110S	P35PGL-S	P45PBL-S	1/bag, 50/case
XL	PXL9110S	P35PGXL-S	P45PBXL-S	1/bag, 50/case
2XL	P2XL9110S	P35PG2XL-S	P45PB2XL-S	1/bag, 50/case
3XL	-	P35PG3XL-S	P45PB3XL-S	1/bag, 50/case



# Warm-Up Jackets



## Reduce costs

During the winter months, you may see an increase in the use of sterile surgical gowns because staff often wear these for extra warmth. By offering warm-up jackets, you can help reduce costs while helping your staff comply with hygiene protocols.



## Keep your staff warm and comfortable

Medline's single-use warm-up jacket range is designed with materials that ensure a high level of comfort while keeping your staff warm. Additionally, they feature an intuitive design with a gentle knit collar, cuffs or elastic wrists, and handy pockets.



Essential long-sleeved warm-up jacket range



Soft long-sleeved warm-up jacket range



Ultimate long-sleeved warm-up jacket range



Essential three-quarter sleeve warm-up jacket range



Essential short-sleeved warm-up jacket range



Soft long-sleeved warm-up jacket

## Disposable Warm-Up Jackets

Essential long-sleeved jacket range			Essential three-quarter sleeve range			Essential short-sleeved range without cuffs		
Item number	Description	Packaging	Item number	Description	Packaging	Item number	Description	Packaging
NONERP600S	Blue, S	10/bag, 30/case	NONERP800S	Blue, S	10/bag, 30/case			
NONERP600M	Blue, M	10/bag, 30/case	NONERP800M	Blue, M	10/bag, 30/case	NONERP700M	Blue, M	10/bag, 30/case
NONERP600L	Blue, L	10/bag, 30/case	NONERP800L	Blue, L	10/bag, 30/case	NONERP700L	Blue, L	10/bag, 30/case
NONERP600XL	Blue, XL	10/bag, 30/case	NONERP800XL	Blue, XL	10/bag, 30/case	NONERP700XL	Blue, XL	10/bag, 30/case
Soft long-sleeved jacket range			Soft long-sleeved warm-up jacket					
			NONERP601S	Purple, S	10/bag, 30/case			
NONERP610M	Elastic at the waist, blue, M	10/bag, 50/case	NONERP601M	Purple, M	10/bag, 30/case			
NONERP610L	Elastic at the waist, blue, L	10/bag, 50/case	NONERP601L	Purple, L	10/bag, 30/case			
NONERP610XL	Elastic at the waist, blue, XL	10/bag, 50/case	NONERP601XL	Purple, XL	10/bag, 30/case			

# Headwear<sup>4</sup>

Covering your head and hair helps reduce the risk of foreign particles getting into surgical wounds and protects patients from potential contamination in operating theatres or ICUs. Medline's disposable headwear products offer a high level of comfort and come in a variety of shapes, sizes and designs that address hospitals' protective and aesthetic needs.

Our broad product range, including surgical hoods and caps, bonnets and bouffant caps, is suitable for a variety of operations and hairstyles. To offer you even more choice and meet your budget, Medline's headwear is available in cluster keyback, viscose and spunbond polypropylene materials that have passed biocompatibility requirements, are low linting and do not contain latex.

## Comfort Line

The *Comfort Line* is Medline's preferred headwear product range for the operating theatre. Made of soft and breathable cluster keyback fabric, this line offers both a high level of comfort and a high level of protection. Comfort Line products are available in blue and green and a large variety of designs. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why our Comfort product range is produced in a WRAP-certified facility.



### Surgical hood

FS61100

100/box, 600/case

Surgical hood, green



### Surgical hood

FS61150A

100/box, 600/case

Surgical hood, blue



### Surgical hood

FS61250

50/box, 300/case

Surgical hood with sweatband, blue



### Bouffant cap

FS61300

100/box, 600/case

Bouffant cap with a latex-free elastic band, green, size 53 cm



### Bouffant cap

FS61350

100/box, 600/case

Bouffant cap with a latex-free elastic band, blue, size 53 cm



### Surgical cap

FS61450

100/box, 600/case

Surgical cap, sailor style, blue



### Surgical cap

FS61500

100/box, 600/case

Surgical cap, elastic at the back, contoured front, green



### Surgical cap

FS61550A

100/box, 600/case

Surgical cap, elastic at the back, contoured front, blue



**Surgical cap**

**FS61750A**

100/box, 600/case

Surgical cap with ties, for long hair, blue



**Surgical cap**

**FS61850A**

100/box, 600/case

Surgical cap with ties, for short hair, blue



**Surgical cap**


**FS61950**

100/box, 600/case

Surgical cap, elastic at the back, flat front, blue

**Colourful Line**

The *Colourful Line* is made of spunlace viscose, which is a soft, breathable and textile-like material. The range comes in a variety of colours and offers good comfort. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why most of our Colourful range is produced in a WRAP-certified facility.



**Over-the-shoulders hood**

**FS60000**

50/box, 300/case

Over-the-shoulders hood with ties and sweatband, blue




**Surgical hood**

**FS60100G**

100/box, 600/case

Surgical hood, green



**Surgical hood**

**FS60150B**

100/box, 600/case

Surgical hood, blue



**Surgical hood**

**FS60250B**

50/box, 300/case

Surgical hood with sweatband, blue



**Surgical cap**

**FS60750B**

100/box, 600/case

Surgical cap with ties, for long hair, blue



**Surgical cap**

**FS60850B**

100/box, 600/case

Surgical cap with ties, for short hair, blue



**Surgical cap**

**FS60950B**

100/box, 600/case

Surgical cap, elastic at the back, flat front, blue



**Surgical cap**

**FS60950G**

100/box, 600/case

Surgical cap, elastic at the back, flat front, green



**Surgical cap**

**FS60950P**

100/box, 600/case

Surgical cap, elastic at the back, flat front, purple



**Surgical cap**

**FS60950PK**

100/box, 600/case

Surgical cap, elastic at the back, flat front, pink



**Surgical cap**

**FS60950W**

100/box, 600/case

Surgical cap, elastic at the back, flat front, white



**Surgical cap**

**FS60950Y**

100/box, 600/case

Surgical cap, elastic at the back, flat front, yellow

## Basic Line

Designed to satisfy the essential needs of the ward at a convenient price, the *Basic Line* is made of a spunbond material with a basis weight from 12 gsm to 30 gsm. Basic Line headwear comes in blue, green, white or printed patterns.



**Surgical cap**

**NONE28625**

100/box, 500/case

Surgical cap with ties, SMS top, absorbent scrim-reinforced material on sides, one size fits most, light blue



**Surgical cap**

**FS67050**

100/box, 600/case

Spunlace bonnet, one size fits most, blue



**Surgical hood**

**FS63100A**

100/box, 600/case

Surgical hood, green



**Surgical hood**

**FS63150A**

100/box, 600/case

Light surgical hood, blue



**Bouffant cap**

**FS63300**

100/box, 600/case

Bouffant cap, size 60 cm, green



**Bouffant cap**

**FS63350**

100/box, 600/case

Bouffant cap, size 60 cm, blue



**Bouffant cap**

**FS63360**

100/box, 600/case

Bouffant cap, size 60 cm, white



**Surgical cap**

**FS63450A**

100/box, 600/case

Surgical cap, sailor style, blue



**Clip cap**

**FS65460**

100/box, 1,200/case

Clip cap, one size fits most, white



**Clip cap**

**FS65450**

100/box, 1,200/case

Clip cap, one size fits most, blue



**Bouffant cap**

**CRIE1000**

100/box, 500/case

Bouffant cap, size 53 cm, green



**Bouffant cap**

**CRIE1002**

100/box, 500/case

Bouffant cap, size 53 cm, white



**Bouffant cap**

**CRIE1004**


100/box, 500/case

Bouffant cap, size 53 cm, blue



# Headwear accessory

## Individual sweatband



Medline's individual sweatband can be placed under any surgical cap and provides comfort to healthcare professionals during long procedures thanks to its soft viscose material.

Item number	Description	Packaging
BAND250	Sweatband with elastic, 18 x 4 cm	50/box, 400/case

# Choose your headwear

## Surgical hoods



Offering a helmet-like cap with ties at the neck, one-size-fits-all surgical hoods provide the most coverage of all products in our headwear range for worry-free performance during surgery.

Some products feature an absorbent sweatband for increased staff comfort, especially for longer procedures.

Item number	Description	Packaging
FS61100	Cluster keyback, green	100/box, 600/case
FS61150A	Cluster keyback, blue	100/box, 600/case
FS61250	Cluster keyback, with sweatband, blue	50/box, 300/case
FS60100G	Viscose, green	100/box, 600/case
FS60150B	Viscose, blue	100/box, 600/case
FS60250B	Viscose, with sweatband, blue	50/box, 300/case
FS63100A	Spunbond polypropylene 30 gsm, green	100/box, 600/case
FS63150A	Spunbond polypropylene 30 gsm, blue	100/box, 600/case

## Surgical caps with ties for long and short hair



Featuring ties at the back of the neck, these adjustable surgical caps are designed for people with long or short hair to safely keep the hair covered.

This surgical cap is ideal for use in the operating theatre or for day surgeries.

Item number	Description	Packaging
<b>For long hair</b>		
FS61750A	Cluster keyback, blue	100/box, 600/case
FS60750B	Viscose, blue	100/box, 600/case
<b>For short hair</b>		
FS61850A	Cluster keyback, blue	100/box, 600/case
FS60850B	Viscose, blue	100/box, 600/case
NONE28625	SMS and scrim-reinforced material, blue	100/box, 500/case



## Surgical caps with elastic at the back



Featuring a broad shape and elastic at the back, these one-size-fits-most surgical caps provide a comfortable and complete hair coverage, meeting the requirements of the most active of days.

This cap design is well-suited for long hairstyles and long periods of wear in the operating theatre, ICU and other laboratory departments.

Item number	Description	Packaging
FS61500	Cluster keyback, green	100/box, 600/case
FS61550A	Cluster keyback, blue	100/box, 600/case
FS61950	Cluster keyback, blue	100/box, 600/case
FS60950B	Viscose, blue	100/box, 600/case
FS60950G	Viscose, green	100/box, 600/case
FS60950P	Viscose, purple	100/box, 600/case
FS60950PK	Viscose, pink	100/box, 600/case
FS60950W	Viscose, white	100/box, 600/case
FS60950Y	Viscose, yellow	100/box, 600/case

## Surgical caps without ties or elastic



These surgical caps do not contain ties or elastic, which makes them comfortable for longer periods of wear thanks to their lack of pressure areas. The caps are designed to fit well around the head due to the natural flexibility of the material.

The caps are ideal for nurses, visitors in the operating theatre and ward staff with short hair. Additionally, the edge of the bonnet (FS67050) can be tuned up to function as a sweatband for even more comfort.

Item number	Description	Packaging
FS61450	Cluster keyback, blue	100/box, 600/case
FS63450A	Spunbond polypropylene 30 gsm, blue	100/box, 600/case
FS67050	Spunlace 35 gsm, blue	100/box, 600/case

## Bouffant caps



These light and roomy caps with a round shape feature elastic all the way around, ensuring all hair is fully and securely covered and that the cap fits comfortably.

Bouffant caps are suitable for long and short hairstyles and provide comfort and breathability for staff in less critical areas in the hospital.

Very similar to the design of the bouffant cap, the compactly packed clip cap is ideal headwear for patients and visitors in less critical areas of the hospital.

Item number	Description	Packaging
FS61300	Cluster keyback, 53 cm, green	100/box, 600/case
FS61350	Cluster keyback, 53 cm, blue	100/box, 600/case
FS63300	SPP 20 gsm, 60 cm, green	100/box, 600/case
FS63350	SPP 20 gsm, 60 cm, blue	100/box, 600/case
FS63360	SPP 20 gsm, 60 cm, white	100/box, 600/case
<b>CRIE1000</b>	SPP 12 gsm, 53 cm, green	100/box, 500/case
CRIE1002	SPP 12 gsm, 53 cm, white	100/box, 500/case
<b>CRIE1004</b>	SPP 12 gsm, 53 cm, blue	100/box, 500/case
<b>FS65450</b>	SPP 14 gsm, one size fits most, blue	100/box, 1,200/case
FS65460	SPP 14 gsm, one size fits most, white	100/box, 1,200/case

SPP = spunbond polypropylene

# Patient apparel

Medline offers a wide variety of designs and sizes to meet the needs of your department. Choose from a selection of patient gowns with front or back openings. Therefore, our patient gowns provide your patients with opacity, comfort and modesty during hospital stays because they are manufactured from a soft, anti-static and low-linting non-woven SMS material.

Item number	Size	Description	Packaging
<b>NONE27046</b> <b>NONE27046XL</b>	L XL	Sleeveless patient gown, blue	50/case
NONE27046SL NONE27046SLXL NONE27046SL2XL	L XL 2XL	Patient gown with short sleeves, blue	50/case
NONE275LK	Suitable for children (3 to 5 years)	Patient gown for children with short sleeves, blue	100/case
CRIE4060K CRIE4062K	L 2XL	Kimono gown with long sleeves, waist ties, blue	60/case
NONE27100M NONE27100 NONE27100XXL	M L 2XL	Patient gown with front opening and short sleeves, waist ties, blue	100/case



NONE27046



NONE27046SL



CRIE4060K

Our apparel provides patients with comfort and dignity during their day surgery or diagnostic procedure because it is manufactured from a soft and opaque polypropylene material. Discover Medline's individual components and choose between colonoscopy shorts, underwear or elastic slippers to offer your patients the very best care and comfort.

Item number	Size	Description	Packaging
NONE27UP NONE27UPXXL	L 2XL	Underwear with elastic, dark blue	200/case
NONECSHL NONECSH2XL	L 2XL	Colonoscopy shorts with back opening and hook-and-loop closures, blue	100/case
CRIE2205 CRIE2210	27 x 10,5 cm 32 x 11 cm	Elastic slippers with rigid sole, white	100/case



NONE27UP



CRIE2210

# Patient accessories

Medline offers convenient accessories to make patient flow and stay in the ambulatory care area more effortless.


Item number	Size	Description	Packaging
PBAG5060	50 x 60 cm	Patient belongings bag, white	250/case
NWP150	106 x 150 cm	Non-woven panel	200/case
TLPC1527	15 x 27 cm	Telemetry pouch for monitoring devices, blue	200/case



PBAG5060

# Patient kits<sup>5</sup>

The kits contain single-use patient apparel, which patients can wear during the diagnostic procedure or day surgery. All components are conveniently packed in a single bag that can be given to the patient upon arrival. Due to its packaging, organising storage and replenishing stock are easy to handle, and you only have to order one item instead of four or five individual components.



Item number	Size	Description	Packaging
KIT1030	L	Kimono, blue Patient gown with short sleeves, blue Underwear, dark blue Clip cap, blue Socks, red	40/case
KIT1035 KIT1035XXL	L 2XL	Kimono, purple Patient gown with short sleeves, purple Underwear, dark blue Clip cap, blue Socks, red	40/case
KIT3020 KIT3020XXL	L 2XL	Patient gown with short sleeves, blue Underwear, dark blue Clip cap, white Shoe covers, white	50/case
NONE27KITK	Suitable for children (3 to 5 years)	Patient gown for child with short sleeves, blue Underwear, dark blue Clip cap, white Shoe covers, white	50/case
NONE35LKIT NONE352XLKIT	L 2XL	Patient gown with short sleeves, blue Scrub suit pants, blue Clip cap, white Shoe covers, textured bottom, white and blue	50/case
KIT1050E	L	Kimono, blue Patient gown with short sleeves, blue Scrub suit pants, blue Underwear, dark blue Shoe covers, textured bottom, white and blue Socks, red	25/case
KIT1440E KIT14402XE	L 2XL	Patient gown with short sleeves, front opening, blue Scrub suit pants, blue Clip cap, blue Shoe covers, textured bottom, white and blue Socks, red	30/case
KIT0060M	M	Patient brassiere; dark-blue underwear	100/case



KIT1030

# Medline also offers a wide range of surgical, examination and respiratory masks.

Check out the facemask brochure for more details.

## Facemasks

Surgical, procedural and respirators



Masks at the ready

For more information, contact your Medline account manager or visit us at: [www.medline.eu/uk](http://www.medline.eu/uk).



### Medline Industries Ltd.

3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237

[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

### Medline Ireland Ltd.

1 Grant's Row  
Mount Street Lower  
Dublin 2  
Ireland  
Tel.: +353 1 691 73 06

[www.medline.eu/ie](http://www.medline.eu/ie)  
[ie-customerservice@medline.com](mailto:ie-customerservice@medline.com)

FOLLOW US  

<sup>1</sup> Cover gowns are class I medical devices according to regulation (EU) 2017/745.



<sup>2</sup> Footwear products are class I medical devices according to Regulation (EU) 2017/745.



<sup>3</sup> Disposable scrub suits are class I medical devices according to Regulation (EU) 2017/745.



<sup>4</sup> Headwear products are class I medical devices according to Regulation (EU) 2017/745.



These products are intended to be used by healthcare professionals.  
Before use, please consult instructions and precautions on the corresponding labelling.

<sup>5</sup> Patient kits are procedure packs according to Article 22 of Regulation (EU) 2017/745.

The procedure pack producer is Medline International France SAS, 5 Rue Charles Lindbergh, 44110 Châteaubriant, France.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093, USA. ML179\_EN 01/2024