



Over **80%**
of all perforations in surgical
gloves are not detected.¹

The longer the gloves are used,
the higher the perforation
risk becomes.²

Glove
perforations
can lead to the
transmission
of infectious
pathogens.

DOUBLE-GLOVING



1. Don a dark-coloured
underglove.



2. Put a light-coloured
exterior glove over it
to create clear contrast.



3. Easily spot perforations
in your glove if they occur.

SensiCare® PI Green MSG92xx

- Made from Isolex polyisoprene, 100% latex-free
- External coating to fit the exterior glove
- Internal aloe vera coating
- 0.21 mm fingertip thickness



DermAssure™ Green MSG65xxl

- Made from neoprene, 100% latex-free
- External coating to fit the exterior glove
- Accelerator-free to reduce the risk of allergic contact dermatitis
- 0.20 mm fingertip thickness



Signature Latex Green MSG55xx

- Made from latex
- External coating to fit the exterior glove
- Internal aloe vera coating
- 0.19 mm fingertip thickness



¹ Thomas, S, Agarwala, M, Mehtab, G. Intraoperative glove perforation—single versus double gloving in protection against skin contamination. Post Graduate Medical Journal. 2001; 77:458-460. Available at: <http://pmj.bmj.com/content/77/909/458.full>. Accessed 18 December 2015.

² Partecke, Lars Ivo, Anna-Maria Goerd, Inga Langner, Bernd Jaeger, Ojan Assadian, Claus-Dieter Heidecke, Axel Kramer and Nils-Olaf Huebner. 'Incidence of Microperforation for Surgical Gloves Depends on Duration of Wear.' Infection Control and Hospital Epidemiology 30.5 (2009): 409-14.

For further information, please contact your Medline
account manager or visit our website: www.medline.eu/uk



Medline Industries Ltd

www.medline.eu/uk
uk-customerservice@medline.com

FOLLOW US 

These gloves are class IIa sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



Medline International France SAS
5, rue Charles Lindbergh
44110 Châteaubriant, France



We reserve the right to correct any errors that may occur in this brochure.

© 2020 Medline Industries, Inc. Medline is a registered trademark of Medline Industries Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML977-EN 01/2020