



Over 80%

of all perforations in surgical gloves are not detected.¹

The longer the gloves are used, the higher the perforation risk becomes.²

Glove perforations can lead to the transmission of infectious pathogens.

Double-Gloving:



1. Don a dark-coloured underglove.

2. Put a light-coloured exterior glove over it to create clear contrast.

3. Easily spot perforations in your glove if they occur.

SensiCare® PI Green MSG92xx

- Made from Isolex polyisoprene, 100% latex-free
- External coating to fit the exterior glove
- Internal aloe vera coating
- 0.21 mm fingertip thickness

DermAssure™ Green MSG65xx

- Made from neoprene, 100% latex-free
- External coating to fit the exterior glove
- Accelerator-free to reduce the risk of allergic contact dermatitis
- 0.21 mm fingertip thickness

Signature Latex Green MSG55xx

- Made from latex
- External coating to fit the exterior glove
- Internal aloe vera coating
- 0.19 mm fingertip thickness



For further information, please contact your Medline account manager or visit our website: www.medline.eu/uk

1. Thomas, S, Agarwala, M, Mehtab, G, Intraoperative glove perforation—single versus double gloving in protection against skin contamination. Post Graduate Medical Journal. 2001; 77:458-460. Available at: <http://pmj.bmj.com/content/77/909/458.full>. Accessed 18 December 2015. 2. Partecke, Lars Ivo, Anna-Maria Goerd, Inga Langner, Bernd Jaeger, Ojan Assadian, Claus-Dieter Heidecke, Axel Kramer and Nils-Olaf Huebner. 'Incidence of Microperforation for Surgical Gloves Depends on Duration of Wear.' Infection Control and Hospital Epidemiology 30.5 (2009): 409-14.

These gloves are class IIa sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



Medline International France SAS
5, rue Charles Lindbergh
44110 Châteaubriant, France



We reserve the right to correct errors that may occur in this brochure.

© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML977_EN 07/2021