



Over 80%

of all perforations in surgical gloves are not detected.¹

The longer the gloves are used, the higher the perforation risk becomes.²

Glove perforations can lead to the transmission of infectious pathogens.

Double-Gloving:



1. Don a dark-coloured underglove.

2. Put a light-coloured exterior glove over it to create clear contrast.

3. Easily spot perforations in your glove if they occur.

SensiCare® PI Green MSG92xx

- Made from Isolex polyisoprene, 100% latex-free
- External coating to fit the exterior glove
- Internal aloe vera coating
- 0.21 mm fingertip thickness

DermAssure™ Green MSG65xx

- Made from neoprene, 100% latex-free
- External coating to fit the exterior glove
- Accelerator-free to reduce the risk of allergic contact dermatitis
- 0.19 mm fingertip thickness

Signature Latex Green MSG55xx

- Made from latex
- External coating to fit the exterior glove
- Internal aloe vera coating
- 0.19 mm fingertip thickness



For further information, please contact your Medline account manager or visit our website: www.medline.eu/uk

1. Thomas, S, Agarwala, M, Mehtab, G, Intraoperative glove perforation—single versus double gloving in protection against skin contamination. Post Graduate Medical Journal. 2001; 77:458-460. Available at: <http://pmj.bmj.com/content/77/909/458.full>. Accessed 18 December 2015. 2. Partecke, Lars Ivo, Anna-Maria Goerd, Inga Langner, Bernd Jaeger, Ojan Assadian, Claus-Dieter Heidecke, Axel Kramer and Nils-Olaf Huebner. 'Incidence of Microperforation for Surgical Gloves Depends on Duration of Wear.' Infection Control and Hospital Epidemiology 30.5 (2009): 409-14.

These gloves are class IIa sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France



We reserve the right to correct errors that may occur in this brochure.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML977_EN 06/2023