

Eco Vinyl



First-generation clear vinyl powder-free exam gloves

These gloves are an economical option for general patient care and other low-risk environments.

Benefits and features

- Material containing traditional polyvinyl chloride (PVC) resin and plasticisers
- Latex-free and accelerator-free
- Acceptable quality level (AQL) for EN455-1: 1.5
- Force at break (FAB) median value: 4.4 newtons
- Personal protective equipment (PPE)-certified

Ordering information

Item number	Size	Packaging
MGE323	X-small	150 ea/bx, 10 bx/cs
MGE324	Small	150 ea/bx, 10 bx/cs
MGE325	Medium	150 ea/bx, 10 bx/cs
MGE326	Large	150 ea/bx, 10 bx/cs
MGE327	X-large	130 ea/bx, 10 bx/cs




Sustainable packaging technology

Flat-layered packaging

Colour-coded sizing

100% recyclable materials



Recommended for use in:



Non-critical care unit



Household cleaning

Regulations and quality standards

- Class 1 non-sterile medical device, per Medical Device Regulation (EU) 2017/745 and the UK Medical Device Regulation
- Compliant with standards EN 455-1, 2, 3, and 4
- Manufactured according to EN ISO 13485
- Meets requirements of the Personal Protective Equipment Regulation (EU) 2016/45 Category III, as brought into UK law and amended
- Compliant with standards EN ISO 21420; EN ISO 374-1; EN ISO 374-2; EN 16523-1; EN ISO 374-4 and EN ISO 374-5
- Tested for viral penetration, in accordance with EN ISO 374-5
- Compliant with Regulation (EU) No. 10/2011 for plastic materials and articles intended to come in contact with food except 'fatty food'



CE 2797 BSI Group The Netherlands B.V.
1066 EP, Amsterdam, Netherlands



Explore how Medline's gloves range, including Eco Vinyl, makes a difference around the world. Contact your local Medline representative or visit our website: www.medline.eu/uk

This product is a class I non-sterile medical device intended for use by healthcare professionals. Before use, consult the instructions and precautions on the corresponding labeling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2024 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1495_EN 06/2024

Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France