


Endobronchial Double-Lumen Tube




Medline's ENDOSID® Endobronchial Double-Lumen Tube—available for both the left and right bronchus—is made from soft atraumatic polyvinyl chloride (PVC) for enhanced patient comfort. This tube features a unique locking mechanism on the elbow connectors, enabling manual opening or closing of the bronchial and tracheal lumens with a simple rotary knob. This design streamlines one-lung ventilation, eliminating the need for Magill forceps during the process. Additional benefits and features include:

- Sterility via ethylene oxide (EO)
- Free from natural rubber latex
- X-ray contrast stripes
- Intubation stylet
- 2 suction catheters
- Soft clear PVC
- Shelf life of 5 years






Manually Lockable Elbow Connectors



Blue Transparent Pivot Balloons



Transparent Tracheal Cuff

Ordering Information

ENDOSID Endobronchial Tube Set Left and Right-Sided

Item number	Size	Design	Pack factor
847028-RZ	28 ch	Right-sided	1
847032-RZ	32 ch	Right-sided	1
847035-RZ	35 ch	Right-sided	1
847037-RZ	37 ch	Right-sided	1
847039-RZ	39 ch	Right-sided	1
847041-RZ	41 ch	Right-sided	1
847028-LZ	28 ch	Left-sided	1
847032-LZ	32 ch	Left-sided	1
847035-LZ	35 ch	Left-sided	1
847037-LZ	37 ch	Left-sided	1
847039-LZ	39 ch	Left-sided	1
847041-LZ	41 ch	Left-sided	1

CH= Charrière



847037-RZ



847037-LZ

We offer a complete range, applicable to the left and right bronchus.

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk



Medline Industries Ltd.
 3rd Floor
 Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington WA4 6HL
 United Kingdom
 Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd.
 1 Grant's Row
 Mount Street Lower
 Dublin 2
 Ireland
 Tel.: +353 1 691 73 06
www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

ENDOSID® Endotracheal Tubes are class IIa sterile medical devices intended to be used by healthcare professionals. Before use, please consult instructions and precautions on the corresponding labelling.



Shanghai International Holding Corp GmbH (Europe)
 Eiffelstraße 80
 20537 Hamburg, Germany



We reserve the right to correct errors that may occur in this brochure.

© 2024 Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1477_EN 04/2024.