



Non-sterile Flexible Enema Kit

A convenient,
comprehensive
and effective
packaged
solution

Medline's Enema Bag Set provides the end user with a 1,500-ml enema bag, outfitted with a spill-proof valve.

Benefits and features

- Facilitates one-handed operation by way of an easy-to-use funnel top with a plastic shut-off clamp
- Allows for safe and easy insertion, as the 150-cm soft vinyl tube has a pre-lubricated, smooth, non-traumatic tip with rounded eyes
- Comes with a sachet of gentle castile soap and an impervious underpad
- Individually packaged in a box; non-sterile



Ordering information

| Item code | Packaging |
|------------|-----------|
| DYNDE70100 | 48 |

Intended medical uses of the enema bag include:

- Bowel emptying prior to surgery (e.g. colonoscopy)
- Constipation relief, if conventional methods are not effective
- Diagnostic tests

Learn more. Contact your Medline account manager or visit: www.medline.eu/uk

This product is a class I non-sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



All rights reserved. We reserve the right to correct errors that may have occurred.
© 2021 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA.
ML1168_EN 07/2021