



# Non-Sterile Flexible Enema Kit

A convenient,  
comprehensive  
and effective  
packaged  
solution

Medline's Enema Bag Set provides the end user with a 1,500-ml enema bag, outfitted with a spill-proof valve.

## Benefits and features

- Designed to facilitate one-handed operation, by way of an easy-to-use funnel top with a plastic shut-off clamp
- Designed to allow for safe and easy insertion, as the 150-cm soft vinyl tube has a pre-lubricated, smooth, non-traumatic tip with rounded eyes
- Comes with a sachet of gentle castile soap and an underpad, designed to be impervious
- Individually packaged in a box; non-sterile



## Ordering information

Item code	Packaging
DYNDE70100A	48

## Intended medical uses of the enema bag include:

- Bowel emptying prior to surgery (e.g. colonoscopy)
- Constipation relief, if conventional methods are not effective
- Diagnostic tests.

**Learn more.** Contact your Medline account manager or visit: [www.medline.eu/uk](http://www.medline.eu/uk)

This product is a class I non-sterile medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



All rights reserved. We reserve the right to correct errors that may have occurred.  
© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA.  
ML1168\_EN 11/2023