


















# Exam Gloves Ordering Guide

Product name	Item numbers	Sizes	Colour	Texture	Thickness			Cuff length	Description	Packaging
					Fingers	Palm	Cuff			
 Restore® Nitrile*	OAT-XS OAT-S OAT-M OAT-L OAT-XL	X-Small Small Medium Large X-Large	White	Textured fingertips	0.09 mm	0.07 mm	0.05 mm	242 mm	Nitrile, non-powdered, non-sterile, with colloidal oatmeal inner coating	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs
 Restore® Sense*	OATN32XS OATN32S OATN32M OATN32L OATN32XL	X-Small Small Medium Large X-Large	Pale mint green	Textured fingertips	0.16 mm	0.12 mm	0.10 mm	250 mm	Nitrile, non-powdered, non-sterile, with colloidal oatmeal inner coating	XS – L: 250 ea/bx, 10 bx/cs XL: 230 ea/bx, 10 bx/cs
 SensiCare® Free*	MDSFREE-XS MDSFREE-S MDSFREE-M MDSFREE-L MDSFREE-XL	X-Small Small Medium Large X-Large	Dark blue	Textured fingertips	0.09 mm	0.07 mm	0.05 mm	250 mm	Nitrile, non-powdered, non-sterile, accelerator-free	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs
 Generation Pink® Nitrile*	PINKE2500 PINKE2501 PINKE2502 PINKE2503 PINKE2504	X-Small Small Medium Large X-Large	Pink	Textured fingertips	0.09 mm	0.07 mm	0.05 mm	240 mm	Nitrile, non-powdered, non-sterile	XS – L: 250 ea/bx, 10 bx/cs XL: 230 ea/bx, 10 bx/cs
 SensiCare® Ice*	486800 486801 486802 486803 486804	X-Small Small Medium Large X-Large	Blue	Textured fingertips	0.09 mm	0.07 mm	0.05 mm	250 mm	Nitrile, non-powdered, non-sterile	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs
 SensiCare® Silk*	MDS9083 MDS9084 MDS9085 MDS9086 MDS9087	X-Small Small Medium Large X-Large	Dark blue	Textured fingertips	0.09 mm	0.07 mm	0.05 mm	243 mm	Nitrile, non-powdered, non-sterile	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs
 MediGuard® Blue*	EM5050 EM5051 EM5052 EM5053 EM5054	X-Small Small Medium Large X-Large	Light blue	Textured fingertips	0.09 mm	0.07 mm	0.05 mm	248 mm	Nitrile, non-powdered, non-sterile	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs
 SensiCare® Pro*	PRO400XS PRO400S PRO400M PRO400L PRO400XL	X-Small Small Medium Large X-Large	Blue	Fully textured	0.23 mm	0.18 mm	0.11 mm	402 mm	Nitrile, non-powdered, non-sterile	XS – L: 50 ea/bx, 10 bx/cs XL: 50 ea/bx, 10 bx/cs

\*These gloves have been tested and approved for use with chemotherapy drugs.

	Product name	Item numbers	Sizes	Colour	Texture	Thickness			Cuff length	Description	Packaging
						Fingers	Palm	Cuff			
	Aloetouch® Long Cuff*	MDS195184 MDS195185 MDS195186 MDS195187	Small Medium Large X-Large	Green	Fully textured	0.16 mm	0.14 mm	0.10 mm	307 mm	Nitrile, non-powdered, non-sterile, with aloe inner coating	<b>S – L:</b> 50 ea/bx, 10 bx/cs <b>XL:</b> 50 ea/bx, 10 bx/cs
	SensiCare® Long Cuff*	MDSEXTXS MDSEXTS MDSEXTM MDSEXTL MDSEXTXL	X-Small Small Medium Large X-Large	Blue	Textured fingertips	0.14 mm	0.10 mm	0.08 mm	284 mm	Nitrile, non-powdered, non-sterile	<b>XS – L:</b> 100 ea/bx, 10 bx/cs <b>XL:</b> 90 ea/bx, 10 bx/cs
	Curad® 3G Stretch Vinyl	CUR8233 CUR8234 CUR8235 CUR8236 CUR8237	X-Small Small Medium Large X-Large	White	Smooth	0.10 mm	0.08 mm	0.06 mm	242 mm	Stretch vinyl, non-powdered non-sterile	<b>XS – L:</b> 100 ea/bx, 10 bx/cs <b>XL:</b> 90 ea/bx, 10 bx/cs
	Accutouch	MDS192073 MDS192074 MDS192075 MDS192076 MDS192077	X-Small Small Medium Large X-Large	Clear	Smooth	0.10 mm	0.08 mm	0.06 mm	254 mm	Synthetic vinyl, non-powdered, non-sterile	<b>XS – L:</b> 100 ea/bx, 10 bx/cs <b>XL:</b> 100 ea/bx, 10 bx/cs
	Eco Vinyl	MGE323 MGE324 MGE325 MGE326 MGE327	X-Small Small Medium Large X-Large	Clear	Smooth	0.09 mm	0.07 mm	0.05 mm	240 mm	Synthetic vinyl, non-powdered, non-sterile	<b>S – L:</b> 150 ea/bx, 10 bx/cs <b>XL:</b> 130 ea/bx, 10 bx/cs
	Curad® Latex	CURE8103 CURE8104 CURE8105 CURE8106 CURE8107	X-Small Small Medium Large X-Large	Beige	Textured fingertips	0.13 mm	0.11 mm	0.09 mm	246 mm	Latex, non-powdered, non-sterile	<b>XS – L:</b> 100 ea/bx, 10 bx/cs <b>XL:</b> 90 ea/bx, 10 bx/cs
	Medline Sterile Procedure Nitrile*	MDS2294 MDS2295 MDS2296 MDS2297	Small Medium Large X-Large	Blue	Textured fingertips	0.14 mm	0.11 mm	0.09 mm	296 mm	Nitrile, non-powdered, sterile	<b>S – L:</b> 50 pairs/bx, 6 bx/cs <b>XL:</b> 50 pairs/bx, 6 bx/cs

\*These gloves have been tested and approved for use with chemotherapy drugs.



**Medline Industries Ltd**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
Fax: +44 844 334 5238

[www.medline.eu/uk](http://www.medline.eu/uk) | [uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

FOLLOW US  

We reserve the right to correct errors that may occur in this brochure.

© 2020 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML313-EN 03/2020.