


















# Guida guanti da esame

Nome Prodotto	Codice	Misura	Colore	Testurizzazione	Spessori			Lunghezza polsino	Descrizione	Confezionamento
					Dito	Palmo	Polsino			
 Restore® Nitrile*	OAT-XS OAT-S OAT-M OAT-L OAT-XL	X-Small Small Medium Large X-Large	Bianco crema	Polpastrello	0.09 mm	0.07 mm	0.05 mm	242 mm	Nitrile Depolverato Non-Sterile Trattato internamente con avena colloidale	<b>XS - L:</b> 200 Scatola, 10 Scatole/Cartone <b>XL:</b> 180 Scatola, 10 Scatole/Cartone
 Restore® Sense*	OATN32XS OATN32S OATN32M OATN32L OATN32XL	X-Small Small Medium Large X-Large	Verde menta chiaro	Polpastrello	0.16 mm	0.12 mm	0.10 mm	250 mm	Nitrile Depolverato Non-Sterile Trattato internamente con avena colloidale	<b>XS - L:</b> 250 Scatola, 10 Scatole/Cartone <b>XL:</b> 230 Scatola, 10 Scatole/Cartone
 SensiCare® Free*	MDSFREE-XS MDSFREE-S MDSFREE-M MDSFREE-L MDSFREE-XL	X-Small Small Medium Large X-Large	Blu scuro	Polpastrello	0.09 mm	0.07 mm	0.05 mm	250 mm	Nitrile Depolverato Non-Sterile Acceleratori free	<b>XS - L:</b> 200 Scatola, 10 Scatole/Cartone <b>XL:</b> 180 Scatola, 10 Scatole/Cartone
 Generation Pink® Nitrile*	PINKE2500 PINKE2501 PINKE2502 PINKE2503 PINKE2504	X-Small Small Medium Large X-Large	Rosa	Polpastrello	0.09 mm	0.07 mm	0.05 mm	240 mm	Nitrile Depolverato Non-Sterile	<b>XS - L:</b> 250 Scatola, 10 Scatole/Cartone <b>XL:</b> 230 Scatola, 10 Scatole/Cartone
 SensiCare® Ice*	486800 486801 486802 486803 486804	X-Small Small Medium Large X-Large	Blu	Polpastrello	0.09 mm	0.07 mm	0.05 mm	250 mm	Nitrile Depolverato Non-Sterile	<b>XS - L:</b> 200 Scatola, 10 Scatole/Cartone <b>XL:</b> 180 Scatola, 10 Scatole/Cartone
 SensiCare® Silk*	MDS9083 MDS9084 MDS9085 MDS9086 MDS9087	X-Small Small Medium Large X-Large	Blu scuro	Polpastrello	0.09 mm	0.07 mm	0.05 mm	243 mm	Nitrile Depolverato Non-Sterile	<b>XS - L:</b> 200 Scatola, 10 Scatole/Cartone <b>XL:</b> 180 Scatola, 10 Scatole/Cartone
 MediGuard® Blue*	EM5050 EM5051 EM5052 EM5053 EM5054	X-Small Small Medium Large X-Large	Azzurro	Polpastrello	0.09 mm	0.07 mm	0.05 mm	248 mm	Nitrile Depolverato Non-Sterile	<b>XS - L:</b> 200 Scatola, 10 Scatole/Cartone <b>XL:</b> 180 Scatola, 10 Scatole/Cartone
 SensiCare® Pro*	PRO400XS PRO400S PRO400M PRO400L PRO400XL	X-Small Small Medium Large X-Large	Blu	Totale	0.23 mm	0.18 mm	0.11 mm	402 mm	Nitrile Depolverato Non-Sterile	<b>XS - L:</b> 50 Scatola, 10 Scatole/Cartone <b>XL:</b> 50 Scatola, 10 Scatole/Cartone

\*I guanti sono stati testati e approvati per essere utilizzati con farmaci chemioterapici.

Nome Prodotto	Codice	Misura	Colore	Testurizzazione	Spessori			Lunghezza polsino	Descrizione	Confezionamento	
					Dito	Palmo	Polsino				
	Aloetouch® Long Cuff*	MDS195184 MDS195185 MDS195186 MDS195187	Small Medium Large X-Large	Verde	Totale	0.16 mm	0.14 mm	0.10 mm	307 mm	Nitrile Depolverato Non-Sterile Con rivestimento interno in Aloe	<b>S - L:</b> 50 Scatola, 10 Scatole/Cartone <b>XL:</b> 50 Scatola, 10 Scatole/Cartone
	SensiCare® Long Cuff*	MDSEXTXS MDSEXTS MDSEXTM MDSEXTL MDSEXTXL	X-Small Small Medium Large X-Large	Blu	Polpastrello	0.14 mm	0.10 mm	0.08 mm	284 mm	Nitrile Depolverato Non-Sterile	<b>XS - L:</b> 100 Scatola, 10 Scatole/Cartone <b>XL:</b> 90 Scatola, 10 Scatole/Cartone
	Curad® 3G Stretch Vinyl	CUR8233 CUR8234 CUR8235 CUR8236 CUR8237	X-Small Small Medium Large X-Large	Bianco	Liscio	0.10 mm	0.08 mm	0.06 mm	242 mm	Vinile elasticizzato Depolverato Non-Sterile	<b>XS - L:</b> 100 Scatola, 10 Scatole/Cartone <b>XL:</b> 90 Scatola, 10 Scatole/Cartone
	Accutouch	MDS192073 MDS192074 MDS192075 MDS192076 MDS192077	X-Small Small Medium Large X-Large	Trasparente	Liscio	0.10 mm	0.08 mm	0.06 mm	254 mm	Vinile sintetico Depolverato Non-Sterile	<b>XS - L:</b> 100 Scatola, 10 Scatole/Cartone <b>XL:</b> 100 Scatola, 10 Scatole/Cartone
	Eco Vinyl	MGE323 MGE324 MGE325 MGE326 MGE327	X-Small Small Medium Large X-Large	Trasparente	Liscio	0.09 mm	0.07 mm	0.05 mm	240 mm	Vinile sintetico Depolverato Non-Sterile	<b>S - L:</b> 150 Scatola, 10 Scatole/Cartone <b>XL:</b> 130 Scatola, 10 Scatole/Cartone
	Curad® Latex	CURE8103 CURE8104 CURE8105 CURE8106 CURE8107	X-Small Small Medium Large X-Large	Beige	Polpastrello	0.13 mm	0.11 mm	0.09 mm	246 mm	Lattice Depolverato Non-Sterile	<b>XS - L:</b> 100 Scatola, 10 Scatole/Cartone <b>XL:</b> 90 Scatola, 10 Scatole/Cartone
	Medline Sterile Procedure Nitrile*	MDS2294 MDS2295 MDS2296 MDS2297	Small Medium Large X-Large	Blu	Polpastrello	0.14 mm	0.11 mm	0.09 mm	296 mm	Nitrile, Depolverato Sterile	<b>S - L:</b> 50 Paia/Scatola, 6 Scatole/Cartone <b>XL:</b> 50 Paia/Scatola, 6 Scatole/Cartone

\*I guanti sono stati testati e approvati per essere utilizzati con farmaci chemioterapici.



**Medline International Italy s.r.l. Unipersonale**  
Piazzale della Resistenza 3  
50018 Scandicci  
Firenze  
Italia  
Tel.: +39 55 776 65 11  
Fax: +39 55 34 01 12  
[www.medline.eu/it](http://www.medline.eu/it)  
[it-customerservice@medline.com](mailto:it-customerservice@medline.com)

**Medline International Switzerland Sàrl**  
1 Place de Longemalle  
1204 Ginevra  
Svizzera  
Tel. +41 848 244 433  
Fax. +41 848 244 100  
[www.medline.eu/it](http://www.medline.eu/it)  
[ch-customerservice@medline.com](mailto:ch-customerservice@medline.com)

FOLLOW US  