

## Standard nitrile



### MediGuard Special Nitrile

<b>Item number</b>	MGSNXS MGSNS MGSNM MGSNL MGSNXL
<b>Sizes</b>	XS – XL
<b>Colour</b>	Violet blue
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	7.0 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	240 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.10; 0.08; 0.07 (± 0.03) mm
<b>Packaging</b>	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



### SensiCare® Silk

<b>Item number</b>	MDS9083 MDS9084 MDS9085 MDS9086 MDS9087
<b>Sizes</b>	XS – XL
<b>Colour</b>	Violet blue
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	6.7 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	244 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.09; 0.07; 0.05 mm
<b>Packaging</b>	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



### SensiCare® Ice

<b>Item number</b>	486800 486801 486802 486803 486804
<b>Sizes</b>	XS – XL
<b>Colour</b>	Violet blue
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	8.1 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	249 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.09; 0.07; 0.05 mm
<b>Packaging</b>	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



# Exam Gloves Ordering Guide

## Special nitrile



### SensiCare® Free Accelerator-Free

<b>Specialty</b>	Accelerator-Free
<b>Item number</b>	MDSFREE-XS MDSFREE-S MDSFREE-M MDSFREE-L MDSFREE-XL
<b>Sizes</b>	XS – XL
<b>Colour</b>	Blue
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	8.2 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	243 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.09; 0.07; 0.05 mm
<b>Packaging</b>	XS – L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs



### Restore® Sense Colloidal Oatmeal Inner Coating

<b>Item number</b>	OATN32XS OATN32S OATN32M OATN32L OATN32XL
<b>Sizes</b>	XS – XL
<b>Colour</b>	Lemon green
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	6.8 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	244 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.05; 0.05; 0.045 mm
<b>Packaging</b>	XS – L: 250 ea/bx, 10 bx/cs XL: 230 ea/bx, 10 bx/cs



### Generation Pink Nitrile Breast cancer awareness

<b>Item number</b>	PINKE2500 PINKE2501 PINKE2502 PINKE2503 PINKE2504
<b>Sizes</b>	XS – XL
<b>Colour</b>	Pink
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	6.7 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	240 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.09; 0.07; 0.05 mm
<b>Packaging</b>	XS – L: 250 ea/bx, 10 bx/cs XL: 230 ea/bx, 10 bx/cs



### Sterile Procedure Nitrile Sterility

<b>Item number</b>	MDS2294 MDS2295 MDS2296 MDS2297
<b>Sizes</b>	S – XL
<b>Colour</b>	Blue
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	7.8 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	302 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.14; 0.011; 0.09 mm
<b>Packaging</b>	S – L: 50 pairs/bx, 6 bx/cs XL: 50 pairs/bx, 6 bx/cs

## Long-cuff nitrile



### SensiCare® Extended Cuff

<b>Item number</b>	MDSEXTXS MDSEXTS MDSEXTM MDSEXTL MDSEXTXL
<b>Sizes</b>	XS – XL
<b>Colour</b>	Blue
<b>Texture</b>	Fingertips
<b>Force at break (Unchallenged)</b>	12.4 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	280 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.14; 0.10; 0.08 mm
<b>Packaging</b>	XS – L: 100 ea/bx, 10 bx/cs XL: 90 ea/bx, 10 bx/cs



### SensiCare® Pro

<b>Item number</b>	PRO400XS PRO400S PRO400M PRO400L PRO400XL
<b>Sizes</b>	XS – XL
<b>Colour</b>	Blue
<b>Texture</b>	Fully textured
<b>Force at break (Unchallenged)</b>	20.4 N
<b>AQL</b>	1.0
<b>Cuff length (Size M)</b>	412 mm
<b>Thickness (Fingers; Palm; Cuff)</b>	0.23; 0.18; 0.11 mm
<b>Packaging</b>	XS – L: 50 ea/bx, 10 bx/cs XL: 50 ea/bx, 10 bx/cs

## Vinyl



### Eco Vinyl

MGE323  
MGE324  
MGE325  
MGE326  
MGE327

#### Item number

#### Sizes

XS - XL

#### Colour

Clear

#### Texture

Smooth

#### Force at break (Unchallenged)

4.4 N

#### AQL

1.5

#### Cuff length (Size M)

245 mm

#### Thickness

(Fingers; Palm; Cuff)

0.11; 0.09; 0.07 (± 0.03) mm

#### Packaging

XS - L: 150 ea/box, 10 bx/cs  
XL: 130 ea/box, 10 bx/cs



### ACCUtouch®

MDS192074  
MDS192075  
MDS192076  
MDS192077

S - XL

Clear

Smooth

4.7 N

1.5

250 mm

0.10; 0.08; 0.06 mm

S - XL: 100 ea/box, 10 bx/cs



### Curad® 3G

CUR8233  
CUR8234  
CUR8235  
CUR8236  
CUR8237

XS - XL

White

Smooth

3.7 N

1.5

248 mm

0.10; 0.08; 0.06 mm

XS - L: 100 ea/box, 10 bx/cs  
XL: 90 ea/box, 10 bx/cs

## Latex



### Curad® Latex

CURE8103  
CURE8104  
CURE8105  
CURE8106  
CURE8107

XS - XL

White

Smooth

7.5 N

1.5

249 mm

0.13; 0.11; 0.09 mm

XS - L: 100 ea/box, 10 bx/cs  
XL: 90 ea/box, 10 bx/cs

We provide  
you the right  
gloves for the  
right task.

## Wire box holders



### Single box holder

MDS191094

#### Item number

#### Material

Coated wire

#### Colour

White

#### Packaging

10 holders per case



### Double box holder

MDS192094B

Coated wire

White

6 holders per case



### Triple box holder

MDS193094

Coated wire

White

10 holders per case

## Plexiglass box holders



### Single box holder

INTL191096

#### Item number

#### Material

Plexiglass

#### Colour

Clear

#### Packaging

12 holders per case



### Triple box holder

INTL193096

Plexiglass

Clear

6 holders per case