

## Standard nitrile



	MediGuard Special Nitrile	SensiCare® Ice
<b>Item number</b>	MGSNXS MGSNS MGSNM MGSNL MGSNXL	486800 486801 486802 486803 486804
<b>Sizes</b>	XS - XL	XS - XL
<b>Colour</b>	Violet blue	Violet blue
<b>Texture</b>	Fingertips	Fingertips
<b>PPE Type</b>	Type B	Type C
<b>Force at break (Unchallenged)</b>	7.0 N	8.1 N
<b>AQL</b>	≤ 1.0	≤ 1.0
<b>Cuff length (Size M)</b>	240 mm	249 mm
<b>Packaging</b>	XS - L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs	XS - L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs

## Long-Cuff nitrile



	SensiCare® 300 Long	SensiCare® Extended Cuff
<b>Item number</b>	MDSMID-S MDSMID-M MDSMID-L MDSMID-XL	MDSEXTXS MDSEXTS MDSEXTM MDSEXTL MDSEXTXL
<b>Sizes</b>	S - XL	XS - XL
<b>Colour</b>	Blue	Blue
<b>Texture</b>	Fingertips	Fingertips
<b>PPE Type</b>	Type B	Type B
<b>Force at break (Unchallenged)</b>	11.0 N	12.4 N
<b>AQL</b>	≤ 1.0	≤ 1.0
<b>Cuff length (Size M)</b>	300 mm	280 mm
<b>Packaging</b>	S - L: 100 ea/bx, 10 bx/cs XL: 90 ea/bx, 10 bx/cs	XS - L: 100 ea/bx, 10 bx/cs XL: 90 ea/bx, 10 bx/cs



# Exam Gloves Ordering Guide

## Special nitrile



Specialty	SensiCare® Free Accelerator-Free	Restore® Sense Colloidal Oatmeal Inner Coating	Generation Pink Nitrile Breast cancer awareness	Sterile Procedure Nitrile Sterility	SensiCare® Pro Extra-long (400mm) for sterilisation and laboratory handling
<b>Item number</b>	MDSFREE-XS MDSFREE-S MDSFREE-M MDSFREE-L MDSFREE-XL	OATN32XS OATN32S OATN32M OATN32L OATN32XL	PINKE2501 PINKE2502 PINKE2503	MDS2294 MDS2295 MDS2296 MDS2297	PRO400XS PRO400S PRO400M PRO400L PRO400XL
<b>Sizes</b>	XS - XL	XS - XL	S, M, L	S - XL	XS - XL
<b>Colour</b>	Blue	Lemon green	Pink	Blue	Blue
<b>Texture</b>	Fingertips	Fingertips	Fingertips	Fingertips	Fully textured
<b>PPE Type</b>	Type B	Type C	Type C	Type B	Type B
<b>Force at break (Unchallenged)</b>	8.2 N	6.8 N	6.7 N	7.8 N	20.4 N
<b>AQL</b>	≤ 1.0	≤ 1.0	≤ 1.0	≤ 1.0	≤ 1.0
<b>Cuff length (Size M)</b>	243 mm	244 mm	240 mm	302 mm	412 mm
<b>Packaging</b>	XS - L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs	XS - L: 250 ea/bx, 10 bx/cs XL: 230 ea/bx, 10 bx/cs	S - L: 250 ea/bx, 10 bx/cs	S - L: 50 pairs/bx, 6 bx/cs XL: 50 pairs/bx, 6 bx/cs	XS - L: 50 ea/bx, 10 bx/cs XL: 50 ea/bx, 10 bx/cs

## Sustainability



SensiCare® Silk Cuff 1 <sup>st</sup>
MDS9083CF MDS9084CF MDS9085CF MDS9086CF MDS9087CF
XS - XL
Violet blue
Fingertips
Type C
6.7 N
≤ 1.0
244 mm
XS - L: 200 ea/bx, 10 bx/cs XL: 180 ea/bx, 10 bx/cs

## Vinyl



	Eco Vinyl	ACCUtouch®	Curad® 3G
<b>Item number</b>	MGE323 MGE324 MGE325 MGE326 MGE327	MDS192074 MDS192075 MDS192076 MDS192077	CUR8233 CUR8234 CUR8235 CUR8236 CUR8237
<b>Sizes</b>	XS – XL	S – XL	XS – XL
<b>Colour</b>	Clear	Clear	White
<b>Texture</b>	Smooth	Smooth	Smooth
<b>PPE Type</b>	Type C	Type C	Type C
<b>Force at break (Unchallenged)</b>	4.4 N	4.7 N	3.7 N
<b>AQL</b>	≤ 1.5	≤ 1.5	1.5
<b>Cuff length (Size M)</b>	245 mm	250 mm	248 mm
<b>Packaging</b>	XS – L: 150 ea/bx, 10 bx/cs XL: 130 ea/bx, 10 bx/cs	S – XL: 100 ea/bx, 10 bx/cs	XS – L: 100 ea/bx, 10 bx/cs XL: 90 ea/bx, 10 bx/cs

## Latex



	Curad® Latex
<b>Item number</b>	CURE8103 CURE8104 CURE8105 CURE8106 CURE8107
<b>Sizes</b>	XS – XL
<b>Colour</b>	White
<b>Texture</b>	Smooth
<b>PPE Type</b>	Type C
<b>Force at break (Unchallenged)</b>	7.5 N
<b>AQL</b>	≤ 1.5
<b>Cuff length (Size M)</b>	249 mm
<b>Packaging</b>	XS – L: 100 ea/bx, 10 bx/cs XL: 90 ea/bx, 10 bx/cs

**We provide you the right gloves for the right task.**

## Wire box holders



	Single box holder	Double box holder	Triple box holder
<b>Item number</b>	MDS191094	MDS192094B	MDS193094
<b>Material</b>	Coated wire	Coated wire	Coated wire
<b>Colour</b>	White	White	White
<b>Packaging</b>	10 holders per case	6 holders per case	10 holders per case

## Plexiglass box holders



	Single box holder	Triple box holder
<b>Item number</b>	INTL191096	INTL193096
<b>Material</b>	Plexiglass	Plexiglass
<b>Colour</b>	Clear	Clear
<b>Packaging</b>	12 holders per case	6 holders per case

## Cuff 1<sup>st</sup> box holders



	Cuff 1 <sup>st</sup> triple box holders and holder bar
<b>Item number</b>	MDSCF3
<b>Material</b>	Plastic
<b>Colour</b>	White
<b>Packaging</b>	2 sets of Cuff 1st box holders and 2 holder bars per case

Medline exam gloves are class I non-sterile medical devices.

Medline International France SAS  
5 rue Charles Lindbergh  
44110 Châteaubriant, France



Sterile Procedure Nitrile gloves are Class I sterile medical devices.



Medline International France SAS  
5 rue Charles Lindbergh  
44110 Châteaubriant, France



We reserve the right to correct errors that may occur in this brochure.

© 2025 Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML313\_EN 01/2025.