



EXUDERM SATIN

Hydrocolloid dressing

Product overview

Exuderm Satin Hydrocolloid Dressing is easy to use with enhanced design features that make both, caregivers and patients happy.

- **Long wear, less change:** 7-day wear time, depending on amount of drainage.
- **Translucent backing:** Allows for wound observation.
- **Satin finish backing:** Low friction for longer wear.
- **Low profile, tapered edge:** Prevents rolling and raised edges.
- **Shaped to fit:** The Exuderm Satin sacrum dressing is anatomically designed to conform to the sacral area, providing secure management of pressure sores and minimizing leakage.



Item Details

Product	Item No.	Size	Order Qty	PIPP/ FP10 Code	NHS SC Code
EXUDERM SATIN hydrocolloid dressing with beveled edge square	MSC5422	5x5 cm	1 box of 20 dressings	398-6981	ELM220
	MSC5444	10x10 cm	1 box of 10 dressings	398-6999	ELM218
	MSC5466	15x15 cm	1 box of 5 dressings	398-7005	ELM219
	MSC5488	20x20 cm	1 box of 5 dressings	398-7013	ELM221
EXUDERM SATIN hydrocolloid dressing with beveled edge sacral	MSC5470	10x9 cm	1 box of 10 dressings	398-7021	
	MSC5475	15x16.5 cm	1 box of 5 dressings	398-7039	

Indications

Pressure injuries, partial and full thickness wounds, leg ulcers, donor sites, wounds with light to moderate drainage, lacerations and abrasions, and first and second degree burns

Contact details

Customer Services

Tel: +44 844 334 5237
uk-customerservice@medline.com

Medline Industries Ltd

3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Fax: +44 844 334 5238
www.medline.eu/uk

FEATURES

- Squared and anatomical shapes for correct wound fit.
- Satin finished backing for low friction and longer wear time.
- Cost effective solution for low to moderate exuding wounds.
- Beveled edges prevent rolling and raised edges.