

Patient Falls Alarms



Let our patient falls alarms help keep your patients safe and secure by alerting you when a patient leaves the chair or bed.



*Image shown may represent a configuration not available in all markets

Set the tone for fast responses

When medical attention is needed, your staff needs to know right away. For patients assessed as having a high falls risk, Medline's patient alarm offers a simple and low cost intervention, available with a two-year warranty.

Alarm Monitor

- **Alert staff and keep your patients safe** with auto-sensing technology that monitors patient movement and triggers an alarm.
- The alarm monitor is **compatible** with most corded sensor pads and nurse call systems.
- **Customisable settings** allow you to set the volume and time delay based on patient and staff needs.



Sensor Pads

- When this **weight-sensitive** pad is connected to a monitor unit, it **activates an alarm** if the patient leaves a chair or bed.
- The pad is **easy to clean**, thanks to durable and high-grade vinyl.
- It is also **comfortable** for patients, thanks to an additional foam layer.
- It is intended for **single-patient use** on a chair or bed (underneath the fitted sheet).



SKU	Description	Size	Packaging
ALARM01	Falls Prevention Alarm Monitor, wired	8.5 x 13 cm	1 ea/box

SKU	Description	Size	Packaging
SENPAD01	Chair Sensor Pad, wired, single-patient use	25.4 x 38.1 cm	5 ea/box
SENPAD02	Bed Sensor Pad, wired, single-patient use	25.4 x 76.2 cm	5 ea/box

For additional information on falls prevention and products, please contact your account manager or visit uk.medline.eu