




**EASY-TO-USE, SAFE
ECONOMICAL SOLUTION.**



Flex-One Surgical Suction

The product for your day-to-day surgical suction needs

The Flex-One system offers a large range of reusable products that are economical, easy-to-use, safe and reliable. Flex-One is the solution for all your medical suction requirements. It can be easily adapted into numerous configurations thanks to its mobile carousel with four canisters. Wall mounting can be completed with a rail mount or by simply mounting it directly to the wall.

Economical

- With just 1 liner available (2 L), Flex-One is a simple and efficient solution.

Easy to use

- 2 ports: 1x patient port and 1x vacuum port
- Colour coded to identify liners that contain a hydroscopic filter
- Colour coded to identify the external canister

Safe

- The hydroscopic filter on each canister with a green lid protects the source from aerosols and infected fluids.
- A non-return valve, placed on the liner's patient port, ensures optimal protection for your staff and the patient by avoiding leaks.

Flex-One single-use products¹

Article No.	Description	Quantity
FO200	2 L liner without filter, white lid	66
FO201	2 L liner with filter, green lid	54
FO201PG	2 L liner with filter, with solidifier, green lid	54

Hardware^{2,3,4}

Article No.	Description	Quantity
CF0201	Green bracket, external liner	15
CF0200	White bracket, external liner	15
DYNDSCCLW	Wall plate	20
LCLC1426CE	Metal rail clamp	1
DYNDSCC	Tall roll stand, 1 to 4 canisters	1
ORBASE	Short roll stand, 1 to 2 canisters	1
TUBVAC	Connecting tube for vacuum canister, 30 m	1



BENEFITS

- A large range of reusable products
- Economical solution
- Easy-to-use
- Safe

For additional information about this product, please contact your account manager or visit us at: www.medline.eu/uk



ALWAYS
ON.

Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk
uk-customerservice@medline.com