



# Wound Care


You're challenged daily to care for wounds, helping patients heal and avoid infection. You need high-quality, high-performance products. We can help.

## Gauze

Medline offers a variety of gauze swabs, available in multiple sizes and plies. Gauze is ideal for wound dressings, wound packing and general wound care.


### Woven gauze

Item number	Size	Ply	Sterility	XR	Packaging
GS10P102012XR	10 x 20 cm	12	Sterile	Yes	10 pcs/pack 60 packs/case
GS10P757512XR	7.5 x 7.5 cm	12	Sterile	Yes	10 pcs/pack 60 packs/case
GS2P101012	10 x 10 cm	12	Sterile	No	2 pcs/pack 800 packs/case
GS5P050508	5 x 5 cm	8	Sterile	No	5 pcs/pack 540 packs/case
GSNS101008	10 x 10 cm	8	Non-sterile	No	100 pcs/pack 10 packs/case



### Non-woven gauze

Item number	Size	Ply	Sterility	XR	Packaging
NWNS30101004	10 x 10 cm	4	Non-sterile	No	100 pcs/pack 10 packs/case






Our XR gauze is sewn with barium sulphate thread, which is X-ray detectable.

## ABD pads

ABD pads are superabsorbent pads with a soft non-woven outer layer and a cellulose centre to prevent pooling. These pads come with a hydrophobic back to impede strikethrough as well as sealed edges to prevent leakage and linting.

Item number	Size	Ply	Sterility	XR	Packaging
ABD1010	10 x 10 cm	n/a	Sterile	No	1 pc/pouch 300 pouches/case
ABD1020	10 x 20 cm	n/a	Sterile	No	1 pc/pouch 200 pouches/case
ABD1525	15 x 25 cm	n/a	Sterile	No	1 pc/pouch 180 pouches/case
ABD2020	20 x 20 cm	n/a	Sterile	No	1 pc/pouch 120 pouches/case
ABD2040	20 x 40 cm	n/a	Sterile	No	1 pc/pouch 80 pouches/case



For more information, please contact your account manager or visit us at [www.medline.eu/uk](http://www.medline.eu/uk)



**Medline Industries Ltd**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

**Medline Ireland Ltd**  
1 Grant's Row  
Mount Street Lower  
Dublin 2  
Ireland  
Tel.: +353 1 691 73 06  
[www.medline.eu/ie](http://www.medline.eu/ie)  
[ie-customerservice@medline.com](mailto:ie-customerservice@medline.com)

FOLLOW US

GS10PXXXX12XR is a class IIa medical device. 2797 (BSI) Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France

GSNS101008 and NWNS30101004 are class I non-sterile medical devices. Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France

GSXPXXXXX and ABDXXXX are class I sterile medical devices. 2797 (BSI) Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France

These medical devices are intended to be used by healthcare professionals. Before use, please consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2022 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093, USA. ML1241\_EN 01/2022