

The background of the advertisement is a close-up photograph of the Gemini Sterilisation Wrap. It shows two layers of the material, one in a vibrant blue and the other in a soft pink. Both layers have a fine, uniform perforated texture. The layers are slightly offset, creating a sense of depth and showing the material's flexibility. The lighting is soft and even, highlighting the texture of the wrap.

**A GEM OF
A WRAP.**

Gemini Sterilisation Wrap

Ensure the Integrity of Your Sterilisation Process

GEMINI STERILISATION WRAP HELPS ENSURE THE INTEGRITY OF YOUR STERILISATION PROCESS.

Drapeable and Stronger

Gemini wrap is resistant to tears and punctures, performing better than standard SMS wraps. It is very soft and easy to fold, which makes it simple to work with.

Low Linting and Better Adhesion Properties

The surface is friendly to sterilisation tape. Naturally low linting, it helps maintain the integrity of sterile packs.

Made with Unique SSMMS Fabric with Proven Barrier Protection

Gemini's strong multiple layers and 100% polypropylene (EN 11607-1 and EN 868-2 tested) construction provide unique high-level barrier protection against fluid and particles.

Versatile

Medline's Gemini Wrap is indicated for use in all major sterilisation cycles, including:

GEMO-GEM5

- > Prevacuum steam cycles
- > Gravity steam cycles

GEM1-GEM5

- > Ethylene oxide (EtO) sterilisation
- > STERRAD® sterilisation
- > STERIS V-PRO® low temperature sterilisation
- > STERIZONE® VP4 sterilisation.

WRAP Certification

Gemini wraps are manufactured in a WRAP-certified facility. Promoting and respecting labour and human rights in our supply chain is part of Medline's corporate social responsibility agenda. More detailed descriptions of WRAP's 12 Principles are available on their website: www.wrapcompliance.org.



GEMINI SINGLE STERILISATION WRAP

Gemini Single Sterilisation Wrap is a great solution for sterilisation units where two sheets are used to wrap a medical device or an instrument set.



	GEMO 40 gr	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr	GEM5 75 gr
46 x 46 cm	GEMO118-EU 1,000 ea/cs	GEM1118-EU 1,000 ea/cs	GEM2118-EU 1,000 ea/cs	GEM3118-EU 600 ea/cs		
51 x 51 cm	GEMO120-EU 1,000 ea/cs	GEM1120-EU 1,000 ea/cs	GEM2120-EU 500 ea/cs	GEM3120-EU 250 ea/cs		
61 x 61 cm	GEMO124-EU 500 ea/cs	GEM1124-EU 500 ea/cs	GEM2124-EU 500 ea/cs	GEM3124-EU 250 ea/cs	GEM4124-EU 250 ea/cs	
76 x 76 cm	GEMO130-EU 300 ea/cs	GEM1130-EU 300 ea/cs	GEM2130-EU 300 ea/cs	GEM3130-EU 200 ea/cs	GEM4130-EU 250 ea/cs	
91 x 91 cm	GEMO136-EU 300 ea/cs	GEM1136-EU 300 ea/cs	GEM2136-EU 300 ea/cs	GEM3136-EU 150 ea/cs	GEM4136-EU 150 ea/cs	GEM5136-EU 144 ea/cs
102 x 102 cm	GEMO140-EU 250 ea/cs	GEM1140-EU 250 ea/cs	GEM2140-EU 250 ea/cs	GEM3140-EU 150 ea/cs	GEM4140-EU 150 ea/cs	GEM5140-EU 96 ea/cs
114 x 114 cm	GEMO145-EU 250 ea/cs	GEM1145-EU 250 ea/cs		GEM3145-EU 100 ea/cs	GEM4145-EU 100 ea/cs	GEM5145-EU 96 ea/cs
122 x 122 cm	GEMO148-EU 250 ea/cs	GEM1148-EU 250 ea/cs	GEM2148-EU 100 ea/cs	GEM3148-EU 50 ea/cs	GEM4148-EU 50 ea/cs	GEM5148-EU 48 ea/cs
137 x 137 cm	GEMO154-EU 100 ea/cs	GEM1154-EU 100 ea/cs	GEM2154-EU 100 ea/cs	GEM3154-EU 50 ea/cs	GEM4154-EU 50 ea/cs	GEM5154-EU 48 ea/cs
137 x 183 cm	GEMO172-EU 100 ea/cs			GEM3172-EU 50 ea/cs	GEM4172-EU 50 ea/cs	

GEMINI INTERLEAVED STERILISATION WRAP

Gemini Interleaved Sheets offer both our multilayer technology and colour differentiation to sterilisation units where sequential wrapping is in place. This solution combines tear resistance and security. It is indicated for use in all major sterilisation cycles.



	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr
61 x 61 cm	GEM1124INT-EU 500 ea/cs	GEM2124INT-EU 500 ea/cs	GEM3124INT-EU 250 ea/cs	GEM4124INT-EU 250 ea/cs
76 x 76 cm	GEM1130INT-EU 300 ea/cs	GEM2130INT-EU 300 ea/cs	GEM3130INT-EU 200 ea/cs	GEM4130INT-EU 250 ea/cs
91 x 91 cm	GEM1136INT-EU 300 ea/cs	GEM2136INT-EU 300 ea/cs	GEM3136INT-EU 150 ea/cs	GEM4136INT-EU 150 ea/cs
102 x 102 cm	GEM1140INT-EU 250 ea/cs	GEM2140INT-EU 250 ea/cs	GEM3140INT-EU 150 ea/cs	GEM4140INT-EU 150 ea/cs
122 x 122 cm	GEM1148INT-EU 250 ea/cs	GEM2148INT-EU 100 ea/cs	GEM3148INT-EU 50 ea/cs	GEM4148INT-EU 50 ea/cs
137 x 137 cm	GEM1154INT-EU 100 ea/cs	GEM2154INT-EU 100 ea/cs	GEM3154INT-EU 50 ea/cs	GEM4154INT-EU 50 ea/cs
137 x 183 cm			GEM3172INT-EU 50 ea/cs	GEM4172INT-EU 50 ea/cs

GEMINI BONDED STERILISATION WRAP

Gemini Bonded Wrap can give you greater confidence in your sterilisation process. Two sheets fused together provide exceptional strength and help improve efficiency. The dual colour facilitates hole detection and helps differentiate between instrument sets.



	GEMO 40 gr	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr	GEM5 75 gr
46 x 46 cm	GEM0118T-EU 500 ea/cs	GEM1118T-EU 500 ea/cs	GEM2118T-EU 500 ea/cs	GEM3118T-EU 300 ea/cs		
51 x 51 cm	GEM0120T-EU 500 ea/cs	GEM1120T-EU 500 ea/cs	GEM2120T-EU 250 ea/cs	GEM3120T-EU 125 ea/cs		
61 x 61 cm	GEM0124T-EU 250 ea/cs	GEM1124T-EU 250 ea/cs	GEM2124T-EU 250 ea/cs	GEM3124T-EU 100 ea/cs	GEM4124T-EU 100 ea/cs	
76 x 76 cm	GEM0130T-EU 150 ea/cs	GEM1130T-EU 150 ea/cs	GEM2130T-EU 150 ea/cs	GEM3130T-EU 75 ea/cs	GEM4130T-EU 100 ea/cs	
91 x 91 cm	GEM0136T-EU 150 ea/cs	GEM1136T-EU 150 ea/cs	GEM2136T-EU 150 ea/cs	GEM3136T-EU 75 ea/cs	GEM4136T-EU 75 ea/cs	GEM5136T-EU 72 ea/cs
102 x 102 cm	GEM0140T-EU 100 ea/cs	GEM1140T-EU 100 ea/cs	GEM2140T-EU 100 ea/cs	GEM3140T-EU 75 ea/cs	GEM4140T-EU 75 ea/cs	GEM5140T-EU 48 ea/cs
114 x 114 cm	GEM0145T-EU 100 ea/cs	GEM1145T-EU 100 ea/cs		GEM3145T-EU 50 ea/cs	GEM4145T-EU 50 ea/cs	GEM5145T-EU 48 ea/cs
122 x 122 cm	GEM0148T-EU 100 ea/cs	GEM1148T-EU 100 ea/cs	GEM2148T-EU 50 ea/cs	GEM3148T-EU 30 ea/cs	GEM4148T-EU 30 ea/cs	GEM5148T-EU 24 ea/cs
137 x 137 cm	GEM0154T-EU 50 ea/cs	GEM1154T-EU 50 ea/cs	GEM2154T-EU 50 ea/cs	GEM3154T-EU 30 ea/cs	GEM4154T-EU 30 ea/cs	GEM5154T-EU 24 ea/cs
137 x 183 cm				GEM3172T-EU 30 ea/cs	GEM4172T-EU 30 ea/cs	GEM5172T-EU 24 ea/cs

WRAP ACCESSORIES

TRAY LINER SHEET

Paper tray liners facilitate the absorption of condensation during sterilisation.



Code	Dimensions	Pack Factor
TLC009	30 x 30 cm	2,000
TLC010	30 x 45 cm	1,000
TLC011	30 x 60 cm	1,000

DUST COVERS

These covers protect instrument trays from dust and environmental agents. They are constructed from strong, quality polyethylene and are self-sealable.

Code	Dimensions	Pack Factor
MDS100815	20 x 30 cm	1,000
MDS100825	30 x 50 cm	500
MDS100830	40 x 55 cm	500
MDS100835	40 x 75 cm	250
MDS100840	60 x 75 cm	250
MDS100845	63 x 94 cm	100



To learn more about Gemini wraps, please contact your Medline account manager or visit our site: www.medline.eu/uk



**ALWAYS
ON.**

Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk
uk-customerservice@medline.com

Gemini wraps are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

Tray liner sheets are class I non-sterile medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct any errors that may occur within this brochure.

© 2019 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc., USA. ML678-EN03/MC 10/2019.



Medline International France SAS
5, rue Charles Lindbergh
44110 Châteaubriant, France



Climpak Ltd
Beech House, Knaves Beech Business Centre, Loudwater
High Wycombe, Bucks, HP10 9SD, UK