

Gemini Sterilisation Wrap

Ensure the integrity of your sterilisation process.



**A gem of
a wrap**



Gemini Sterilisation Wrap

Helps ensure the integrity of your sterilisation process

Drapeable and stronger

Gemini wrap is made to be resistant to tears and punctures, performing better than standard SMS wraps. It is very soft and easy to fold, which makes it simple to work with.

Low linting and better adhesion properties

The surface is designed to be friendly to sterilisation tape. Naturally low linting, it helps maintain the integrity of sterile packs.

Made with unique SSMMS fabric with proven barrier protection

Gemini's strong multiple layers and 100% polypropylene (EN 11607-1 and EN 868-2 tested) construction provide unique high-level barrier protection against fluid and particles.

Versatile

Medline's Gemini Wrap is indicated for use in all major sterilisation cycles, including:

GEM0-GEM5

- Prevacuum steam cycles
- Gravity steam cycles

GEM1-GEM5

- Ethylene oxide (EtO) sterilisation
- STERRAD® sterilisation
- STERIS V-PRO® low-temperature sterilisation
- STERIZONE® VP4 sterilisation.

WRAP Certification

Gemini wraps are manufactured in a WRAP-certified facility. Promoting and respecting labour and human rights in our supply chain is part of Medline's corporate social responsibility agenda. More detailed descriptions of WRAP's 12 Principles are available on their website: www.wrapcompliance.org.



Gemini Single Sterilisation Wrap

Gemini Single Sterilisation Wrap is tailored to be an ideal solution for use in sterilisation units where two sheets are used to wrap a medical device or an instrument set.




	GEM0 40 gr	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr	GEM5 75 gr
46 x 46 cm	GEM0118-EU 1,000 ea/cs	GEM1118-EU 1,000 ea/cs	GEM2118-EU 1,000 ea/cs	GEM3118-EU 600 ea/cs		
61 x 61 cm	GEM0124-EU 500 ea/cs	GEM1124-EU 500 ea/cs	GEM2124-EU 500 ea/cs	GEM3124-EU 250 ea/cs	GEM4124-EU 250 ea/cs	
76 x 76 cm	GEM0130-EU 300 ea/cs	GEM1130-EU 300 ea/cs	GEM2130-EU 300 ea/cs	GEM3130-EU 200 ea/cs	GEM4130-EU 250 ea/cs	
91 x 91 cm	GEM0136-EU 300 ea/c	GEM1136-EU 300 ea/cs	GEM2136-EU 300 ea/cs	GEM3136-EU 150 ea/cs	GEM4136-EU 150 ea/cs	
102 x 102 cm	GEM0140-EU 250 ea/cs	GEM1140-EU 250 ea/cs	GEM2140-EU 250 ea/cs	GEM3140-EU 150 ea/cs	GEM4140-EU 150 ea/cs	GEM5140-EU 96 ea/cs
114 x 114 cm					GEM4145-EU 100 ea/cs	
122 x 122 cm	GEM0148-EU 250 ea/c	GEM1148-EU 250 ea/cs	GEM2148-EU 100 ea/cs	GEM3148-EU 50 ea/cs	GEM4148-EU 50 ea/cs	GEM5148-EU 48 ea/cs
137 x 137 cm	GEM0154-EU 100 ea/cs	GEM1154-EU 100 ea/cs	GEM2154-EU 100 ea/cs	GEM3154-EU 50 ea/cs	GEM4154-EU 50 ea/cs	
137 x 183 cm	GEM0172-EU 100 ea/cs					

Gemini Interleaved Sterilisation Wrap

Gemini Interleaved Sheets offer both our multilayer technology and colour differentiation to sterilisation units where sequential wrapping is in place. This solution combines tear resistance and security. It is indicated for use in all major sterilisation cycles.



	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr
61 x 61 cm	GEM1124INT-EU 500 ea/cs	GEM2124INT-EU 500 ea/cs	GEM3124INT-EU 250 ea/cs	GEM4124INT-EU 250 ea/cs
76 x 76 cm	GEM1130INT-EU 300 ea/cs	GEM2130INT-EU 300 ea/cs	GEM3130INT-EU 200 ea/cs	GEM4130INT-EU 250 ea/cs
91 x 91 cm	GEM1136INT-EU 300 ea/cs	GEM2136INT-EU 300 ea/cs	GEM3136INT-EU 150 ea/cs	GEM4136INT-EU 150 ea/cs
102 x 102 cm	GEM1140INT-EU 250 ea/cs	GEM2140INT-EU 250 ea/cs	GEM3140INT-EU 150 ea/cs	GEM4140INT-EU 150 ea/cs
122 x 122 cm	GEM1148INT-EU 250 ea/cs	GEM2148INT-EU 100 ea/cs	GEM3148INT-EU 50 ea/cs	GEM4148INT-EU 50 ea/cs
137 x 137 cm			GEM3154INT-EU 50 ea/cs	GEM4154INT-EU 50 ea/cs
137 x 183 cm				GEM4172INT-EU 50 ea/cs



Gemini Bonded Sterilisation Wrap

Gemini Bonded Wrap can give you greater confidence in your sterilisation process. Two sheets fused together provide exceptional strength and help improve efficiency. The dual colour facilitates hole detection and helps differentiate between instrument sets.




	GEM0 40 gr	GEM1 47 gr	GEM2 54 gr	GEM3 61 gr	GEM4 71 gr	GEM5 75 gr
46 x 46 cm	GEM0118T-EU 500 ea/cs	GEM1118T-EU 500 ea/cs	GEM2118T-EU 500 ea/cs	GEM3118T-EU 300 ea/cs		
61 x 61 cm	GEM0124T-EU 250 ea/cs	GEM1124T-EU 250 ea/cs	GEM2124T-EU 250 ea/cs	GEM3124T-EU 100 ea/cs	GEM4124T-EU 100 ea/cs	
76 x 76 cm	GEM0130T-EU 150 ea/cs	GEM1130T-EU 150 ea/cs	GEM2130T-EU 150 ea/cs	GEM3130T-EU 75 ea/cs	GEM4130T-EU 100 ea/cs	
91 x 91 cm	GEM0136T-EU 150 ea/cs	GEM1136T-EU 150 ea/cs	GEM2136T-EU 150 ea/cs	GEM3136T-EU 75 ea/cs	GEM4136T-EU 75 ea/cs	GEM5136T-EU 72 ea/cs
102 x 102 cm	GEM0140T-EU 100 ea/cs	GEM1140T-EU 100 ea/cs	GEM2140T-EU 100 ea/cs	GEM3140T-EU 75 ea/cs	GEM4140T-EU 75 ea/cs	GEM5140T-EU 48 ea/cs
114 x 114 cm		GEM1145T-EU 100 ea/cs		GEM3145T-EU 50 ea/cs	GEM4145T-EU 50 ea/cs	GEM5145T-EU 48 ea/cs
122 x 122 cm	GEM0148T-EU 100 ea/cs	GEM1148T-EU 100 ea/cs	GEM2148T-EU 50 ea/cs	GEM3148T-EU 30 ea/cs	GEM4148T-EU 30 ea/cs	GEM5148T-EU 24 ea/cs
137 x 137 cm	GEM0154T-EU 50 ea/cs	GEM1154T-EU 50 ea/cs	GEM2154T-EU 50 ea/cs	GEM3154T-EU 30 ea/cs	GEM4154T-EU 30 ea/cs	GEM5154T-EU 24 ea/cs
137 x 183 cm				GEM3172T-EU 30 ea/cs	GEM4172T-EU 30 ea/cs	

Strength Table

Support the weight of every wrapped set by choosing the correct GEM series.

Medline GEM Series	Recommended Tray Weight Limit per Medline's IFU
Super lightweight GEM0XXX	0 – 2.7 kg
Lightweight GEM1XXX	2.7 kg
Regular weight GEM2XXX	4 kg
Medium weight GEM3XXX	5.9 kg
Heavyweight GEM4XXX	7.7 kg
Super heavyweight GEM5XXX	11.3 kg

The weight limits are in place for the pre-vacuum, gravity steam and EO sterilisation modalities. For other modalities, please refer to the IFU.



Wrap Accessories

Dust Cover

These covers are designed to protect instrument trays from dust and environmental agents. They are constructed from strong, quality polyethylene and are self-sealable.

Item number	Description	Packaging
MDS100815	20 x 30 cm	1,000
MDS100825	30 x 50 cm	500
MDS100830	40 x 55 cm	500
MDS100835	40 x 75 cm	250
MDS100840	60 x 75 cm	250
MDS100845	63 x 94 cm	100



Tray Liner Sheet

Paper tray liners facilitate the absorption of condensation during sterilisation.

Item number	Description	Packaging
TLC009	30 x 30 cm	2,000
TLC010	30 x 45 cm	1,000
TLC011	30 x 60 cm	1,000



Transport Tray

- Rigid plastic is designed to permit air flow
- Helps reduce exposed wire shelving
- Helps protect instrument tray from holes or tears during sterilisation, storage and transport


Item number	Description	Packaging
GEMTRAY3CE	31 x 56 cm	10
GEMTRAY4CE	43 x 69 cm	8



Aluminium Sterilisation Plate

The Aluminium Sterilisation Plate offers a reusable solid aluminium sheet, intended for use under sterilisation trays or baskets, to create a heat sink, leading to improved evaporation of steam to ward off water during steam sterilisation.

Item number	Description	Packaging
LCLC1850CE	52 x 24 cm	1
LCLC1851CE	48 x 24 cm	1




Silicone Tubing Roll

This non-sterile silicone tubing roll can be cut and used to protect the sharp tips of surgical instruments in storage, handling and sterilisation.

Item number	Inner diameter (ID)	Outer diameter (OD)	Packaging
510204	2	4	16
510305	3	5	16
510406	4	6	16
510407	4	7	10
510508	5	8	7
510509	5	9	7
510609	6	9	4
510612	6	12	4
510710	7	10	4
510812	8	12	4
511014	10	14	4

One 15-m roll is easy to cut to the required length, depending on clinical needs



To learn more about Gemini wraps, please contact your Medline account manager or visit our site: www.medline.eu/uk



Medline Industries Ltd.
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237

www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd.
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

Gemini wraps and transport trays are class I non-sterile medical devices intended to be used by healthcare professionals.



Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France



Tray liner sheets are class I non-sterile medical devices intended to be used by healthcare professionals.



Medical Device Safety Service (MDSS)
Schiffgraben 41
30175 Hannover, Germany



Aluminium sterilisation plates are class I non-sterile medical devices intended to be used by healthcare professionals.



PROSIMED S.r.l.
Via della Magliana, 295
00146 Roma, Italia



Silicone tubing rolls are class I non-sterile medical devices intended to be used by healthcare professionals.



Shanghai International Holding Corp. GmbH
Eiffelstraße 80
20537 Hamburg, Germany



Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur in this brochure.

© 2024 Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML678_EN 02/2024