

Disposable Headwear



**Comfort
and variety**



Headwear

Covering your head and hair helps reduce the risk of foreign particles getting into surgical wounds and protects patients from potential contamination in operating theatres or ICUs. Medline's disposable headwear products offer a high level of comfort and come in a variety of shapes, sizes and designs that address hospitals' protective and aesthetical needs.

Our broad product range, including surgical hoods and caps, bonnets and bouffant caps, is suitable for a variety of operations and hairstyles. To offer you even more choice and meet your budget, Medline's headwear is available in cluster keyback, viscose and spunbond polypropylene materials that have passed biocompatibility requirements, are low linting and do not contain latex.

Comfort Line

The *Comfort Line* is Medline's preferred headwear product range for the operating theatre. Made of soft and breathable cluster keyback fabric, this line offers both a high level of comfort and a high level of protection. Comfort Line products are available in blue and green and a large variety of designs. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why our Comfort product range is produced in a WRAP-certified facility.



Surgical hood

FS61100

100/box, 600/case

Surgical hood, green



Surgical hood

FS61150A

100/box, 600/case

Surgical hood, blue



Surgical hood

FS61250

50/box, 300/case

Surgical hood with sweatband, blue



Bouffant cap

FS61300

100/box, 600/case

Bouffant cap with a latex-free elastic band, green, size 53 cm



Bouffant cap

FS61350

100/box, 600/case

Bouffant cap with a latex-free elastic band, blue, size 53 cm



Surgical cap

FS61450

100/box, 600/case

Surgical cap, sailor style, blue



Surgical cap

FS61500

100/box, 600/case

Surgical cap, elastic at the back, contoured front, green



Surgical cap

FS61550A

100/box, 600/case

Surgical cap, elastic at the back, contoured front, blue



Surgical cap

FS61750A

100/box, 600/case

Surgical cap with ties, for long hair, blue



Surgical cap

FS61850A

100/box, 600/case

Surgical cap with ties, for short hair, blue



Surgical cap


FS61950

100/box, 600/case

Surgical cap, elastic at the back, flat front, blue

Colourful Line

The *Colourful Line* is made of spunlace viscose, which is a soft, breathable and textile-like material. The range comes in a variety of colours and offers good comfort. Medline is committed to sourcing products from ethical and socially responsible factories around the world. That's why most of our Colourful range is produced in a WRAP-certified facility.




Over-the-shoulder hood

FS60000

50/box, 300/case

Over-the-shoulder hood with ties and sweatband, blue




Surgical hood

FS60100G

100/box, 600/case

Surgical hood, green




Surgical hood

FS60150B

100/box, 600/case

Surgical hood, blue



Surgical hood

FS60250B

50/box, 300/case

Surgical hood with sweatband, blue



Surgical cap

FS60750B

100/box, 600/case

Surgical cap with ties, for long hair, blue



Surgical cap

FS60850B

100/box, 600/case

Surgical cap with ties, for short hair, blue



Surgical cap

FS60950B

100/box, 600/case

Surgical cap, elastic at the back, flat front, blue



Surgical cap

FS60950G

100/box, 600/case

Surgical cap, elastic at the back, flat front, green



Surgical cap

FS60950P

100/box, 600/case

Surgical cap, elastic at the back, flat front, purple



Surgical cap

FS60950PK

100/box, 600/case

Surgical cap, elastic at the back, flat front, pink



Surgical cap

FS60950W

100/box, 600/case

Surgical cap, elastic at the back, flat front, white



Surgical cap

FS60950Y

100/box, 600/case

Surgical cap, elastic at the back, flat front, yellow

Basic Line

Designed to satisfy the essential needs of the ward at a convenient price, the *Basic Line* is made of a spunbond material with a basis weight from 12 gsm to 30 gsm. *Basic Line* headwear comes in blue, green, white or printed patterns.



Surgical cap

NONE28625

100/box, 500/case

Surgical cap with ties, SMS top, absorbent scrim-reinforced material on sides, one size fits most, light blue



Surgical cap

FS67050

100/box, 600/case

Spunlace bonnet, one size fits most, blue



Surgical hood

FS63100A

100/box, 600/case

Surgical hood, green



Surgical hood

FS63150A

100/box, 600/case

Light surgical hood, blue

**Bouffant cap****FS63300**

100/box, 600/case

Bouffant cap, size 60 cm, green

**Bouffant cap****FS63350**

100/box, 600/case

Bouffant cap, size 60 cm, blue

**Bouffant cap****FS63360**

100/box, 600/case

Bouffant cap, size 60 cm, white

**Surgical cap****FS63450A**

100/box, 600/case

Surgical cap, sailor style, blue

**Clip cap****FS65460**

100/box, 1,200/case

Clip cap, one size fits most, white

**Clip cap****FS65450**

100/box, 1,200/case

Clip cap, one size fits most, blue

**Bouffant cap****CRIE1000**

100/box, 500/case

Bouffant cap, size 53 cm, green

**Bouffant cap****CRIE1002**

100/box, 500/case

Bouffant cap, size 53 cm, white


**Bouffant cap****CRIE1004**

100/box, 500/case

Bouffant cap, size 60 cm, blue

Headwear accessory

Individual sweatband



Medline's individual sweatband can be placed under any surgical cap and provides comfort to healthcare professionals during long procedures thanks to its soft viscose material.

Item number	Description	Packaging
BAND250*	Sweatband with elastic, 18 x 4 cm	50/box, 400/case

Choose your headwear

Surgical hoods



Offering a helmet-like cap with ties at the neck, one-size-fits-all surgical hoods provide the most coverage of all products in our headwear range for worry-free performance during surgery.

Some products feature an absorbent sweatband for increased staff comfort, especially for longer procedures.

Item number	Description	Packaging
FS61100	Cluster keyback, green	100/box, 600/case
FS61150A	Cluster keyback, blue	100/box, 600/case
FS61250	Cluster keyback, with sweatband, blue	50/box, 300/case
FS60100G	Viscose, green	100/box, 600/case
FS60150B	Viscose, blue	100/box, 600/case
FS60250B	Viscose, with sweatband, blue	50/box, 300/case
FS63100A	Spunbond polypropylene 30 gsm, green	100/box, 600/case
FS63150A	Spunbond polypropylene 30 gsm, blue	100/box, 600/case

Surgical caps with ties for long and short hair



Featuring ties at the back of the neck, these adjustable surgical caps are designed for people with long or short hair to safely keep the hair covered.

This surgical cap is ideal for use in the operating theatre or for day surgeries.

Item number	Description	Packaging
For long hair		
FS61750A	Cluster keyback, blue	100/box, 600/case
FS60750B	Viscose, blue	100/box, 600/case
For short hair		
FS61850A	Cluster keyback, blue	100/box, 600/case
FS60850B	Viscose, blue	100/box, 600/case
NONE28625	SMS and scrim-reinforced material, blue	100/box, 500/case

Surgical caps with elastic at the back



Featuring a broad shape and elastic at the back, these one-size-fits-most surgical caps provide a comfortable and complete hair coverage, meeting the requirements of the most active of days.

This cap design is well-suited for long hairstyles and long periods of wear in the operating theatre, ICU and other laboratory departments.

Item number	Description	Packaging
FS61500	Cluster keyback, green	100/box, 600/case
FS61550A	Cluster keyback, blue	100/box, 600/case
FS61950	Cluster keyback, blue	100/box, 600/case
FS60950B	Viscose, blue	100/box, 600/case
FS60950G	Viscose, green	100/box, 600/case
FS60950P	Viscose, purple	100/box, 600/case
FS60950PK	Viscose, pink	100/box, 600/case
FS60950W	Viscose, white	100/box, 600/case
FS60950Y	Viscose, yellow	100/box, 600/case

Surgical caps without ties or elastic



These surgical caps do not contain ties or elastic, which makes them comfortable for longer periods of wear thanks to their lack of pressure areas. The caps are designed to fit well around the head due to the natural flexibility of the material.

The caps are ideal for nurses, visitors in the operating theatre and ward staff with short hair. Additionally, the edge of the bonnet (FS67050) can be tuned up to function as a sweatband for even more comfort.

Item number	Description	Packaging
FS61450	Cluster keyback, blue	100/box, 600/case
FS63450A	Spunbond polypropylene 30 gsm, blue	100/box, 600/case
FS67050	Spunlace 35 gsm, blue	100/box, 600/case

Bouffant caps



These light and roomy caps with a round shape feature elastic all the way around, ensuring all hair is fully and securely covered and that the cap fits comfortably.

Bouffant caps are suitable for long and short hairstyles and provide comfort and breathability for staff in less critical areas in the hospital.

Very similar to the design of the bouffant cap, the compactly packed clip cap is ideal headwear for patients and visitors in less critical areas of the hospital.

Item number	Description	Packaging
FS61300	Cluster keyback, 53 cm, green	100/box, 600/case
FS61350	Cluster keyback, 53 cm, blue	100/box, 600/case
FS63300	SPP 20 gsm, 60 cm, green	100/box, 600/case
FS63350	SPP 20 gsm, 60 cm, blue	100/box, 600/case
FS63360	SPP 20 gsm, 60 cm, white	100/box, 600/case
CR1E1000	SPP 12 gsm, 53 cm, green	100/box, 500/case
CR1E1002	SPP 12 gsm, 53 cm, white	100/box, 500/case
CR1E1004	SPP 12 gsm, 60 cm, blue	100/box, 500/case
FS65450	SPP 14 gsm, one size fits most, blue	100/box, 1,200/case
FS65460	SPP 14 gsm, one size fits most, white	100/box, 1,200/case

SPP = spunbond polypropylene



Internationally
recognised social
compliance certification
body that verifies
whether factories
respect workers' labour
and human rights

**Medline's Comfort
and Colourful
product ranges
are produced in a
WRAP-certified
facility.**

For more information, contact your Medline
account manager or visit us at: www.medline.eu/uk.



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237

www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd
1 Grant's Row
Mount Street Lower
Dublin 2
Ireland
Tel.: +353 1 691 73 06

www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

These products are class I non-sterile medical devices intended to be used by healthcare professionals.
Before use, please consult instructions and precautions on the corresponding labelling.

*BAND250 is a commodity.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2022 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093, USA. ML832_EN 12/2022

