




CENTURION
MEDICAL PRODUCTS

HubGuard®
Catheter securement

**Protects skin
and provides
comfort**




Benefits for caregivers, hospitals and patients

- Adhesive-backed foam securely anchors the catheter at the infusion site
- Enhances stability to reduce 'pistoning' of the catheter in the vein, reducing catheter migration, dislodgement and unscheduled IV restarts
- Cushion protects the skin and provides comfort from rigid catheter hubs
- Non-irritating adhesive is gentle on skin, yet stays in place for 72 hours or until clinically indicated, per hospital protocol
- Available in 2 materials that do not contain natural rubber latex:
 - For maximum securement and cushioning, it integrates 1.6 mm thick PE foam

Available sizes

Item code	Size	Material
HGL4	1.9 x 8.6 cm	1.6 mm thick PE foam
HGL6	1.6 x 6.0 cm 0.95 x 2.2 cm	
HGL7	1.6 x 4.0 cm 0.95 x 1.1 cm	
HGL8	1.6 x 5.0 cm 0.95 x 2.2 cm	
HGDBSXT	1.6 x 4.0 cm 0.95 x 2.2 cm	
HG2	0.95 x 3.8 cm	
HG05	1.6 x 6.0 cm	



Other sizes and styles may be available. Contact your Medline account manager to learn about additional catheter securement products.



HGL4



Medline Industries Ltd
 3rd Floor
 Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington WA4 6HL
 United Kingdom
 Tel.: +44 844 334 5237
www.medline.eu/uk
uk-customerservice@medline.com

Medline Ireland Ltd
 1 Grant's Row
 Mount Street Lower
 Dublin 2
 Ireland
 Tel.: +353 1 691 73 06
www.medline.eu/ie
ie-customerservice@medline.com

FOLLOW US  

Securement devices are class I sterile medical devices intended to be used by healthcare professionals. Before use, please consult instructions and precautions on the corresponding labelling.



Medline International France SAS
 5 Rue Charles Lindbergh,
 44110 Châteaubriant, France



All rights reserved. We reserve the right to correct errors that may have occurred.

© 2022 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. ML946_EN 04/2022
 © 2022 Centurion Medical Products Corporation LIT026V7