



**PROTECTION IN THE
TOUGHEST SITUATION.**



INTENSIVE CARE SOLUTIONS
Comprehensive Range of Products for Patients and Staff

FROM PATIENT CARE TO STAFF PROTECTION, WE'RE HERE TO HELP.

Every day, hospitals look for innovative solutions to help care for patients and protect staff in an efficient manner and with the highest quality products.

With patient care solutions (such as patient hygiene, moisture management, pressure ulcer prevention and fluid management) and staff protection (like protective apparel and exam gloves), Medline has what you need for your intensive care unit (ICU).

Medline products are designed to support your standard of practice, standardise care and improve patient outcomes while protecting your staff.

TABLE OF CONTENTS

Patient Hygiene	3
Moisture Management	4
Pressure Ulcer Prevention	4
Suction	5
Protective Apparel	6
Exam Gloves	7

PATIENT HYGIENE

EasyBath® Patient Bathing

The EasyBath Patient Bathing solution does not require additional water, rinsing or drying.

With ultra-soft and hypoallergenic materials, EasyBath provides skin protection and can be warmed for increased patient comfort.

By using a packaged solution, caregivers can give a bath in less time than with the traditional basin methods. EasyBath also assists with moisture management and, ultimately, the prevention of pressure ulcers.

Impregnated¹

Item	Description	Wipes/pk
MSCE095103A	EasyBath, cloth, fragrance free	8
MSCE095104	EasyBath+, scented	8
MSCE095105	EasyBath+, fragrance free	8
MSCE095309	EasyBath lightweight cloth, fragrance free	8
MSCE095230	EasyBath+, shampoo cap, scented	1
MSCE095231	EasyBath+, shampoo cap, fragrance free	1
MSCE263654	Aloetouch, wipe, scented	48
MSCE263810	ReadyFlush, scented	24
MSCE995300WM	EasyBath, wash gloves, scented	8
MSCE995301WM	EasyBath, wash gloves, fragrance free	8
MSCE995302WM	EasyBath, wash gloves, scented	5
MSCE995310WM	EasyBath, wash gloves, scented, antibacterial	8
MSCE995311WM	EasyBath, wash gloves, fragrance, antibacterial	8

Non-impregnated

Item	Description	Wipes/pk
MEDGLOVE01	Dry wash gloves, 16 x 23 cm, 75 g	100
MEDGLOVE02	Dry wash gloves, 14.5 x 20 cm, 65 g	100
MEDGLOVE03	Dry wash gloves, 14.5 x 20 cm, 75 g	100

VAPprevent® Oral Care²

For ventilated patients in the ICU, VAPprevent products offer a convenient, all-in-one solution for cleaning patients' oral cavities and reducing the risk factors correlated with ventilator-associated pneumonia (VAP). These tools and components are designed for ease of use and help provide consistent oral care to all patients.

Item	Description
MDS096013MEU	VAPprevent standard kit, 2 swabs, rinse and gel
MDS096202EU	DenTips, oral swab, blue untreated, including wrap
MDS096502EU	DenTips, oral swab, green treated, including wrap
MDS09651300EU	VAPprevent suction kit, 1 swab, 2 swabs
MDS096513MEU	VAPprevent suction kit, 1 swab, 2 swabs, rinse, gel
MDS096525EU	Suction swab, treated, including wrap
MDS09657200EU	VAPprevent toothbrush kit, 1 swab
MDS096572MEU	VAPprevent toothbrush kit, 1 swab, rinse, gel
MDS096575EU	Suction toothbrush, treated, including wrap



Dry wash gloves used with cosmetic products of choice.



MOISTURE MANAGEMENT

Ultrasorbs^{®3}

Ultrasorbs are unique, breathable drypads developed to manage moisture and help prevent moisture-associated skin damage (MASD) in the most challenging care settings.

Whether it is ongoing fluid loss from wound exudates or drains, oedematous legs, incontinence or increased perspiration, Ultrasorbs will absorb and lock away any bodily moisture from patients' skin, keeping it dry, while also protecting the surface below.

Ultrasorbs are particularly useful in high-dependency wards like the ICU where caring for patients is difficult due to their fragile condition.



Item	Size	Packaging
ULTRASRBE2540	24 x 40 cm	100/cs, 10/bg
ULTRASRBE4561	45 x 61 cm	60/cs, 10/bg
ULTRASRBE6191	61 x 91 cm	70/cs, 10/bg
MSCE4090PUPS	101.6 x 228.6 cm	30/cs
MSCE4090PUPSB	101.6 x 228.6 cm	1/pk, 30/cs
ULTRSRB4561ST sterile	45 x 61 cm	1/pk, 160/cs
ULTRSRB6191STA sterile	61 x 91 cm	1/pk, 50/cs
MSCE4090PPSTA sterile	101.6 x 228 cm	1/pk, 25/cs
USAPE6191ES	61 x 91 cm	60/cs, 10/bg
USAP4090ES	101.6 x 228 cm	1/pk, 25/cs

PRESSURE ULCER PREVENTION

Heelmedix[®] and Heel Raiser[®] Pro⁴

Heelmedix and Heel Raiser Pro are boots that help prevent and manage heel pressure ulcers for patients with limited mobility. Especially in the ICU, these boots aid healthcare workers in delivering a high quality of care by helping patients avoid heel pressure ulcers through an offloading device.

Item	Description	Packaging
MDTE823330P	Heelmedix, petite, heel protector	1/ea
MDTE823330SW	Heelmedix, with wedge, heel protector	1/ea
MDTE823330XL	Heelmedix, XL, with wedge, heel protector	1/ea
MDTE823200CS	Heel Raiser Pro, heel protector	4/cs



SUCTION

Suction Tubing, Canisters and Yankauers

Medline offers a full range of suction circuits with reusable canisters and disposable liners, tubing and yankauers.

Canisters and Liners⁵

MED-SOFT disposable flexible liners, with lids permanently attached to the liners, minimise the risk of infectious medical waste leaking and are engineered for clinician safety.

Capacity	Reusable Canisters	Disposable Liners
1 l	DYNDCL01000	DYNDSC1000
1.5 l	DYNDCL01500	DYNDSC1500
3 l	DYNDCL03000	DYNDSC3000



Tubing^{6,7}

ORNEX VC tubing⁶ is built with a vacuum control connector, which provides a good fastening point for bronchial-tracheal catheters.

Item	Type	ID	Length	Packaging
ORNEX56VC	F/F+MVC	5 mm	1.80 m	50/cs
ORNEX510VC	F/F+MVC		3 m	50/cs
ORNEX66VC	F/MVC	6 mm	1.80 m	50/cs
ORNEX610VC	F/MVC		3 m	50/cs



ORNEX Roll Tube⁷ is perfect for the ICU because you can cut the tubing to the size needed for the room or procedure.

Item	ID	Length	Packaging
ORNEX5100	5 mm	30 m roll = 1.8 m x 16 suction tubes	1/cs
ORNEX6100	6 mm		1/cs
ORNEX7100	7 mm		1/cs



Yankauers⁸

Medline offers multiple sizes of yankauers to assist in oral care procedures, available with or without vacuum control. The rounded tip is designed to minimise tissue trauma, even on the most delicate areas.

Item	Description	Length	Diameter	Holes	Vacuum Control	Packaging
Y1200	Fine tip for paediatric, dental and ENT use	275 mm	12 CH	2	No	50/cs
Y1201		275 mm	12 CH	2	Yes	50/cs
Y1500		195 mm	15 CH	0	No	100/cs
Y1501		195 mm	15 CH	0	Yes	100/cs
Y1800	Regular tip for pharyngeal use	250 mm	18 CH	6	No	50/cs
Y1801		250 mm	18 CH	6	Yes	50/cs
Y240FS	Straight tip	297 mm	24 CH	4	Yes	50/cs
Y240FB	Bulb tip	297 mm	24 CH	4	Yes	50/cs



PROTECTIVE APPAREL

Medline protective apparel helps safeguard healthcare professionals and ICU visitors and contributes to better infection control in the ICU environment. We use high-quality materials that offer great protection and comfort.

Isolation Gowns⁹



NONE386CF

The examination isolation gown's front and sleeves are made of medium-weight SMS material, and the back is made of spunbond material, for additional comfort and breathability. The examination gown is equipped with ties at the neck and waist, along with knit cuffs at the wrists. The gown is breathable, comfortable and generously sized for full coverage and flexibility.



CRIE4040C

The isolation gown's front and sleeves are made of an impervious bilaminate material and the back of breathable SPP material. This design offers a great balance between comfort and protection. The impervious isolation gown and is equipped with knit cuffs and ties at the neck and waist, which ensure that the gown remains securely in place.

Scrub Suits¹⁰



Medline disposable scrub suits are designed in collaboration with healthcare professionals and represent a premium alternative to reusable scrub suits. The Soft scrub suits range is manufactured from a comfortable and soft 35 gsm SMS material, which has increased opacity to promote confident daily activity in the ward.


The shirt has a V-neck, three pockets and a bell-shaped design for a perfect fit, and the pants are equipped with ties at the waist. The range is available in blue and green.

Blue scrub suits set, L size: P35PBL
Green scrub suit set, L size: P35PGL

Headwear¹¹



FS61750A
for long hair



FS61850A
for short hair

Headwear is commonly worn by healthcare workers and ICU visitors.

Facemasks¹²



NONE27378



NON27412EL

All Medline procedure masks are manufactured from carefully selected, high-quality materials that comply with biocompatibility and performance standards according to EN14683. All facemasks are latex-free.

EXAM GLOVES

Medline exam gloves come in a variety of materials and lengths and help prevent the spread of infection.

Restore™ Nitrile Exam Gloves¹³



Restore Exam Gloves provide optimal user protection and are dual-certified nitrile exam gloves that conform to both the Medical Device Directive (MDD) and Personal Protective Equipment directive (PPE).

Restore gloves also have an inner coating of colloidal oatmeal to promote healthy skin and drive proper hand hygiene compliance. Restore provides an ideal mix of strength and sensitivity plus resistance to a variety of chemicals and chemotherapy drugs.



Additionally, the packaging is equipped with MicroStop film, which helps reduce waste and improve efficiency when dispensing gloves.

X-Small	Small	Medium	Large	X-Large
OAT-XS	OAT-S	OAT-M	OAT-L	OAT-XL

SensiCare® Extended Cuff Exam Gloves¹³



SensiCare Extended Cuff Exam Gloves are also MDD and PPE dual-certified. SensiCare Extended Cuff gloves offer additional thickness and length for situations that require extra protection.

They are resistant to a variety of chemicals and chemotherapy drugs.

X-Small	Small	Medium	Large	X-Large
MDSEXTXS	MDSEXTS	MDSEXTM	MDSEXTL	MDSEXTXL

Medline Sterile Procedure Nitrile Exam Gloves¹⁴



Medline Sterile Procedure Nitrile Exam Gloves are sterile, have an extended cuff and come packed in pairs. They are perfect for situations where increased protection and sterility are required.

Small	Medium	Large	X-Large
MDS2294	MDS2295	MDS2296	MDS2297


PATIENT HYGIENE



MOISTURE MANAGEMENT



SUCTION



PRESSURE ULCER PREVENTION

PROTECTIVE APPAREL



EXAM GLOVES



For Medline's full range of products, please visit www.medline.eu/uk or speak with your local account manager.



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk
uk-customerservice@medline.com

ALWAYS ON.

¹EasyBath products, except MSC095309, are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

²VAPrevent products are class I and IIa non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

³Non-sterile Ultrasorbs are class I medical devices. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France. Sterile Ultrasorbs are class I medical devices intended to be used by healthcare professionals. B. Braun Medical GmbH, 38110 Helmstedt, Germany.

⁴Heelmedix and Heel Raiser Pro are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

⁵Canisters and liners are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

⁶ORNEX VC tubing is a class I sterile medical device intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

⁷ORNEX Roll Tubes are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

⁸Yankauers are class IIa sterile medical devices intended to be used by healthcare professionals. WINGSMEAD SERVICE LIMITED, 19 Mezzanine Floor, 19-21 Crawford Street, London, W1H 1PL, UK.

⁹Isolation gowns are in Personal Protective Equipment Category I according to the 89/686/EEC Directive.

¹⁰Scrub suits are in Personal Protective Equipment Category I according to the 89/686/EEC Directive.

¹¹Headwear products are class I non-sterile medical devices and are intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

¹²NON27378 and NON27412EL facemasks are class I non-sterile medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

¹³Non-sterile exam gloves are class I medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

¹⁴Sterile exam gloves are class I medical devices intended to be used by healthcare professionals. Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France.

Before use, consult instructions and precautions on the corresponding labelling.
We reserve the right to correct any errors that may occur within this brochure.
©2019 Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML491-EN02/CM 02/2019.