



# EFFECTIVELY REDUCE INFLAMMATION.



**Cold Therapy Solutions**  
Safe and Easy, Specialised Cold Therapy

# ACCU-THERM™ ICE BAGS

Accu-Therm ice bags are used in cold therapy treatments to reduce inflammation and swelling as well as to help relieve muscle spasms by lowering skin temperature.



Medline ice bags are designed with an impervious inner layer that prevents moisture from forming on the outside and they feature a spunbond outer layer to help prevent ice burn.

We offer a **wide range of products in different sizes and shapes** to meet all your needs.

Each ice bag is designed for specialised cold therapy. Reusable closures make the products easy to refill for repeated use on a single patient.

				
<b>Tie-on bag</b>	<b>Hook-and-loop bag</b>	<b>EENT bag</b>	<b>Perineal bag</b>	<b>Maxillo-facial bag</b>
A universal ice bag in medium and large that can be attached directly to the affected area with 4 ties	A universal ice bag with 2 hook-and-loop straps	Ice bag designed for patients recovering from eye or nasal surgery	Anatomically shaped for great comfort to help women recover after delivery	Specially designed for cold therapy on the face and head; secured with hook-and-loop strap

Item	Description	Size	Packaging
NON4410E NON4420E	Tie-on bag	35.7 x 13.3 cm 41.3 x 17.9 cm	25 ea/bx 2 bx/cs
NON4430E	Maxillo-facial bag	35.5 x 13 cm	12 ea/bx 2 bx/cs
NON4440E	EENT bag	28.5 x 19.9 cm	15 ea/bx 2 bx/cs
NON4460E	Hook-and-loop bag	41.3 x 17.9 cm	15 ea/bx 2 bx/cs
NON4480E	Perineal bag	35.7 x 11 cm	50 ea/bx 2 bx/cs

Accu-Therm ice bags are medical devices in class I non sterile intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



# INSTANT COLD PACKS

- » Soft-touch non-woven (fabric-like) lining helps prevent the pack from sweating and it can be placed directly on the skin.
- » Economic solutions that eliminate the need for freezing and for holders or towels.
- » Non-toxic urea-blend-based formula
- » Fold-to-activate design enables easier activation.
- » Products are available in two sizes.

Item	Description	Size	Target Temperature	Packaging
CPACKS	Instant cold pack - small size	12.7 x 19 cm (5 x 7.5 in)	0°C/32°F	24/cs
CPACKL	Instant cold pack - large size	16 x 24 cm (6.25 x 9.5 in)	0°C/32°F	24/cs




# PERINEAL INSTANT COLD PACKS

- » Perineal instant cold pack is a solution that combines post-delivery absorbency with soothing cool relief.
- » Economic solution that eliminates the need for freezing and for pads and packs.
- » Adhesive strip for secure placement
- » Non-toxic urea-blend-based formula
- » Fold-to-activate design enables easier activation.
- » Universal size

Item	Description	Size	Absorbency	Adhesive	Packaging
CPACKP	Perineal instant cold pack	11.4 x 36.2 cm (4.5 x 14.25 in)	Heavy	Strip	24/cs



Instant cold packs are medical devices in class IIa intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



Medline offers a great range of different products and services for you.



**Prolance<sup>®</sup>**  
Safety lancets for capillary blood sampling.

Prolance safety lancets are medical devices in class IIa intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



HTL Strefa S.A ul. Adamówek 7,  
95-035, Ozorkow, Poland

CE0344 (DEKRA Certification)



**Disposable Tourniquets**  
Cost Effective Alternative to Reusable Tourniquets

Tourniquets are medical devices in class I non sterile intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



Medline International France SAS,  
5 rue Charles Lindbergh,  
44110 Châteaubriant, France

CE



**ALWAYS  
ON.**

**Medline Industries Ltd**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
Fax: +44 844 334 5238

[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)