

# Invisishield™ Drapes

Unique products to help improve safety and patient care



**Protect and  
see clearly**

# Overview

Medline's Invisishield specialty drapes are designed to make the clinician's job easier, while improving safety for both patients and staff.

## Invisishield patient drapes

Medline's Invisishield drapes are clear so that clinicians can easily see theatre equipment and the incision area. The impervious, tear-resistant material provides exceptional protection and offers superb drapeability.

## Clear, impervious and drapeable

Configurations for all major procedures and applications can be found in each section of our different Invisishield patient drapes.

- Transparent – Clinicians can easily view the area around incision sites.
- Impervious – Made of tear-resistant polyethylene (PE) and polyurethane films, Invisishield drapes provide superb protection from fluids.
- Drapeable – Invisishield drapes have just the right amount of flexibility and softness for superb drapeability.

## Useful features

Invisishield drapes have handy features that help make the clinician's job easier. Our fluid pouches feature handy tabs on the adhesive liner that make removal easy. Capped, closable ports eliminate the need to cut or pierce the pouch for drainage, improving fluid control and reducing set-up times.

# Table of contents

|  |    |
|--|----|
| General accessories .....                      | 3  |
| Orthopaedic drapes .....                       | 7  |
| Ophthalmic drapes .....                        | 9  |
| Other drapes .....                             | 12 |
| Dermadine® and Dermadrape® Incise drapes ..... | 14 |
| Equipment drapes .....                         | 15 |
| Stockinettes .....                             | 22 |
| Cross reference .....                          | 23 |

# General accessories

## Towel drapes

### Description

- Strong adhesive with no residue
- Non-linting, tear-resistant, poly-release liner
- Strong, drapeable, low-glare matte-finish film

### Use

Prepping, tourniquet and instrument cover, squaring, isolating sterile from non-sterile areas, fluid barrier, non-linting surface

| Item number | Dimensions                               | Packaging          |
|-------------|--|--------------------|
| IDE1000     | 45 x 30 cm<br>adhesive only: 45 x 5.1 cm | 10 ea/bx; 40 ea/cs |
| IDE1010     | 45 x 61 cm<br>adhesive only: 45 x 5.1 cm | 10 ea/bx; 40 ea/cs |



Easy-peel liner: less force is needed to separate it from the adhesive section. The reinforced adhesive section reduces the chance of wrinkling.

## Fluid pouches


### Description

- Strong adhesive with no residue
- Non-linting, tear-resistant, waxy paper release liner
- Mouldable strip to enhance pouch performance

### Use

ENT, general

| Item number | Dimensions                                     | Packaging          |
|-------------|--|--------------------|
| IDE2114     | 30.5 x 27.9 cm<br>Adhesive only: 30.5 x 2.5 cm | 10 ea/bx; 40 ea/cs |



## Isolation bags

### Description

- Soft, strong, clear film with seamless bottom
- Cotton drawstrings for easy closure
- Cuffed opening to improve sterile delivery

### Use

Isolate foot, hand, non-linting surface





| Item number | Dimensions | Packaging           |
|-------------|------------|---------------------|
| ICE1003     | 50 x 50 cm | 10 ea/bx; 40 ea/cs  |
| 586560-1    | 60 x 65 cm | 30 ea/bx; 180 ea/cs |

## Fluid pouches with drainage ports

### Description

- Wide adhesive area for better adhesion
- Strong adhesive with no residue
- Non-linting, tear-resistant, waxy paper release liner
- Mouldable strip to enhance pouch performance
- Micro screen for improved filtering system
- Suction port with cap for improved fluid control
- Heavy, low-glare matte-finish film



Extra-long, mouldable strip allows for customisation.

### Use

CV, orthopaedic, neurology, urology, general

| Item number | Dimensions                              | Packaging |
|-------------|---|-----------|
| IDE1016     | 50 x 50 cm<br>Adhesive only: 50 x 15 cm | 40 ea/cs  |
| 4638CE      | 61 x 51 cm<br>Adhesive only: 51 x 13 cm | 40 ea/cs  |
| DYJPEFP020  | 51 x 56 cm<br>Adhesive only: 46 x 5 cm  | 50 ea/cs  |


## Instrument pouches

### Description

- 2 compartments for cauteries or catheters
- New high-tack, aggressive adhesive for better adhesion to linen or non-woven materials
- Heavy, low-glare matte-finish film

### Use

CV, orthopaedic, neurology



| Item number | Dimensions                                   | Packaging          |
|-------------|--|--------------------|
| IDE1018     | 17 x 29 cm<br>Adhesive only (2): 17 x 3.8 cm | 10 ea/bx; 40 ea/cs |


## Long instrument pouches

### Description

- 3 long compartments for cauteries or catheters
- New high-tack, aggressive adhesive for better adhesion to linen or non-woven materials
- Heavy, low-glare matte-finish film

### Use

CV, orthopaedic, neurology



| Item number | Dimensions                                 | Packaging          |
|-------------|--|--------------------|
| IDE1018L    | 25 x 47 cm<br>Adhesive only (2): 25 x 5 cm | 10 ea/bx; 40 ea/cs |

## Instrument pouches

### Description

- 1 large compartment for instruments
- New high-tack adhesive tape
- Clear, strong polyethylene film

### Use




CV, orthopaedic, neurology, general

| Item number  | Dimensions                                | Packaging |
|--------------|---|-----------|
| DYJPEISTPOM1 | 33 x 41 cm<br>Adhesive only: 41 x 5 cm    | 100 ea/cs |
| DYJPEISTPOM2 | 33 x 41 cm<br>Adhesive only: 41 x 5 cm    | 100 ea/cs |
| DYJPEISTPOM3 | 43 x 48 cm<br>Adhesive (2): 43 x 5 cm     | 60 ea/cs  |
| 04A040CE     | 33 x 51 cm<br>Adhesive only: 51 x 5 cm    | 180 ea/cs |
| 96860CE      | 30.5 x 29 cm<br>Adhesive (2): 30.5 x 5 cm | 150 ea/cs |
| 96890CE      | 43 x 49 cm<br>Adhesive (2): 49 x 5 cm     | 250 ea/cs |
| 9686CE*      | 27 x 29 cm<br>Adhesive (2): 29 x 5 cm     | 150 ea/cs |
| 9687CE*      | 36 x 86 cm<br>Adhesive (2): 36 x 5 cm     | 120 ea/cs |
| 9689CE*      | 43.5 x 49 cm<br>Adhesive (2): 49 x 5 cm   | 150 ea/cs |

\*Made of spunlace material

## Set-up covers

| Item number | Dimensions           | Packaging |
|-------------|----------------------|-----------|
| ICE2158     | 46 x 61 cm, adhesive | 40 ea/cs  |
| ICE2008     | 137 x 229 cm         | 20 ea/cs  |
| ICE35975    | 150 x 190 cm         | 40 ea/cs  |



## Orthopaedic drapes

### Polyethylene U-drape



#### Description

- 'U' fenestration with wide adhesive area
- Advanced adhesive, providing excellent adhesion without residue
- Strong, drapeable, low-glare matte-finish film
- Split, wax paper liner for ease of application

#### Use

Prepping procedures: arthroscopies, hip, knee, shoulder, back, cranial, isolation and applications

| Item number | Dimensions   | Packaging          |
|-------------|--|--------------------|
| IDE1015     | 121 x 132 cm<br>Adhesive only: 25 x 48 cm                                      | 10 ea/bx; 40 ea/cs |
| IDE1019     | 121 x 180 cm<br>Adhesive only: 25 x 48 cm<br>'U' fenestration only: 15 x 48 cm | 10 ea/bx; 40 ea/cs |



These drapes include courtesy tabs for easy removal.

### Shoulder pouches


#### Description

- Oval aperture with slit and adhesive areas
- Strong adhesive with no residue
- Heavy, low-glare matte-finish film for additional strength
- 2 mouldable strips and 1 foam strip to enhance pouch performance
- 2 suction ports with caps for improved fluid control

#### Use

Shoulder arthroscopy, implant applications

| Item number | Dimensions  | Packaging |
|-------------|---|-----------|
| IDE1067     | 91 x 73 cm<br>Adhesive only: 25 x 48 cm<br>Aperture: 19 x 25 cm | 24 ea/cs  |



The capped, closable port eliminates the step of cutting or piercing for drainage and makes it easy to stop the flow at the end of the procedure.


## U fluid collection pouches

### Description

- 'U' aperture with adhesive area
- Strong adhesive with no residue
- 3 mouldable strips to enhance pouch performance
- Strong, drapeable, clear polyethylene film

### Use

Shoulder arthroscopy, hip surgery, orthopaedic



| Item number | Dimensions   | Packaging |
|-------------|--|-----------|
| DYJPEFP090  | 76 x 112 cm<br>Adhesive only: 25 x 51 cm<br>Aperture: 15 x 46 cm | 25 ea/cs  |


## Patient isolation drapes

### Description

- Less aggressive adhesive strip for easy removal
- Aggressive incise adhesive with no residue
- Clear, strong polyethylene film

### Use

Hip pinning



| Item number | Dimensions   | Packaging         |
|-------------|--|-------------------|
| IDE1014     | 317 x 213 cm<br>Adhesive only: 49 x 25 cm<br>Adhesive strip only: 317 x 1.3 cm | 5 ea/bx; 15 ea/cs |

## Patient isolation drapes with pouch


### Description

- Strong incise adhesive with no residue
- Clear, strong polyethylene film
- 2 cord holders and instrument pouches for better organisation
- Large pouch with micro screen, suction port and cap for improved fluid control

### Use

Hip pinning, nephroscopies, biliaries

| Item number | Dimensions   | Packaging |
|-------------|--|-----------|
| IDE1017     | 317 x 213 cm<br>Adhesive only: 49 x 25 cm<br>Adhesive strip only: 317 x 1.3 cm<br>Pouch only: 51 x 51 cm       | 15 ea/cs  |
| IDE1017W    | 317 x 254 cm<br>Adhesive only: 85.1 x 27.9 cm<br>Adhesive strip only: 317 x 1.3 cm<br>Pouch only: 69.9 x 94 cm | 15 ea/cs  |
| 5801CE      | 330 x 270 cm<br>Adhesive only: 78 x 26 cm<br>Adhesive strip only: 330 x 1.9 cm<br>Pouch only: 25.5 x 28 cm     | 16 ea/cs  |



These drapes have a superb filter with deep placement for quick drainage and numerous small perforations.

## Ophthalmic drapes

### Aperture drapes



### Description

- Circular aperture with adhesive area
- Non-linting, tear-resistant, wax paper release liner
- Thin, soft, drapeable, low-glare matte-finish film

### Use

Cataract, corneal, refractive, lid/plastic, minor

| Item number | Dimensions  | Packaging          |
|-------------|---|--------------------|
| IDE1020     | 40 x 40 cm<br>Adhesive only: 12 x 12 cm<br>Aperture only: 6 cm circle | 10 ea/bx; 40 ea/cs |



Easy-peel liner: less force is needed to separate it from the adhesive section.

## Ophthalmic drapes

### Description

- Incise adhesive area optimised for skin
- Non-linting, tear-resistant, wax paper release liner
- 2 adhesive tabs for additional support
- Mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film

### Use

Cataract, vitreoretinal, corneal



| Item number | Dimensions  | Packaging          |
|-------------|---|--------------------|
| IDE1023     | 48 x 41 cm<br>Adhesive only: 12 x 10 cm<br>Adhesive tabs only: 5 x 5 cm | 10 ea/bx; 40 ea/cs |


## Ophthalmic aperture pouches

### Description

- Oval aperture with adhesive optimised for skin
- Non-linting, tear-resistant, wax paper release liner
- 2 adhesive tabs for additional support
- Mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film

### Use

Cataract, vitreoretinal, corneal



| Item number | Dimensions   | Packaging          |
|-------------|--|--------------------|
| IDE1024     | 48 x 41 cm<br>Adhesive only: 6 x 3.5 cm<br>Adhesive tabs only: 5 x 5 cm<br>Aperture only: 6 x 3.5 cm | 10 ea/bx; 40 ea/cs |


## Aperture drapes

### Description

- Oval aperture with adhesive area
- Non-linting, tear-resistant, wax paper release liner
- Soft, drapeable, low-glare matte-finish film

### Use

Cataract, vitreoretinal, corneal, lid/plastic, ENT



| Item number | Dimensions  | Packaging          |
|-------------|---|--------------------|
| IDE1030     | 83 x 129 cm<br>Adhesive only: 12 x 11 cm<br>Aperture only: 7 x 5 cm | 10 ea/bx; 40 ea/cs |


## Ophthalmic incise drapes with fluid pouch

### Description


- Incise adhesive area
- Non-linting, tear-resistant, wax paper release liner
- Fluid pouch with mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film

### Use

Cataract, vitreoretinal



| Item number | Dimensions  | Packaging          |
|-------------|---|--------------------|
| IDE1061     | 121 x 129 cm<br>Adhesive only: 12 x 11 cm<br>Pouch only: 30 x 27 cm | 10 ea/bx; 40 ea/cs |



# Other drapes

## Aperture pouches


### Description

- Oval aperture with large adhesive area
- Strong adhesive with no residue
- Mouldable strip for additional fluid control
- Heavy, low-glare matte-finish film

### Use

Vaginal, ophthalmic

| Item number | Dimensions  | Packaging          |
|-------------|---|--------------------|
| IDE1072     | 45 x 22 cm<br>Adhesive only: 30 x 22 cm<br>Aperture only: 11 x 9 cm | 10 ea/bx; 40 ea/cs |



## Clear adhesive drape

| Item number | Dimensions   | Packaging |
|-------------|--------------|-----------|
| 5810CE      | 150 x 300 cm | 17 ea/cs  |

## Urological drapes


### Description

- Strong adhesive with no residue
- Latex-free finger cot
- Strong, low-glare matte-finish film
- Non-linting, tear-resistant, wax paper release liner

### Use


TUR, urology, gynaecology

| Item number | Dimensions  | Packaging          |
|-------------|---|--------------------|
| IDE1071     | 60 x 91 cm<br>Adhesive only: 60 x 5 cm<br>Finger cot adhesive only: O.D. 10 x I.D. 3 cm | 10 ea/bx; 40 ea/cs |



## Clear epidural drapes

| Item number | Dimensions  | Packaging |
|-------------|---|-----------|
| 8461CE      | 60 x 100 cm, diameter 9 cm<br>Adhesive fenestration                                   | 150 ea/cs |
| 8465CE      | 80 x 140 cm, 6 x 10 cm<br>Adhesive oval fenestration<br>Pre-cut line for easy removal | 150 ea/cs |



## Cesarean pouches


### Description

- High-moisture, permeable incise film for extended adhesion
- Strong adhesive with no residue
- Large-cell foam to provide fluid control
- 360° fluid pouch for complete fluid containment
- 2 cord holders for line organisation
- Heavy, low-glare matte-finish film for additional pouch strength

### Use

C-section, trauma, general

| Item number | Dimensions                              | Packaging |
|-------------|---|-----------|
| IDE2057     | 91 x 76 cm<br>Adhesive only: 38 x 36 cm | 10 ea/cs  |



## Wound ring protectors


### Description

- Strong, drapeable, low-glare matte-finish film
- 5 ring diameters to meet different surgical needs
- Application with adhesive tape

### Use

General, laparotomy, cardiothoracic, colectomy

| Item number | Dimensions                       | Packaging |
|-------------|----------------------------------|-----------|
| IDE3707     | 90 x 90 cm, ring diameter: 7 cm  | 40 ea/cs  |
| IDE3712     | 90 x 90 cm, ring diameter: 12 cm | 40 ea/cs  |
| IDE3718     | 90 x 90 cm, ring diameter: 18 cm | 40 ea/cs  |
| IDE3723     | 90 x 90 cm, ring diameter: 23 cm | 40 ea/cs  |
| IDE3728     | 90 x 90 cm, ring diameter: 28 cm | 40 ea/cs  |




The wound ring protectors have a flexible but rigid ring at the fenestration site to help keep the drape in the correct position. The frosted clear plastic drape around the fenestration helps maintain the sterile barrier.

# Dermadine® and Dermadrape® incise drapes

## Description

- Thin, conformable, low-glare matte-finish film
- Strong adhesive with no residue
- Breathable, water-repellent and waterproof
- Comfortable and enveloping
- Long-lasting incise drapes
- Dermadine incise drapes: impregnated with an antibacterial solution
- Dermadrape drapes: clear incise drapes



| Dermadine item number |   | Dimensions                         | Packaging             | Dermadrape item number |   |
|-----------------------|---|------------------------------------|-----------------------|------------------------|---|
| 42020                 |   | 20 x 20 cm<br>Adhesive: 15 x 20 cm | 20 ea/bx<br>120 ea/cs | 32020                  |   |
| n/a                   |   | 35 x 25 cm<br>Adhesive: 30 x 25 cm | 20 ea/bx<br>120 ea/cs | 33525                  |  |
| 44040                 |  | 40 x 40 cm<br>Adhesive: 34 x 40 cm | 15 ea/bx<br>90 ea/cs  | 34040                  |  |
| n/a                   |   | 50 x 48 cm<br>Adhesive: 44 x 48 cm | 15 ea/bx<br>90 ea/cs  | 35050                  |  |
| 47045                 |  | 70 x 45 cm<br>Adhesive: 60 x 45 cm | 15 ea/bx<br>90 ea/cs  | 37045                  |  |
| 48060                 |  | 80 x 60 cm<br>Adhesive: 60 x 60 cm | 10 ea/bx<br>60 ea/cs  | 38060                  |  |
| 48080                 |  | 80 x 80 cm<br>Adhesive: 70 x 80 cm | 10 ea/bx<br>60 ea/cs  | 38080                  |  |


# Equipment drapes

## Circular banded bags

### Description

- Clear, drapeable polyethylene film
- Latex-free elastic bands
- Versatile, often used to cover intensifiers and other equipment

| Item number | Dimensions  | Packaging |
|-------------|-------------|-----------|
| ICE60020S   | 25 x 55 cm  | 25 ea/cs  |
| ICE60028S   | 35 x 70 cm  | 25 ea/cs  |
| ICE60030S   | 40 x 80 cm  | 25 ea/cs  |
| ICE60040S   | 50 x 100 cm | 25 ea/cs  |
| ICE60050S   | 60 x 120 cm | 25 ea/cs  |
| ICEBBG070   | 65 x 130 cm | 50 ea/cs  |
| ICE60060S   | 70 x 165 cm | 10 ea/cs  |




## Rectangular banded bags

### Description

- Clear, drapeable polyethylene film
- Latex-free elastic bands
- Versatile, often used to cover intensifiers and other equipment

| Item number | Dimensions   | Packaging |
|-------------|--------------|-----------|
| ICE63036R   | 76 x 91 cm   | 25 ea/cs  |
| ICE63628R   | 91 x 71 cm   | 25 ea/cs  |
| ICE63636R   | 91 x 91 cm   | 25 ea/cs  |
| ICEBBG095   | 95 x 85 cm   | 25 ea/cs  |
| ICE63654S   | 95 x 130 cm  | 25 ea/cs  |
| ICE64848S   | 122 x 122 cm | 25 ea/cs  |
| ICE65656S   | 142 x 142 cm | 25 ea/cs  |





## C-arm set Philips

### Description

- Clear banded bags for X-ray tank/intensifier
- Clear polyethylene drape to cover the C-arm


DYNJE4526 for BV29/300

DYNJE4527 for BV300

DYNJE4530 for BV25/26

DYNJE4569 for Pulsera/Endura

| Item number | Dimensions   | Packaging             |
|-------------|--|-----------------------|
| DYNJE4526   | Banded bags: 60 x 115 cm   70 x 130 cm<br>C-arm cover: 30 x 198 cm       | 10 ea/bx;<br>20 ea/cs |
| DYNJE4527   | Banded bags: 60 x 115 cm   rect.: 95 x 85 cm<br>C-arm cover: 30 x 198 cm | 10 ea/bx;<br>20 ea/cs |
| DYNJE4530   | Banded bags: 60 x 115 cm<br>C-arm cover: 30 x 198 cm                     | 10 ea/bx;<br>20 ea/cs |
| DYNJE4569   | Banded bags: 60 x 115 cm   70 x 140 cm<br>C-arm cover: 30 x 198 cm       | 10 ea/bx;<br>20 ea/cs |




## C-arm set Ziehm

### Description

- Clear banded bags for X-ray tank/intensifier
- Clear polyethylene drape to cover the C-arm

For starter set SPTE-SS634

| Item number | Dimensions   | Packaging             |
|-------------|--|-----------------------|
| DYNJE4531   | Banded bags: 60 x 115 cm   70 x 130 cm<br>C-arm cover: 35 x 185 cm | 10 ea/bx;<br>20 ea/cs |



## C-arm set Siemens

### Description



- Clear banded bags for X-ray tank/intensifier
- Clear polyethylene drape to cover the C-arm

DYNJE4532 for SIREMOBIL N4

- Clear polyethylene drape with clips for starter set SPTE-SS628

DYNJE4556 for ISO-C/Arcadis/Orbic, 2 parts

| Item number | Dimensions   | Packaging             |
|-------------|--|-----------------------|
| DYNJE4532   | Banded bags: 60 x 115 cm   60 x 115 cm<br>C-arm cover: 35 x 183 cm | 10 ea/bx;<br>20 ea/cs |
| DYNJE4547   | Banded bags: 70 x 130 cm   60 x 115 cm<br>C-arm cover: 45 x 183 cm | 20 ea/cs              |
| DYNJE4556   | Top: 51 – 82 x 130 cm<br>Bottom: 53 – 46 x 134 cm                  | 10 ea/bx;<br>20 ea/cs |




## Fluoroscope drapes

### Description

- 2 adhesive strips for attachment, with 2 adhesive tabs for proper fit
- Clear, drapeable polyethylene film

| Item number | Dimensions   | Packaging             |
|-------------|--|-----------------------|
| ICE1012     | 91 x 111 cm<br>Adhesive only (2): 91 x 1.2 cm<br>Adhesive tabs only: 5 x 10 cm | 10 ea/bx;<br>40 ea/cs |



## Full-size C-arm drapes

### Description

- 3 adhesive polyethylene straps for securement
- Clear, drapeable polyethylene film

ICE4400R: has 2 latex-free elastic bands for additional securement

| Item number | Dimensions   | Packaging |
|-------------|--------------|-----------|
| ICE4400     | 104 x 188 cm | 20 ea/cs  |
| ICE4400R    | 104 x 188 cm | 20 ea/cs  |
| ICE4405     | 90 x 225 cm  | 75 ea/cs  |




## Contoured C-arm drapes

### Description

- 4 adhesive polyethylene straps for securement
- Clear, drapeable polyethylene film
- Contoured to fit full C-arm

| Item number | Dimensions   | Packaging |
|-------------|--------------|-----------|
| ICE4430     | 152 x 267 cm | 20 ea/cs  |



## Mini C-arm drapes


### Description

- Held in place with an elastic band on the back
- 3 adhesive polyethylene straps for additional securement
- Clear, drapeable polyethylene film

ICE4410: fits Fluoroscan® mini C-arm models

ICE4420: fits Fluoroscan® and OEC6600/6800 mini C-arm models

| Item number | Dimensions   | Packaging |
|-------------|--------------|-----------|
| ICE4410     | 137 x 160 cm | 20 ea/cs  |
| ICE4420     | 137 x 216 cm | 20 ea/cs  |




## Full-size C-arm kits

### Description

- Clear, drapeable polyethylene film
- Latex-free bands

| Item number | Dimensions  | Packaging |
|-------------|---|-----------|
| ICE4500     | Intensifier drape (1): 91 x 91 cm<br>C-Arm drape (1): 41 x 218 cm<br>Footswitch drape (1): 46 x 41 cm<br>Banded bag (1): 71 x 31 cm | 20 ea/cs  |
| ICE8130     | C-Arm drape (1): 10/14 x 225 cm<br>Banded bag (2): 75 x 90 cm   | 20 ea/cs  |
| ICE8132     | C-Arm drape (1): 10/14 x 225 cm<br>Banded bag (1): 75 x 90 cm   | 20 ea/cs  |
| ICE8230     | C-Arm drape (1): 8/10 x 205 cm<br>Banded bag (2): 75 x 75 cm  | 20 ea/cs  |
| ICE8232     | C-Arm drape (1): 8/10 x 205 cm<br>Banded bag (1): 75 x 75 cm  | 20 ea/cs  |
| ICE8330     | C-Arm drape (1): 8/10 x 180 cm<br>Banded bag (2): 75 x 75 cm  | 20 ea/cs  |
| ICE8332     | C-Arm drape (1): 8/10 x 180 cm<br>Banded bag (1): 75 x 75 cm  | 20 ea/cs  |



## Remote control covers

### Description

- Clear polyethylene film bag
- Pinch-top closure

| Item number | Dimensions | Packaging              |
|-------------|------------|------------------------|
| DYNJE4635   | 10 x 28 cm | 50 ea/bx;<br>200 ea/cs |
| DYNJE4641   | 18 x 35 cm | 50 ea/bx;<br>200 ea/cs |




## Video camera covers

### Description

- Adhesive polyethylene strap for holding it in place
- Adhesive strip at base for additional securement
- Clear, drapeable polyethylene film
- Cable cover

| Item number | Dimensions                      | Packaging |
|-------------|---------------------------------|-----------|
| ICE4250     | 13 x 244 cm reverse             | 25 ea/cs  |
| ICE4300     | 23 x 244 cm reverse             | 25 ea/cs  |
| ICE4240     | 13 x 240 cm telescopic with tab | 25 ea/cs  |
| ICE4260     | 15 x 240 cm with ring           | 25 ea/cs  |
| ICE4270     | 15 x 298 cm concertina          | 25 ea/cs  |
| ICE4200     | 18 x 244 cm elastic end         | 25 ea/cs  |
| ICE4100     | 8 x 200 cm telescopic with tab  | 25 ea/cs  |
| ICE4245     | 18 x 250 cm telescopic with tab | 25 ea/cs  |



## Microscope drapes

### Description

- All fitted with clear, drapeable polyethylene film

#### ICE5010, ICE5040

- Fits 48 mm lens (Zeiss® OPML)
- 3 adhesive PE straps for securement
- Fits 1 pair of oculars

#### ICE5200, ICE5210, ICE5030

- Fits 65 mm lens (Zeiss MD)
- 3 adhesive PE straps for securement
- Fits 3 pairs of oculars


#### ICE5220

- Fits 65 mm lens (Zeiss MD)
- 5 adhesive PE straps for securement
- Fits 3 pairs of oculars

#### ICE5300, ICE5330

- Fits 68 mm lens (Leica®)
- 3 adhesive PE straps for securement
- Fits 3 pairs of oculars

| Item number                         | Dimensions   | Packaging |
|-------------------------------------|--------------|-----------|
| <b>Fits 48 mm lens (Zeiss OPML)</b> |              |           |
| ICE5010                             | 51 x 163 cm  | 20 ea/cs  |
| ICE5040                             | 51 x 267 cm  | 40 ea/cs  |
| ICE5020                             | 104 x 305 cm | 10 ea/cs  |
| <b>Fits 65 mm lens (Zeiss MD)</b>   |              |           |
| ICE5200                             | 117 x 203 cm | 10 ea/cs  |
| ICE5030                             | 117 x 305 cm | 10 ea/cs  |
| ICE5210                             | 117 x 267 cm | 10 ea/cs  |
| ICE5220                             | 137 x 381 cm | 10 ea/cs  |
| <b>Fits 68 mm lens (Leica)</b>      |              |           |
| ICE5300                             | 104 x 163 cm | 10 ea/cs  |
| ICE5320                             | 137 x 267 cm | 10 ea/cs  |
| ICE5330                             | 137 x 381 cm | 10 ea/cs  |




## X-ray cassette covers

### Description

- Clear, drapeable polyethylene film

ICE4010: release liner for facilitating sterile technique

| Item number | Dimensions                      | Packaging |
|-------------|---------------------------------|-----------|
| ICE4000     | Large: 61 x 102 cm              | 20 ea/cs  |
| ICE4010     | Standard with liner: 61 x 64 cm | 48 ea/cs  |
| ICE4050     | Standard: 61 x 64 cm            | 50 ea/cs  |



## Probe cover kits

### Description

- 1 probe cover: soft and drapeable polyethylene or polyurethane film
- 1 ultrasound gel packet: high-quality, sterile ultrasound gel for transducer
- 2 elastic bands for securement

| Item number              | Dimensions  | Packaging |
|--------------------------|-------------|-----------|
| <b>Kits with gel</b>     |             |           |
| ICE5900                  | 10 x 122 cm | 20 ea/cs  |
| ICE5910                  | 15 x 61 cm  | 20 ea/cs  |
| ICE5920                  | 15 x 122 cm | 20 ea/cs  |
| ICE5930                  | 15 x 244 cm | 20 ea/cs  |
| <b>Cover without gel</b> |             |           |
| ICE18120                 | 18 x 120 cm | 25 ea/cs  |
| ICE1830                  | 18 x 30 cm  | 25 ea/cs  |
| ICE1880                  | 18 x 80 cm  | 25 ea/cs  |




High-quality polyurethane is soft and drapeable to conform to the transducer and has a smooth finish for easy application and glide. The slightly matte surface helps prevent glare.

## Thermo cover

### Description

- Clear polyethylene film
- Aluminium coated foil
- Temperature resistance up to 60° C

| Item number | Dimensions | Packaging |
|-------------|------------|-----------|
| DYNJE7000   | 68 x 86 cm | 40 ea/cs  |



# Stockinettes

## Impervious stockinettes


### Description

- Offers exceptional, high-level barrier properties and optimal fluid control, thanks to the use of an impervious thermoplastic elastomer material
- Gives the patient comfort by way of a soft, knitted jersey fabric on the inside.
- Available in 6 sizes

### Use

Orthopaedic, arthroplasty, as well as hand, knee, ankle, foot, limb and vascular surgeries

| Item number | Width 'A' | Width 'B' | Packaging |
|-------------|-----------|-----------|-----------|
| 22520CE     | 15 cm     | 76 cm     | 20 ea/cs  |
| 22525CE     | 20 cm     | 97 cm     | 20 ea/cs  |
| 22530CE     | 23 cm     | 91 cm     | 20 ea/cs  |
| 22534CE     | 23 cm     | 122 cm    | 16 ea/cs  |
| 22535CE     | 30 cm     | 122 cm    | 12 ea/cs  |
| 22540CE     | 36 cm     | 122 cm    | 12 ea/cs  |



## Impervious stretchable stockinettes


### Description

- Allows for an ideal fit over the patient's anatomy, due to the use of a stretchable, bilaminate material
- Facilitates set-up with the help of 2 pulling ties
- Both lower and upper extremities can be covered, as it is available in 2 sizes

### Use

Orthopaedic, arthroplasty, as well as hand, knee, ankle, foot, limb and vascular surgeries

| Item number | Width 'A' | Width 'B' | Packaging |
|-------------|-----------|-----------|-----------|
| 22620CE     | 17 cm     | 76 cm     | 18 ea/cs  |
| 22634CE     | 24 cm     | 118 cm    | 10 ea/cs  |



# Invisishield cross reference

| Page | Description                              | Product number Medline | Product number 3M® | Product number Microtek® |
|------|--|------------------------|--------------------|--------------------------|
| 3    | Towel drape                              | IDE1000                | 1000               | ISO1000                  |
| 4    | Isolation bag                            | ICE1003                | 1003               | ISO1003                  |
| 3    | Towel drape                              | IDE1010                | 1010               | ISO1010                  |
| 8    | Patient isolation drape                  | DYNJSD1014             | 1014               | n/a                      |
| 7    | Polyethylene U drape                     | IDE1015                | 1015               | ISO1015                  |
| 4    | Fluid pouch                              | IDE1016                | 1016               | ISO1016                  |
| 9    | Patient isolation drape                  | DYNJSD1017             | 1017               | 5960                     |
| 9    | Wide patient isolation drape             | DYNJSD1017W            | n/a                | 6960                     |
| 5    | Instrument pouch                         | IDE1018                | 1018               | ISO1018                  |
| 5    | Long instrument pouch                    | IDE1018L               | 1018L              | n/a                      |
| 7    | Polyethylene U drape                     | IDE1019                | 1019               | ISO1019                  |
| 10   | Aperture drape                           | IDE1020                | 1020               | ISO1020                  |
| 10   | Ophthalmic incise drape                  | IDE1023                | 1023               | n/a                      |
| 11   | Ophthalmic aperture pouch                | IDE1024                | 1024               | ISO1024                  |
| 11   | Aperture drape                           | IDE1030                | 1030               | ISO1030                  |
| 14   | Incise drape                             | 32020                  | 1035               | ISO1035                  |
| 14   | Incise drape                             | 34040                  | 1040               | ISO1040                  |
| 14   | Incise drape                             | 37045                  | 1050               | ISO1050                  |
| 14   | Incise drape                             | 38080                  | 1051               | ISO1051                  |
| 12   | Ophthalmic incise drape with fluid pouch | IDE1061                | 1061               | ISO1061                  |
| 7    | Shoulder pouch                           | DYNJSD1067             | 1067               | ISO1067                  |
| 12   | Urological drape                         | IDE1071                | 1071               | ISO1071                  |
| 13   | Aperture pouch                           | IDE1072                | 1072               | ISO1072                  |
| 8    | U fluid collection pouch                 | DYJPEFP090             | n/a                | n/a                      |
| 6    | Instrument pouch                         | DYJPEISTPOM1           | n/a                | n/a                      |
| 6    | Instrument pouch                         | DYJPEISTPOM2           | n/a                | n/a                      |
| 6    | Instrument pouch                         | DYJPEISTPOM3           | n/a                | n/a                      |

 All products are latex free and sterile, unless indicated in the description.

 Adhesive
  Matte film
  Polyethylene film

# Dermadine cross reference

| Type of film       | Description  | Product number Medline | Product number 3M® |
|--------------------|--|------------------------|--------------------|
| PU                 | Incise, sterile drape 15 x 20 cm; incise area 10 x 20 cm, Iodophor impregnated | 42020                  | 6635               |
| PU                 | Incise, sterile drape 44 x 35 cm; incise area 34 x 35 cm, Iodophor impregnated | 44040                  | 6640               |
| PU                 | Incise, sterile drape 66 x 60 cm; incise area 56 x 60 cm, Iodophor impregnated | 48060                  | 6648               |
| PU                 | Incise, sterile drape 66 x 45 cm; incise area 56 x 45 cm, Iodophor impregnated | 47045                  | 6650               |
| PU                 | Incise, sterile drape 66 x 85 cm; incise area 56 x 85 cm, Iodophor impregnated | 48080                  | 6651               |
| PU with applicator | Incise, sterile drape 30 x 20 cm; incise area 26 x 20 cm, Iodophor impregnated | 42020                  | 6661EZ             |
| PU with applicator | Incise, sterile drape 60 x 35 cm; incise area 35 x 35 cm, Iodophor impregnated | 44040                  | 6640EZ             |
| PU with applicator | Incise, sterile drape 90 x 60 cm; incise area 60 x 60 cm, Iodophor impregnated | 48060                  | 6648EZ             |
| PU with applicator | Incise, sterile drape 90 x 45 cm; incise area 60 x 45 cm, Iodophor impregnated | 47045                  | 6650EZ             |
| PU with applicator | Incise, sterile drape 90 x 85 cm; incise area 60 x 85 cm, Iodophor impregnated | 48080                  | 6651EZ             |

# Dermadrape cross reference

| Type of film | Description  | Product number Medline | Product number 3M® |
|--------------|--|------------------------|--------------------|
| PE           | Incise, sterile drape 15 x 20 cm; incise area 10 x 20 cm | 32020                  | 1035               |
| PE           | Incise, sterile drape 29 x 25 cm; incise area 19 x 25 cm | 32020                  | 1037               |
| PE           | Incise, sterile drape 60 x 35 cm; incise area 35 x 35 cm | 34040                  | 1040               |
| PE           | Incise, sterile drape 82 x 45 cm; incise area 60 x 45 cm | 37045                  | 1050               |
| PE           | Incise, sterile drape 82 x 85 cm; incise area 60 x 85 cm | 38080                  | 1051               |
| PU           | Incise, sterile drape 20 x 20 cm; incise area 10 x 20 cm | 32020                  | 2035               |
| PU           | Incise, sterile drape 38 x 25 cm; incise area 28 x 25 cm | 33525                  | 2037               |
| PU           | Incise, sterile drape 38 x 41 cm; incise area 28 x 41 cm | 34040                  | 2040               |
| PU           | Incise, sterile drape 60 x 45 cm; incise area 50 x 45 cm | 35050                  | 2045               |
| PU           | Incise, sterile drape 90 x 60 cm; incise area 60 x 60 cm | 38060                  | 2050               |
| PU           | Incise, sterile drape 90 x 85 cm; incise area 60 x 85 cm | 38080                  | 2051               |



**Medline Industries Ltd**  
 3rd Floor  
 Quayside Wilderspool Business Park  
 Greenalls Avenue  
 Warrington WA4 6HL  
 United Kingdom  
 Tel.: +44 844 334 5237  
 Fax: +44 844 334 5238  
[www.medline.eu/uk](http://www.medline.eu/uk) | [uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

FOLLOW US  

Dermadine® and Dermadrape® are class IIa sterile medical devices.

Wound ring protectors are class IIa sterile medical devices.

All other products listed are class I sterile medical devices.

These products are intended to be used by healthcare professionals.

Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2020 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093, USA. ML164-EN 06/2020.

