

# Invisishield™ Drapes

Unique products to help improve safety and patient care



**Protect and  
see clearly**



**Clear,  
impervious  
and drapeable**

Medline's Invisishield drapes are clear so that clinicians can easily see theatre equipment and incision areas. The impervious, tear-resistant material provides exceptional protection and offers high-level drapeability.

Invisishield drapes have handy features that help make clinicians' job easier. Our fluid pouches, for example, feature handy tabs on an adhesive liner that make removal easy. Capped, closable ports also eliminate the need to cut or pierce the pouch for drainage, thereby improving fluid control and reducing set-up times.



### **Transparent**

Clinicians can easily view the area around incision sites.



### **Impervious**

Made of tear-resistant polyethylene (PE) and polyurethane films, Invisishield drapes provide a high level of protection from fluids.



### **Drapeable**

Invisishield drapes have just the right amount of flexibility and softness for a higher level of drapeability.







Clear polyethylene



Clear polyethylene




Clear, frosted  
polyethylene




Clear, frosted  
polyethylene



Clear polyurethane



Clear polyurethane




Clear polyurethane  
with iodine




Clear polyurethane  
with iodine



SMS




SMS



Spunlace




Spunlace




Stretchable  
bilaminate



Stretchable  
bilaminate



Thermoplastic  
elastomer



Thermoplastic  
elastomer


## Invisishield Raw Materials



# Table of Contents

<b>General Accessories</b>	<b>5</b>	<b>Video Camera Covers</b>	<b>24</b>
Towel drapes	5	Video camera covers	24
Isolation bags	5		
Fluid pouches with drainage ports	5	<b>Microscope Covers</b>	<b>25</b>
Instrument pouches	6	Microscope covers	25
Long instrument pouches	6		
Set-up covers	7	<b>Probe Covers</b>	<b>27</b>
Sterile sponge counter bags	7	Probe cover kits	27
Non-sterile sponge counter bags	7	Probe cover kits with adhesive	27
		Endocavity probe covers	28
<b>Orthopaedic Drapes</b>	<b>8</b>		
Polyethylene U-drapes	8	<b>Other Covers</b>	<b>29</b>
Shoulder pouch	8	Remote control covers	29
U-fluid collection pouch	9	Thermo cover	29
Patient isolation drapes	9	X-ray cassette covers	29
<b>Ophthalmic Drapes</b>	<b>10</b>	<b>Stockinettes</b>	<b>30</b>
Circular/oval aperture drapes	10	Stockinettes	30
Ophthalmic aperture pouches	10		
Ophthalmic incise drapes with fluid pouch	11	<b>References</b>	<b>31</b>
		Invisishield Cross-Reference Guide	31
<b>Other Drapes</b>	<b>12</b>	Dermadine Plus Cross-Reference Guide	31
Aperture drape	12	Dermadrape Cross-Reference Guide	32
Urology drape	12	Stretchable Stockinettes Cross-Reference Guide	32
Clear epidural drapes	12		
Cesarean pouch	13		
Wound ring protectors	13		
<b>Incise Drapes</b>	<b>14</b>		
Dermadine Plus and Dermadrape	15		
<b>Invisishield in the Cath Lab</b>	<b>17</b>		
<b>C-arm</b>			
Customised C-arm covers	18		
Fluoroscope drape	19		
Full-size C-arm drapes	20		
Contoured C-arm drape	20		
Mini C-arm drapes	20		
Full-size C-arm kits	21		
<b>Banded Bags</b>			
Banded bags	22		
Banded bags for X-ray shields	23		
Cover for thyroid collars	23		

## Towel drapes



**Easy-peel liner:** Less force is needed to separate it from the adhesive section. The reinforced adhesive section reduces the chance of wrinkling.

### Benefits and features

- Drapeable with strong residue-free adhesive
- Non-linting, tear-resistant, release liner
- Low-glare matte-finish film

**Use:** Prepping, tourniquet and instrument cover, squaring, isolating sterile from non-sterile areas, fluid barrier, non-linting surface

Item number	Dimensions	Packaging
IDE1000	45 x 30 cm adhesive only: 45 x 5.1 cm	10 ea/bx; 40 ea/cs
IDE1010	45 x 61 cm adhesive only: 45 x 5.1 cm	10 ea/bx; 40 ea/cs

## Isolation bags




### Benefits and features

- Easy closure facilitated by cotton drawstrings
- Helps improve sterile delivery via cuffed opening
- Soft, strong, clear film with seamless bottom

**Use:** Isolate foot, hand, non-linting surface

Item number	Dimensions	Packaging
ICE1003	50 x 50 cm	10 ea/bx; 40 ea/cs
586560-1	60 x 65 cm	30 ea/bx; 180 ea/cs

## Fluid pouches with drainage ports



### Benefits and features

- Better adhesion made possible by wide adhesive area
- Enhanced pouch performance, thanks to mouldable strip
- Improved fluid control facilitated by suction port with cap
- Non-linting, tear-resistant, paper release liner
- Heavy, low-glare matte-finish film

**Use:** CV, orthopaedic, neurology, urology, general










Item number	Dimensions	Packaging
IDE1016	50 x 50 cm Adhesive only: 50 x 15 cm	40 ea/cs
4638CE	61 x 51 cm Adhesive only: 51 x 13 cm	40 ea/cs
DYJPEFP020	51 x 56 cm Adhesive only: 46 x 5 cm	50 ea/cs

## Instrument pouches

### Benefits and features

- Improved visibility and durability via clear, strong polyethylene film
- Facilitates holding instruments via large compartment
- New high-tack adhesive tape

**Use:** CV, orthopaedic, neurology, general

1 compartment				2 compartments			
Item number		Dimensions	Packaging	Item number		Dimensions	Packaging
9686CE		27 x 29 cm Adhesive (2): 29 x 5 cm	150 ea/cs	DYJPEISTPOM1		33 x 41 cm Adhesive only: 41 x 5 cm	100 ea/cs
96860CE		30.5 x 29 cm Adhesive (2): 30.5 x 5 cm	150 ea/cs	DYJPEISTPOM3		43 x 48 cm Adhesive (2): 43 x 5 cm	60 ea/cs
DYJPEISTPOM2		33 x 45 cm Adhesive only: 41 x 5 cm	100 ea/cs	96890CE		43 x 49 cm Adhesive (2): 49 x 5 cm	250 ea/cs
04A040CE		33 x 51 cm Adhesive only: 51 x 5 cm	180 ea/cs	9687CE		36 x 86 cm Adhesive (2): 36 x 5 cm	120 ea/cs
				9689CE		43.5 x 49 cm Adhesive (2): 49 x 5 cm	150 ea/cs



## Long instrument pouches

### Benefits and features

- Better grip to linen or non-woven materials by way of new high-tack, strong adhesive
- Heavy, low-glare matte-finish film


**Use:** CV, orthopaedic, neurology

**Use:** CV, laparoscopy, gynaecology

2 compartment				3 compartments			
Item number		Dimensions	Packaging	Item number		Dimensions	Packaging
IDE1018		17 x 29 cm Adhesive only (2): 17 x 3.8 cm	10 ea/bx; 40 ea/cs	IDE1018L		25 x 47 cm Adhesive only (2): 25 x 5 cm	10 ea/bx; 40 ea/cs




## Set-up covers



Item number	Dimensions	Packaging
ICE2158	41 x 61 cm (adhesive)	40 ea/cs
ICE2008	137 x 229 cm	20 ea/cs
5810CE	150 x 300 cm	17 ea/cs

## Sterile sponge counter bags



### Benefits and features

- Clear, drapeable polyethylene film
- 1 adhesive strip for attachment
- 3 sizes
- 3 compartment options

Sterile sponge counter bags allow for safe, effective and convenient counting, handling and disposal of contaminated swabs. They minimise the risk of pathogen inhalation and keep the surgical area clean and organised.

Item number	Dimensions	Number of compartments	Packaging
IBE5242	52 x 42 cm	40	50 ea/cs
IBE5637	56 x 37 cm	6	50 ea/cs
IBE5652	56 x 52 cm	20	50 ea/cs

## Non-sterile sponge counter bags





### Benefits and features

- Clear, drapeable polyethylene film
- Available with both a transparent and blue backing
- Can be hung from sponge counter rack NON50511, which attaches to most IV poles

Non-sterile sponge counter bags can be purchased in blue or transparent for different customer preferences. The sponge counter bags are hung on a sponge counter rack in the operating theatre and are used after surgical procedures to quickly count lap sponges.

Item number	Dimensions	Colour	Number of compartments	Packaging
DYNJE1000	25 x 94 cm	Clear	10	50 ea/cs
DYNJE1001	25 x 94 cm	Blue	10	50 ea/cs
NON50511	Sponge counter rack	White	Can hold up to 4 sets of bags	1 ea/cs

## Polyethylene U-drapes



These drapes include courtesy tabs for easy removal.


### Benefits and features

- Provides excellent adhesion without residue
- Facilitates easier application via split, paper liner
- 'U' fenestration with wide adhesive area
- Strong, drapeable, low-glare matte-finish film

**Use:** Prepping procedures: arthroscopies, hip, knee, shoulder, back, cranial, isolation, and applications

Item number	Dimensions	Packaging
IDE1015	121 x 132 cm Adhesive only: 25 x 48 cm	10 ea/bx; 40 ea/cs
IDE1019	121 x 180 cm Adhesive only: 25 x 48 cm 'U' fenestration only: 15 x 48 cm	10 ea/bx; 40 ea/cs

## Shoulder pouch



The capped, closable port eliminates the step of cutting or piercing for drainage and makes it easy to stop the flow at the end of the procedure.

### Benefits and features


- Offers stronger adhesive with no residue
- Oval aperture with slit and adhesive areas
- Provides additional strength due to heavy, low-glare matte-finish film
- Enhanced pouch performance via 2 mouldable strips and 1 foam strip
- Improved fluid control, thanks to 2 suction ports with caps

**Use:** Shoulder arthroscopy, implant applications

Item number	Dimensions	Packaging
IDE1067	91 x 73 cm Adhesive only: 25 x 48 cm Aperture: 19 x 25 cm	24 ea/cs



## U-fluid collection pouch






### Benefits and features

- Enhanced performance via 3 mouldable strips
- 'U' aperture with adhesive area
- Strong adhesive with no residue
- Drapeable, clear polyethylene film

**Use:** Shoulder arthroscopy, hip surgery, orthopaedic

Item number	Dimensions	Packaging
DYJPEFP090	76 x 112 cm Adhesive only: 25 x 51 cm Aperture: 15 x 46 cm	25 ea/cs

## Patient isolation drapes



These drapes have a highly effective filter with deep placement for quick drainage and numerous small perforations.




### Benefits and features

- Better organisation, thanks to 2 cord holders and instrument pouches
- Improved fluid control thanks to large pouch with micro screen, suction port and cap
- Strong incise adhesive with no residue
- Clear, strong polyethylene film

**Use:** Hip pinning, nephroscopies, biliaries

Item number	Dimensions	Packaging
IDE1014	317 x 213 cm Adhesive only: 49 x 25 cm Adhesive strip only: 317 x 1.3 cm	5 ea/bx; 15 ea/cs
IDE1017	317 x 213 cm Adhesive only: 49 x 25 cm Adhesive strip only: 317 x 1.3 cm Pouch only: 51 x 51 cm	15 ea/cs
IDE1017W	317 x 254 cm Adhesive only: 85.1 x 27.9 cm Adhesive strip only: 317 x 1.3 cm Pouch only: 69.9 x 94 cm	15 ea/cs
5801CE	330 x 270 cm Adhesive only: 78 x 26 cm Adhesive strip only: 330 x 1.9 cm Pouch only: 25.5 x 28 cm	16 ea/cs

## Circular/Oval aperture drapes



**Easy-peel liner:** Less force is needed to separate it from the adhesive section.



**Benefits and features**

- Non-linting and tear-resistant
- Paper release liner
- Circular/oval aperture with adhesive area
- Thin, soft, drapeable, low-glare matte-finish film

**Use:** Cataract, corneal, refractive, lid/plastic, minor

Item number	Dimensions	Packaging
IDE1020	40 x 40 cm Adhesive only: 12 x 12 cm Aperture only: 6 cm circle	10 ea/bx; 40 ea/cs
IDE1030	83 x 129 cm Adhesive only: 12 x 11 cm Aperture only: 7 x 5 cm	10 ea/bx; 40 ea/cs
IDE1033	83 x 130 cm Adhesive only: 23 x 11 cm Aperture only: 15 x 5 cm	10 ea/bx; 40 ea/cs

## Ophthalmic aperture pouches



**Benefits and features**

- Incise adhesive area (with oval aperture option) optimised for skin
- Non-linting, tear-resistant, paper release liner
- 2 Adhesive tabs for additional support
- Mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film

**Use:** Cataract, vitreoretinal, corneal

Item number	Dimensions	Packaging
IDE1023	48 x 41 cm Adhesive only: 12 x 10 cm Adhesive tabs only: 5 x 5 cm	10 ea/bx; 40 ea/cs
IDE1024	48 x 41 cm Adhesive only: 6 x 3.5 cm Adhesive tabs only: 5 x 5 cm Aperture only: 6 x 3.5 cm	10 ea/bx; 40 ea/cs

Ophthalmic incise drapes with fluid pouch



Benefits and features


- Non-linting and tear-resistant
- Paper release liner
- Incise adhesive area
- Fluid pouch with mouldable strip for additional fluid control
- Soft, drapeable, low-glare matte-finish film (IDE1061)
- Repellent/low-linting SMS fabric (IDE1062)

Item number	Dimensions	Packaging
IDE1061	121 x 129 cm Adhesive only: 12 x 11 cm Pouch only: 30 x 27 cm	10 ea/bx; 40 ea/cs
IDE1062	122 x 129.5 cm Adhesive only: 13 x 11.5 cm Pouch only: 31 x 28 cm	10 ea/cs





### Aperture drape




#### Benefits and features

- Allows for additional fluid control, thanks to mouldable strip
- Oval aperture with large adhesive area
- Strong adhesive with no residue
- Heavy, low-glare matte-finish film

**Use:** Vaginal, ophthalmic

Item number	Dimensions	Packaging
IDE1072	47 x 23 cm Adhesive only: 30 x 22 cm Aperture only: 11 x 9 cm	10 ea/bx; 40 ea/cs

### Urology drape




#### Benefits and features

- Non-linting and tear-resistant
- Paper release liner
- Strong adhesive with no residue
- Latex-free finger cot
- Strong, low-glare matte-finish film

**Use:** TUR, urology, gynaecology


Item number	Dimensions	Packaging
IDE1071	60 x 91 cm Adhesive only: 60 x 5 cm Finger cot adhesive only: O.D. 10 x I.D. 3 cm	10 ea/bx; 40 ea/cs

### Clear epidural drapes



Item number	Dimensions	Packaging
8461CE	60 x 100 cm, diameter 9 cm Adhesive fenestration	150 ea/cs
8465CE	80 x 140 cm, 6 x 10 cm Adhesive oval fenestration Pre-cut line for easy removal	150 ea/cs

## Cesarean pouch




### Benefits and features

- Provides strong adhesive with no residue
- Promotes improved fluid control via large-cell foam and 2 suction ports with cap
- Complete fluid containment is possible with 360° fluid pouch
- Additional pouch strength via heavy, low-glare matte-finish film
- Line organisation through 2 cord holders
- High-moisture, permeable incise film for extended adhesion

**Use:** C-section, trauma, general

Item number	Dimensions	Packaging
IDE2057	92 x 76 cm Adhesive only: 38 x 36 cm	10 ea/cs

## Wound ring protectors



The wound ring protectors have a flexible but rigid ring at the fenestration site, to help keep the drape in the correct position. The frosted clear plastic drape around the fenestration helps maintain the sterile barrier.


### Benefits and features

- Meets different surgical needs as it features 5 different ring diameters
- Strong, drapeable, low-glare matte-finish film
- Application with adhesive tape

**Use:** General, laparotomy, cardiothoracic, colectomy

Item number	Dimensions	Packaging
IDE3707	90 x 90 cm, ring diameter: 7 cm	40 ea/cs
IDE3712	90 x 90 cm, ring diameter: 12 cm	40 ea/cs
IDE3718	90 x 90 cm, ring diameter: 18 cm	40 ea/cs
IDE3723	90 x 90 cm, ring diameter: 23 cm	40 ea/cs
IDE3728	90 x 90 cm, ring diameter: 28 cm	40 ea/cs





## Class III Medical Devices

### Incise Drapes

Cleaning in preparation for surgery and skin disinfection are necessary steps in the initial reduction of surface bacteria. But with our Dermadine Plus Incise Drape, you are guaranteed to maintain a sterile field, by way of the characteristic broad-spectrum antimicrobial activity against a range of microorganisms present on the skin.




Breathable, water-repellent  
and waterproof



Resistant and anti-reflective for  
good operating comfort



Leaves no residue, thanks to  
strong adhesive




Thin, low-glare polyurethane  
matte-finish film



Impregnated with an antibacterial  
solution



# Dermadine Plus and Dermadrape



Dermadine Plus item number		Dimensions	Packaging	Dermadrape item number	
42020P		20 x 20 cm Adhesive: 15 x 20 cm	20 ea/bx 120 ea/cs	32020	
n/a		35 x 25 cm Adhesive: 30 x 25 cm	20 ea/bx 120 ea/cs	33525	
44040P		40 x 40 cm Adhesive: 34 x 40 cm	15 ea/bx 90 ea/cs	34040	
n/a		50 x 48 cm Adhesive: 44 x 48 cm	15 ea/bx 90 ea/cs	35050	
47045P		70 x 45 cm Adhesive: 60 x 45 cm	15 ea/bx 90 ea/cs	37045	
47060P		70 x 60 cm Adhesive: 60 x 60 cm	10 ea/bx 60 ea/cs	n/a	
47080P		70 x 80 cm Adhesive: 60 x 80 cm	10 ea/bx 60 ea/cs	n/a	
48060P		80 x 60 cm Adhesive: 60 x 60 cm	10 ea/bx 60 ea/cs	38060	
48080P		80 x 80 cm Adhesive: 70 x 80 cm	10 ea/bx 60 ea/cs	38080	

## Equipment Covers

Medline offers a complete range of medical equipment covers that can help provide the right protection, combined with easy handling, to cover most of the equipment in your operating theatre. These include covers for video cameras, C-arms, probes, microscopes and other essential medical equipment.






## Invisishield in the Cath Lab

A catheterisation laboratory (cath lab) is an examination room in a hospital or clinic with diagnostic imaging equipment, used to visualise the arteries and chambers of the heart, and treat stenosis or abnormalities.



## Customised C-arm covers

### Philips C-arm Set



#### Benefits and features


- Promotes visibility even when covered via clear polyethylene material
- Provides complete banded coverage bags for X-ray tank/intensifier
- Provides completed coverage for the C-arm

#### Fits with

- DYNJE4526 for BV29/300
- DYNJE4527 for BV300
- DYNJE4530 for BV25/26
- DYNJE4569 for Pulsera/Endura

Item number	Dimensions	Packaging
DYNJE4526	Banded bags: 60 x 115 cm   70 x 130 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs
DYNJE4527	Banded bags: 60 x 115 cm   rect.: 95 x 85 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs
DYNJE4530	Banded bags: 60 x 115 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs
DYNJE4569	Banded bags: 60 x 115 cm   70 x 140 cm C-arm cover: 30 x 198 cm	10 ea/bx; 20 ea/cs

### Ziehm C-arm Set



#### Benefits and features


- Promotes visibility via clear polyethylene drape material

#### Fits with

- For starter set SPTE-SS634

Item number	Dimensions	Packaging
DYNJE4531	Banded bags: 60 x 115 cm   70 x 130 cm C-arm cover: 35 x 185 cm	10 ea/bx; 20 ea/cs

## Siemens C-arm Set



### Benefits and features


- Promotes visibility via clear polyethylene drape material
- Comes with clips for starter set SPTE-SS628

### Fits with

- DYNJE4532 for SIREMOBIL N4
- DYNJE4556 for ISO-C/Arcadis/Orbic, 2 parts

Item number	Dimensions	Packaging
DYNJE4532	Banded bags: 60 x 115 cm   60 x 115 cm C-arm cover: 35 x 183 cm	10 ea/bx; 20 ea/cs
DYNJE4547	Banded bags: 70 x 130 cm   60 x 115 cm C-arm cover: 45 x 183 cm	20 ea/cs
DYNJE4556	Top: 51 – 82 x 130 cm Bottom: 53 – 46 x 134 cm	10 ea/bx; 20 ea/cs

## Fluoroscope drape



### Benefits and features

- Encourages proper fit via 2 adhesive strips for attachment, with 2 adhesive tabs
- Promotes visibility via clear, drapeable polyethylene film

Item number	Dimensions	Packaging
ICE1012	91 x 111 cm Adhesive only (2): 91 x 1.2 cm Adhesive tabs only: 5 x 10 cm	10 ea/bx; 40 ea/cs

## Full-size C-arm drapes




### Benefits and features

- Provides securement via 3 adhesive polyethylene straps
- Provides addition securement via 2 latex-free elastic bands on ICE4400R
- Clear, drapeable polyethylene film

Item number	Dimensions	Packaging
ICE4400	104 x 188 cm	20 ea/cs
ICE4400R	104 x 188 cm	20 ea/cs
ICE4405	90 x 225 cm	75 ea/cs

## Contoured C-arm drape




### Benefits and features

- Contoured to fit full C-arm
- Offers securement by way of 4 adhesive polyethylene straps
- Clear, drapeable polyethylene film

Item number	Dimensions	Packaging
ICE4430	152 x 267 cm	20 ea/cs

## Mini C-arm drapes



### Benefits and features


- Held in place with an elastic band on the back
- 3 Adhesive polyethylene straps for additional securement
- Clear, drapeable polyethylene film

### Fits with

- ICE4410: fits Fluoroscanner® mini C-arm models
- ICE4420: fits Fluoroscanner® and OEC6600/6800 mini C-arm models

Item number	Dimensions	Packaging
ICE4410	137 x 160 cm	20 ea/cs
ICE4420	137 x 216 cm	20 ea/cs

## Full-size C-arm kits



### Benefits and features

- Available with latex-free bands
- Clear, drapeable polyethylene film

Item number	Dimensions	Packaging
ICE4500	Intensifier drape (1): 91 x 91 cm C-arm drape (1): 41 x 218 cm Footswitch drape (1): 46 x 41 cm Banded bag (1): 71 x 31 cm	20 ea/cs
ICE8130	C-arm drape (1): 10/14 x 225 cm Banded bag (2): 75 x 90 cm	20 ea/cs
ICE8132	C-arm drape (1): 10/14 x 225 cm Banded bag (1): 75 x 90 cm	20 ea/cs
ICE8230	C-arm drape (1): 8/10 x 205 cm Banded bag (2): 75 x 75 cm	20 ea/cs
ICE8232	C-arm drape (1): 8/10 x 205 cm Banded bag (1): 75 x 75 cm	20 ea/cs
ICE8330	C-arm drape (1): 8/10 x 180 cm Banded bag (2): 75 x 75 cm	20 ea/cs
ICE8332	C-arm drape (1): 8/10 x 180 cm Banded bag (1): 75 x 75 cm	20 ea/cs




## Banded bags

## Benefits and features


- Versatile; often used to cover intensifiers and other equipment
- Clear, drapeable polyethylene film
- Latex-free elastic bands

## Circular banded bags



Item number	Dimensions	Stretched width	Packaging
ICE60020S	30 x <b>25</b> cm	<b>55</b> cm	25 ea/cs
ICE60028S	45 x <b>35</b> cm	<b>70</b> cm	25 ea/cs
ICE60030S	56 x <b>40</b> cm	<b>80</b> cm	25 ea/cs
ICE60040S	58 x <b>50</b> cm	<b>100</b> cm	25 ea/cs
ICE60050S	72 x <b>60</b> cm	<b>120</b> cm	25 ea/cs
ICEBBG070	75 x <b>65</b> cm	<b>130</b> cm	50 ea/cs
ICE60060S	140 x <b>70</b> cm	<b>165</b> cm	10 ea/cs


## Rectangular banded bags



Item number	Dimensions	Elastic band- thread color	Packaging
ICE63036R	76 x 91 cm	Brown	25 ea/cs
ICE63628R	91 x 71 cm	White	25 ea/cs
ICE63636R	91 x 91 cm	White	25 ea/cs
ICEBBG095	95 x 85 cm	White	25 ea/cs
ICE63654S	95 x 130 cm	White	25 ea/cs
ICE64848S	122 x 122 cm	Green	25 ea/cs
ICE65656S	142 x 142 cm	White	25 ea/cs

[illegible]

## Banded bags for x-ray shields




Medline's new banded bags help healthcare professionals cover protective shields, which are generally used in angiography procedures. These items provide a sterile barrier between the equipment and the patient undergoing the procedure.

### Benefits and features

- Provides a design suitable for most X-ray panels
- Latex-free elastic bands
- Clear, drapeable polyethylene film

Item number	Diameter 'A' cm	Height 'B' cm	Stretched width 'C' cm	Packaging
ICE8090	30	90	80	25 ea/cs
ICE8090	23	120	90	25 ea/cs

## Cover for thyroid collar



During procedures requiring fluoroscopy, an X-ray protective thyroid collar is worn to shield the wearer from the harmful effects of radiation. Medline's disposable thyroid collars were designed to wrap around and universally fit any thyroid collar shield. Our disposable thyroid collar cover will reduce or prevent the build-up of bacteria, make-up and other debris.

### Benefits and features

- Keeps thyroid collar sanitary and clean
- Comfortable spunlace fabric
- 2 Adhesive tape strips for closing and securement
- Fits standard thyroid shields
- Easy to use and dispose of after use


**Use:** Medical and dental X-rays, orthopaedic, radiology, CT radiation, mammography

Item number	Dimensions	Packaging
SPTE-ICE9000	42 x 21.5 cm	25 ea/bx; 200 ea/cs

Video camera covers

Benefits and features

- Provides stability as it features an adhesive polyethylene strap to hold it in place
- Offers additional securement via an adhesive strip at the base
- Clear, drapeable polyethylene film
- Cable cover



Telescopic with tab

Item number	Dimensions	Packaging
ICE4240	13 x 240 cm	25 ea/cs

ICE4100	8 x 200 cm	25 ea/cs
---------	------------	----------

ICE4245	18 x 250 cm	25 ea/cs
---------	-------------	----------

Concertina

ICE4270	15 x 298 cm	25 ea/cs
---------	-------------	----------

Reverse

ICE4250	13 x 244 cm	25 ea/cs
---------	-------------	----------

ICE4300	23 x 244 cm	25 ea/cs
---------	-------------	----------

With ring

ICE4260	15 x 240 cm	25 ea/cs
---------	-------------	----------

Elastic end

ICE4200	18 x 244 cm	25 ea/cs
---------	-------------	----------

With closed system

ICE4290	14 x 250 cm	25 ea/cs
---------	-------------	----------

ICE4280	18 x 250 cm	25 ea/cs
---------	-------------	----------

Allows for quick telescope changes, saving time in the operation theater

Eliminate the need for several covers




Prevent telescope rotation by way of locking system

Reduces the cost of camera system maintenance , cleaning and inventory

## Microscope covers

### Zeiss OPMI 48mm




#### Benefits and features

- Fits 48-mm lens (Zeiss® OPMI)
- Provides securement via 3 adhesive PE straps
- Fits 1 pair of oculars
- All fitted with clear, drapeable, polyethylene film

Item number	Dimensions	Packaging
ICE5010	51 x 163 cm	20 ea/cs
ICE5040	51 x 267 cm	40 ea/cs

### Zeiss MD 65mm




#### Benefits and features

- Fits 65-mm lens (Zeiss MD)
- Fits 3 pairs of oculars
- 3 Adhesive PE straps for securement (5 for ICE5220)

Item number	Dimensions	Packaging
ICE5200	117 x 203 cm	10 ea/cs
ICE5030	117 x 305 cm	10 ea/cs
ICE5210	117 x 267 cm	10 ea/cs
ICE5220	137 x 381 cm	10 ea/cs



Leica (68mm)




Benefits and features

- Fits 68-mm lens (Leica®)
- Fits 3 pairs of oculars
- 3 Adhesive PE straps for securement

Item number	Dimensions	Packaging
ICE5300	104 x 163 cm	10 ea/cs
ICE5320	137 x 267 cm	10 ea/cs
ICE5330	137 x 381 cm	10 ea/cs

Microscope handle cover



Benefits and features

- Provides complete coverage for microscope handles

Item number	Dimensions	Packaging
ICE5400	15 x 35 cm	20 ea/cs

## Probe covers kits



The high-quality polyurethane material is soft and drapeable, allowing it to conform to the transducer. It also has a smooth finish for easy application and glide. The slightly matte surface helps prevent glare.

### Benefits and features

- 1 Probe cover: soft and drapeable polyurethane film
- 1 Ultrasound gel packet: high-quality, sterile ultrasound gel for transducer
- 2 Elastic bands for securement

Item number	Dimensions	Packaging
<b>Kits with gel</b>		
ICE5900	10 x 122 cm	20 ea/cs
ICE5910	15 x 61 cm	20 ea/cs
ICE5920	15 x 122 cm	20 ea/cs
ICE5930	15 x 244 cm	20 ea/cs
<b>Kits without gel</b>		
ICE18120	18 x 120 cm	25 ea/cs
ICE1830	18 x 30 cm	25 ea/cs
ICE1880	18 x 80 cm	25 ea/cs

## Probe covers kits with adhesive



Medline's SonoGuard Plus Sterile Ultrasound Probe Cover eliminates the need for contact gel use. The cover comes with adhesive ClearView film for attaching the ultrasound transducer head to the inside of the cover. This reduces setup time. And as a result of its non-slip properties, it is optional to use a second elastic band or adhesive strip for additional securement.


### Benefits and features

- 1 Probe cover: clear polyethylene film with adhesive tip
- 1 Ultrasound gel packet; high-quality, sterile ultrasound gel for transducer
- 2 Elastic bands and 2 adhesive tapes for additional securement

Item number	Dimensions	Packaging
<b>Kits with gel</b>		
US.20.180.040-G	18 x 40 cm	50 ea/bx; 200 ea/cs
US.20.180.080-G	18 x 80 cm	50 ea/bx; 200 ea/cs
US.20.180.120-G	18 x 120 cm	50 ea/bx; 200 ea/cs
US.20.180.160-G	18 x 160 cm	50 ea/bx; 200 ea/cs
<b>Kits without gel</b>		
US.20.180.040	18 x 40 cm	100 ea/bx; 200 ea/cs
US.20.180.080	18 x 80 cm	100 ea/bx; 200 ea/cs
US.20.180.120	18 x 120 cm	100 ea/bx; 200 ea/cs
US.20.180.160	18 x 160 cm	100 ea/bx; 200 ea/cs




Traditional ultrasound probe cover



**SonoGuard Plus**  
ultrasound probe cover

Endocavity probe covers



Benefits and features

- Compatible with most ultrasound transducers
- Available as sterile or non-sterile ultrasound covers
- Available in both latex and polyisoprene

Use: Gynaecology, urology and proctology

	Item number	Dimensions	Packaging	Sterile
Latex	NS-ICE28295L	28 x 295 mm	144 ea/bx; 2,880 ea/cs	No
	NS-ICE33230L	33 x 195 mm	144 ea/bx; 2,880 ea/cs	No
	NS-ICE33230LL*	33 x 195 mm	144 ea/bx; 2,880 ea/cs	No
	ICE28295L	28 x 295 mm	100 ea/bx; 2,000 ea/cs	Yes
	ICE33230L	33 x 195 mm	100 ea/bx; 2,000 ea/cs	Yes
Polyisoprene	NS-ICE28295P	28 x 295 mm	100 ea/bx; 2,000 ea/cs	No
	ICE28295P	28 x 295 mm	100 ea/bx; 2,000 ea/cs	Yes

\*With lubricant

## Remote control covers




### Benefits and features

- Promotes visibility via clear polyethylene film
- Bag with pinch-top closure

Item number	Dimensions	Packaging
DYNJE4635	10 x 28 cm	50 ea/bx; 200 ea/cs
DYNJE4641	18 x 35 cm	50 ea/bx; 200 ea/cs

## Thermo cover



### Benefits and features



- Provides temperature resistance up to 60° C
- Clear polyethylene film
- Aluminium coated foil

Draping  
instruction  
video:



Item number	Dimensions	Packaging
DYNJE7000	68 x 86 cm	40 ea/cs

## X-ray cassette covers



### Benefits and features

- Clear, drapeable polyethylene film
- ICE4010: release liner for facilitating sterile technique

Item number	Dimensions	Packaging
ICE4000	Large: 61 x 102 cm	20 ea/cs
ICE4010	Standard with liner: 61 x 64 cm	48 ea/cs
ICE4050	Standard: 61 x 64 cm	50 ea/cs
8720CEA	Large: 54 x 91 cm	48 ea/cs

## Stockinettes

### Impervious Stockinettes




#### Benefits and features

- Offers exceptional, high-level barrier properties and optimal fluid control, thanks to the use of an impervious thermoplastic elastomer material
- Gives the patient comfort by way of a soft, knitted jersey fabric on the inside
- Available in 6 sizes

**Use:** Orthopaedic, arthroplasty, as well as hand, knee, ankle, foot, limb and vascular surgeries

Item number	Width 'A'	Width 'B'	Packaging
22520CE	15 cm	76 cm	20 ea/cs
22525CE	20 cm	97 cm	20 ea/cs
22530CE	23 cm	91 cm	20 ea/cs
22534CE	23 cm	122 cm	16 ea/cs
22535CE	30 cm	122 cm	12 ea/cs
22540CE	36 cm	122 cm	12 ea/cs

### Impervious stretchable stockinettes



#### Benefits and features

- Allows for an ideal fit over the patient's anatomy, due to the use of a stretchable, bilaminate material
- Facilitates set-up with the help of 2 pulling ties
- Both lower and upper extremities can be covered, as it is available in 2 sizes

**Use:** Orthopaedic, arthroplasty, as well as hand, knee, ankle, foot, limb and vascular surgeries

Item number	Width 'A'	Width 'B'	Packaging
22620CE	17 cm	76 cm	18 ea/cs
22634CE	24 cm	118 cm	10 ea/cs



Impervious stretchable  
stockinettes instruction  
video:



## Invisishield Cross-Reference Guide

Page	Description	Product number Medline	Product number 3M®	Product number Microtek®
5	Towel drape	IDE1000	1000	ISO1000
5	Isolation bag	ICE1003	1003	ISO1003
5	Towel drape	IDE1010	1010	ISO1010
9	Patient isolation drape	IDE1014	1014	n/a
9	Polyethylene U-drape	IDE1015	1015	ISO1015
5	Fluid pouch	IDE1016	1016	ISO1016
9	Patient isolation drape	IDE1017	1017	5960
9	Wide patient isolation drape	IDE1017W	n/a	6960
6	Instrument pouch	IDE1018	1018	ISO1018
6	Long instrument pouch	IDE1018L	1018L	n/a
8	Polyethylene U drape	IDE1019	1019	ISO1019
10	Aperture drape	IDE1020	1020	ISO1020
10	Ophthalmic incise drape	IDE1023	1023	n /a
10	Ophthalmic aperture pouch	IDE1024	1024	ISO1024
10	Aperture drape	IDE1030	1030	ISO1030
10	Aperture drape	IDE1033	1033	n/a
15	Incise drape	32020	1035	ISO1035
15	Incise drape	34040	1040	ISO1040
15	Incise drape	37045	1050	ISO1050
15	Incise drape	38080	1051	ISO1051
11	Ophthalmic incise drape with fluid pouch	IDE1061	1061	ISO1061
8	Shoulder pouch	IDE1067	1067	ISO1067
12	Urological drape	IDE1071	1071	ISO1071
12	Aperture pouch	IDE1072	1072	ISO1072
9	U-fenestrated drape with fluid collection pouch	DYJPEFP090	n/a	n/a
6	Instrument pouch	DYJPEISTPOM1	n/a	n/a
6	Instrument pouch	DYJPEISTPOM2	n/a	n/a
6	Instrument pouch	DYJPEISTPOM3	n/a	n/a

## DermaDine Plus Cross-Reference Guide

Type of film	Description	Product number Medline	Product number 3M®
PU	Incise, sterile drape 15 x 20 cm; incise area 10 x 20 cm, Iodophor impregnated	42020P	6635
PU	Incise, sterile drape 44 x 35 cm; incise area 34 x 35 cm, Iodophor impregnated	44040P	6640
PU	Incise, sterile drape 66 x 60 cm; incise area 56 x 60 cm, Iodophor impregnated	47060P	6648
PU	Incise, sterile drape 66 x 45 cm; incise area 56 x 45 cm, Iodophor impregnated	47045P	6650
PU	Incise, sterile drape 66 x 85 cm; incise area 56 x 85 cm, Iodophor impregnated	47080P	6651
PU with applicator	Incise, sterile drape 30 x 20 cm; incise area 26 x 20 cm, Iodophor impregnated	42020P	6661EZ
PU with applicator	Incise, sterile drape 60 x 35 cm; incise area 35 x 35 cm, Iodophor impregnated	44040P	6640EZ
PU with applicator	Incise, sterile drape 90 x 60 cm; incise area 60 x 60 cm, Iodophor impregnated	47060P	6648EZ
PU with applicator	Incise, sterile drape 90 x 45 cm; incise area 60 x 45 cm, Iodophor impregnated	47045P	6650EZ
PU with applicator	Incise, sterile drape 90 x 85 cm; incise area 60 x 85 cm, Iodophor impregnated	48080P	6651EZ

## DermaDrape Cross-Reference Guide

Type of film	Description	Product number Medline	Product number 3M®
PE	Incise, sterile drape 15 x 20 cm; incise area 10 x 20 cm	32020	1035
PE	Incise, sterile drape 29 x 25 cm; incise area 19 x 25 cm	32020	1037
PE	Incise, sterile drape 60 x 35 cm; incise area 35 x 35 cm	34040	1040
PE	Incise, sterile drape 82 x 45 cm; incise area 60 x 45 cm	37045	1050
PE	Incise, sterile drape 82 x 85 cm; incise area 60 x 85 cm	38080	1051
PU	Incise, sterile drape 20 x 20 cm; incise area 10 x 20 cm	32020	2035
PU	Incise, sterile drape 38 x 25 cm; incise area 28 x 25 cm	33525	2037
PU	Incise, sterile drape 38 x 41 cm; incise area 28 x 41 cm	34040	2040
PU	Incise, sterile drape 60 x 45 cm; incise area 50 x 45 cm	35050	2045
PU	Incise, sterile drape 90 x 60 cm; incise area 60 x 60 cm	38060	2050
PU	Incise, sterile drape 90 x 85 cm; incise area 60 x 85 cm	38080	2051

## Stretchable Stockinettes Cross-Reference Guide

Product number Medline	Description	Product number Halyard
22620CE	Impervious, stretchable stockinette; 17 x 76 cm	89811
22634CE	Impervious, stretchable stockinette; 24 x 118 cm	89861




**Medline Industries Ltd**  
3rd Floor  
Quayside Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL  
United Kingdom  
Tel.: +44 844 334 5237  
[www.medline.eu/uk](http://www.medline.eu/uk)  
[uk-customerservice@medline.com](mailto:uk-customerservice@medline.com)

**Medline Ireland Ltd**  
1 Grant's Row  
Mount Street Lower  
Dublin 2  
Ireland  
Tel.: +353 1 691 73 06

[www.medline.eu/ie](http://www.medline.eu/ie)  
[ie-customerservice@medline.com](mailto:ie-customerservice@medline.com)

FOLLOW US  

Isolation bag 586560-1 on page 4 is a class I sterile medical device.  SENGEWALD Klinikprodukte GmbH  
Adlerstr. 2  
D-83101 Rohrdorf-Thansau, Germany

 0123

DermaDrape Plus is a class III medical device, and DermaDrape is a class IIa sterile medical device. 

TIASET di Lorenzo Carresi  
SP loc Botriolo Km 5100  
52026 Castelfranco-Piandisco (AR), Italy

 0373 Istituto Superiore di Sanità

Endocavity probe covers are class IIa medical devices (sterile or non-sterile).  Advena Limited Tower Business  
Centre 2nd Floor, Tower Street,  
Swatara, BKR 4013 Malta

 0123

Wound ring protectors are class IIa sterile medical devices. 

Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France

 2797 (BSI)

C-arm covers ICE4405, ICE81xx, ICE82xx and ICE83xx on pages 19 and 20 and camera covers ICE4280 and ICE4290 are class I sterile medical devices.

 Medbar Tıbbi Malzemeler Turizm San. ve Tic. A.Ş.  
1142 SK. N:35 Sarıı, Gaziantep, Turkey

 1984

Probe cover kits with gel ICE59xx on page 27 are procedure packs according to Article 12 of the MDD 93/42/EEC or Article 22 of the MDR 2017/745.

Procedure pack producer: Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France

Probe cover kits with adhesive on page 27 are class I sterile medical devices intended to be used by healthcare professionals.


Kits without gel:  0537

Promedcon GmbH  
Kreuzerstraße 10  
20095 Hamburg, Germany

Kits with gel: Procedure pack according to Article 12 of the MDD 93/42/EEC

 Probe cover:  
Promedcon GmbH  
Kreuzerstraße 10  
20095 Hamburg, Germany




 Gel:  
Türkmen Medikal  
Akçaburgaz, Muhsin Yazıcıoğlu Cd. No.45/5  
34522 Esenyurt/Istanbul, Turkey

 2292

 PPP Procedure pack producer:  
Promedcon GmbH  
Kreuzerstraße 10  
20095 Hamburg, Germany

All other products listed are class I sterile medical devices, except the cover for thyroid collar on page 23.

 Medline International France SAS  
5 Rue Charles Lindbergh  
44110 Châteaubriant, France

 2797 (BSI)

These products are intended to be used by healthcare professionals.

Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct errors that may occur within this brochure.

© 2023 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML164\_EN 01/2023.